



 **promedics**[®]

Product Catalogue

A wide range of Orthopaedic products



www.promedics.co.uk

**Scottish
NDC Codes inside**

Ordering Information

By Telephone: 01475 746 400

Our customer service department is open from 8.30am to 5.00pm, Monday to Friday. There is a voicemail service available at all other times for orders & queries.

By Fax: 01475 746 410

Fax details of your order, including product codes and a delivery & invoice address.

By Post:

Post details of your order, including product codes and a delivery & invoice address.

Promedics Orthopaedics Ltd.
Port Glasgow Industrial Estate
Gareloch Road, Port Glasgow
Scotland PA14 5XH

By E-mail:

Orders can be e-mailed to:
orders@promedics.co.uk

Online:

See online for our complete range of products, training courses, videos and web shop:
www.promedics.co.uk

Enquiries:

If you have an order or product enquiry, or would like to order a product catalogue, telephone our customer service department on: **01475 746 400**

Delivery (UK mainland only)

We endeavour to deliver all orders within 48 hours, with the majority being despatched on the same day that purchase orders are received.

Delivery dates for made-to-measure or specially requested products will be advised on receipt of an official purchase order. Please note, made-to-measure orders must be received in writing by fax or post and once acknowledged cannot be cancelled and are subject to additional charges.

Promedics reserves the right to make partial deliveries unless otherwise instructed. All orders are subject to a standard post and packaging charge of only £5.95, which is guaranteed next day service. Saturday delivery and timed delivery incur additional costs.

Payment

All prices are exclusive of VAT. If you have an account, payment should be made within 30 days of date of invoice. Payment can be made by credit card, most major cards are accepted. Payment can also be made by Switch or Delta. Cheques should be made payable to 'Promedics Orthopaedics Ltd.'

Return of Goods Policy:

All goods for return must have authorisation. Please call or e-mail our customer service department for a returns note: All returns must be notified within 28 days of receipt of goods. A restocking fee may be applied. All returned goods must be as new and in original packaging.

Contents

Introduction	2 - 4
Hand Therapy	5 - 34
Cervical & Upper Limb	35 - 46
Lower Limb	47 - 64
Ankle & Footwear	65 - 88
Spinal Support	89 - 94
Maternity Support	95 - 98
Fracture Bracing & Splinting	99 - 108
Cold Therapy	109 - 110
Paediatric	111 - 120
Index	121 - 123

Welcome to the Promedics Orthopaedic Product Catalogue

We are very proud to present to you our new catalogue which is packed with new and innovative products. The majority have been designed and manufactured by us in our 40,000 square foot facility in Scotland.

Promedics remain absolutely committed to using Scottish employees to manufacture our products. We employ over 180 people in Scotland to ensure we continue to deliver to you, high quality products that have been designed specifically for the Scottish medical market and your patients needs.

The vast majority of our products are despatched on the same day that we receive your order. Our commitment to manufacturing our goods in Scotland means that we continue to offer unrivalled levels of service to our customers in the UK and we are consistently a Top 5 supplier according to NHS Supply Chains own statistics.

In addition to our own products, Promedics are proud to be the exclusive UK distributor of world leading BREG products from the United States. We offer the same level of customer service for the entire BREG range and these products are available from stock in our Scottish facility.

Your business is very important to us, and for over 30 years we have remained committed to offering exceptional service, quality and educational support to healthcare professionals. Again this year Promedics will run a full training and education workshop programme for over 50 events throughout the UK on a wide range of our orthopaedic, orthotic and therapy products.

If you would like any further information or help, then please get in touch with me, our Customer Service Team, your local Commercial Manager or contact us through our website at www.promedics.co.uk

Yours sincerely

David Baxendale
CEO



promedics

Manufacturing, Partnership & Development

Promedics constantly listen to the needs of our customers, not only in Scotland but also internationally where we sell to over 40 countries. For over 30 years we have invested in Scottish manufacturing, design and production and our facility in Port Glasgow continues to grow and be a leader in European Orthotics and Orthopaedics.

Manufacturing in Scotland not only gives us the ability to react quickly to the changing demands and needs of our customers, but it also allows us to control the raw materials and final quality of our products to ensure that we adhere to the highest possible standards. Our facility is fully accredited to the Medical Device Directive ISO 13485:2003 and is fully audited by British Standards Institute.

We are committed to Research & Development in Scotland. Our in-house team work independently and in-association with our partners and the NHS, to constantly strive to integrate the latest textiles, designs and benefits to ensure our products are designed specifically to the needs of our customers in each country we work in. Everyone is different and it is our job to meet your needs.

Promedics are very proud to represent on an exclusive basis, the world class companies of BREG, North Coast Medical and Orfit.



airstep
walker



Our factory in Port Glasgow

Customer Service

Tel: 01475 746 400

Fax: 01475 746 410

Email: enquiries@promedics.co.uk

twitter  [@promedics](https://twitter.com/promedics)

Promedics operate throughout the UK and are able to supply every hospital, department and Trust in the UK on a next day basis if required. We are a key supplier to NHS Supply Chains in England, Scotland, Northern Ireland and Wales. So whether you prefer to order direct or via contract, we are there to help.

We have a full dedicated Customer Service department to ensure we provide you with outstanding, first class service every time. Whatever your needs our team are available to help, answer your call, and to ensure you get the product you are looking for. We have over 15,000 product lines in stock and our team are fully trained and waiting to help you.

We also have a full national team of local Commercial Managers throughout the UK. To get in touch with yours please contact our Customer Service Team.

Export

Promedics export to over 40 countries and we can offer you a full design and manufacture service to give you a range of products that work for your country.

We recognise that each country has its own unique demands and for many years have been a leading supplier throughout USA, Germany, Scandinavia and Western Europe.

We would be happy to work with you in Great Britain, to design and manufacture the Orthopaedic range that's right for you.

Hand Therapy





 Scotland
NDC Codes

Standard Wrist Brace

Product Code: **SWB**

An 18 cm wrist brace made of quality woven elastic for comfort and durability, enclosing an aluminium palmar bar which can be easily shaped to fit the wrist. The unique overlapped edges offer excellent comfort especially to minimise chafing at the thumb web and palmar crease. Also available in black.

Contains Latex



Size	Left (Beige)	NDC Code Left (Beige)	Right (Beige)	NDC Code Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12	090815	*/22	082209	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13	066094	*/23	040339	*/13/BK	*/23/BK	14 - 15 cm
M	*/15	066100	*/25	066124	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17	066117	*/27	066131	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18	040322	*/28	040360	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code

(Supplied in boxes of 6)



Wrightington Wrist Brace

Product Code: **WWB**

An 18cm wrist brace made of quality woven elastic, with easy open overstitched hook/loop tabs designed for the patient with limited dexterity. The unique overlapped edges offer excellent comfort especially to minimise chafing at the thumb web and palmar crease.

The aluminium palmar bar is adjustable to any wrist position.

Also available in black.

Contains Latex

Size	Left (Beige)	Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12	*/22	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13	*/23	*/13/BK	*/23/BK	14 - 15 cm
M	*/15	*/25	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17	*/27	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18	*/28	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code

(Supplied in boxes of 6)



 Scotland
NDC Codes

Beta Wrist Brace

Product Code: **BWB**

An 18cm wrist brace made of a dual compression elastic which conforms to the wrist. The vinyl pocket encloses an aluminium palmar bar which can be shaped to the wrist. Lycra binding finishes the edges for comfort and prevents chafing at the thumb web.

Also available in beige.



Size	Left (Beige)	NDC Code Left (Beige)	Right (Beige)	NDC Code Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12	-	*/22	-	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13	027071	*/23	027088	*/13/BK	*/23/BK	14 - 15 cm
M	*/15	027057	*/25	027064	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17	027033	*/27	027040	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18	-	*/28	-	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code

(Supplied in boxes of 6)

Nuform Wrist Brace

Product Code: **NUF**

An 18cm wrist brace made of a ventilated elastic for cool comfort. The vinyl pocket encloses an aluminium palmar bar which can be shaped to the wrist.

Lykra binding finishes the splint edges for a comfortable fit and prevents chafing at the thumb web.

Also available in black.



Size	Left (Beige)	Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12	*/22	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13	*/23	*/13/BK	*/23/BK	14 - 15 cm
M	*/15	*/25	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17	*/27	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18	*/28	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code

(Supplied in boxes of 6)

Deltaform Wrist Brace

Product Code: **DFS**

An 18cm woven elastic/cotton support, with an aluminium stay inside a vinyl pocket and a longitudinal seam for excellent conformity to the wrist and forearm.

Also available in beige.

Contains Latex



Size	Left (Beige)	NDC Code Left (Beige)	Right (Beige)	NDC Code Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12	-	*/22	-	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13	021697	*/23	019731	*/13/BK	*/23/BK	14 - 15 cm
M	*/15	019717	*/25	019748	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17	019724	*/27	019755	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18	090792	*/28	090808	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code

(Supplied in boxes of 6)



Deltaform PLUS Wrist Brace

Product Code: **DFP**

A longer, 23cm woven elastic/cotton support, with an aluminium stay inside a vinyl pocket and a longitudinal seam for excellent conformity to the wrist and forearm.

Extends to 2/3 of the forearm for radial/ulnar stability.

Also available in black.

Contains Latex



Size	Left (Beige)	NDC Code Left (Beige)	Right (Beige)	NDC Code Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12	-	*/22	-	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13	027149	*/23	-	*/13/BK	*/23/BK	14 - 15 cm
M	*/15	027163	*/25	027156	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17	027187	*/27	027170	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18	027200	*/28	027194	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code





Freedom Wrist Brace

Product Code: **FWB**

An 18cm long brace with cut-away design to allow flexion at the palmar crease and opposition of the thumb to the fingers. The soft Lycra binding prevents chafing of the thumb web.

Longitudinal seam gives good conformity and support to the wrist and forearm.

Also available in black.

Contains Latex

Size	Left (Beige)	Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12	*/22	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13	*/23	*/13/BK	*/23/BK	14 - 15 cm
M	*/15	*/25	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17	*/27	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18	*/28	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code

(Supplied in boxes of 6)



 Freedom

Freedom SX Wrist Brace

Product Code: **FSX**

Using latex free elastic, the Freedom SX brace is manufactured with a full length seam assembly that shapes the brace to fit the wrist and forearm perfectly. This newly developed elastic feels soft and kind to the skin. 18cm long with reinforced casing and aluminium palmar bar.

Available in grey with black (Pictured) Also available in black.

Size	Left (Grey)	Right (Grey)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12/GB	*/22/GB	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13/GB	*/23/GB	*/13/BK	*/23/BK	14 - 15 cm
M	*/15/GB	*/25/GB	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17/GB	*/27/GB	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18/GB	*/28/GB	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code

(Supplied in boxes of 6)



 JURA Medical

Jura Latex Free Wrist Brace

Product Code: **JLFW**

A latex free elastic brace with horizontal wrist seam. This, together with tri-weave compression elastic, provides improved fit and compression. 18cm long with vinyl casing and aluminium palmar bar

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code

poroflex

Poroflex is a new three dimensional fabric offering excellent properties in the design of orthopaedic braces. Poroflex offers improved patient comfort due to breathability and temperature control.



Poroflex Lace

Product Code: **PFLB**

- Simple "slip-on" design with one strap closure. Excellent for patients with limited dexterity
- Adjustable aluminium palmar bar with low profile dorsal stays
- Excellent wrist compression and stability

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Poroflex Pull

Product Code: **PFWB**

- Simple "slip-on" design with two strap D ring closures. Excellent for patients with limited dexterity
- Adjustable aluminium palmar bar with flexible dorsal stays
- Excellent wrist compression and stability

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code




Poroflex Pull & Thumb

Product Code: **PFTB**

- Simple "slip-on" design with two strap D ring closures. Excellent for patients with limited dexterity
- Thumb extension for MCP & CMC control
- Adjustable aluminium thumb stay
- Adjustable aluminium palmar bar with flexible dorsal stays
- Excellent wrist compression and stability

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Lycrafleece™

Lycrafleece™ is a Lycra® Spandex™ 3 ply laminate, designed for the ultimate in patient comfort. The unique construction makes this material perfect for orthopaedic bracing. It is our most comfortable and durable material which makes this range of products the elite in medical bracing.

Lycrafleece™ has unique properties which keep the patient warm when it's cold and cool when it's warm for optimum comfort without sacrificing performance.

Lycra® Spandex™ gives the material unique stretch and recovery properties whilst a polyester lining material wicks moisture away from the skin to allow greater patient compliance.

Lycrafleece Thumb Spica

Product Code: **BLTS**

Supports the wrist and thumb, leaving the palm completely free. An 11cm removable flexible stay restricts movement at the CMC and MCP joints, leaving the IP free. The angled hook/loop attachment at the wrist and thumb ensures a perfect fit. Use for patients with Rheumatoid and Osteoarthritis.



Lycrafleece™

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code

Lycrafleece Wrist Brace

Product Code: **BLWB**

An 18cm three-piece wrist brace featuring longitudinal seams to give greater conformity at the wrist. The malleable palmar bar has a spoon shaped end which conforms to the palm of the hand for a comfortable fit.



Lycrafleece™

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Pro-rheuma Wrist Brace

Product Codes: **PRHW** (Beige) **PRHWB** (Black)

Made from a breathable high compression material with microfleece lining. Easy 'slip-on' design with D-ring closure straps. This brace is very comfortable and excellent for patients with limited dexterity, such as osteoarthritis and wrist instability. Supplied with medium stability plastic palmar support which can be heat moulded.

Available in black & beige.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	13 - 14 cm
S	*/13	*/23	14 - 16 cm
M	*/15	*/25	16 - 19 cm
L	*/17	*/27	19 - 20 cm
XL	*/18	*/28	20 - 22 cm
XXL	*/19	*/29	22 - 24 cm

* Insert product code



Pro-rheuma Wrist Thumb Brace

Product Code: **PRHT** (Beige) **PRHTB** (Black)

Made from a breathable high compression material with microfleece lining. Easy 'slip-on' design with D-ring closure straps. This brace is very comfortable and excellent for patients with limited dexterity, such as osteoarthritis and wrist instability. Supplied with medium stability plastic palmar support which can be heat moulded. Semi-flexible stay provides thumb abduction and support to CMC and MCP joints. Available in black & beige.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	13 - 14 cm
S	*/13	*/23	14 - 16 cm
M	*/15	*/25	16 - 19 cm
L	*/17	*/27	19 - 20 cm
XL	*/18	*/28	20 - 22 cm
XXL	*/19	*/29	22 - 24 cm

* Insert product code



Economy Standard Neobrace

Product Code: **ESNB**

An 18cm wrist brace made in Neoprene which retains heat to soothe painful joints, the brace conforms well and can be cut to fit swollen arthritic joints without fraying.

Made in thinner Neoprene for light support.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Neobrace Wrist Brace

Product Code: **SNB** (Standard) **LNB** (Long)

Made in Neoprene which retains heat to soothe painful joints, this brace conforms well and can be cut to fit swollen, arthritic joints without fraying. Standard: 18cm length. Long: 23cm length.

Size	Left	NDC Code	Right	NDC Code	Wrist Circumference
XS	SNB/12	-	SNB/22	-	11.5 - 13 cm
S	SNB/13	074174	SNB/23	052592	14 - 15 cm
M	SNB/15	052578	SNB/25	052585	16.5 - 18 cm
L	SNB/17	052653	SNB/27	052646	20 - 21.5 cm
XL	SNB/18	052707	SNB/28	052608	23 - 24 cm
XXL	SNB/19	-	SNB/29	-	25 - 28 cm

(SNB Supplied in boxes of 5)

Size	Left	NDC Code	Right	NDC Code	Wrist Circumference
XS	LNB/12	-	LNB/22	-	11.5 - 13 cm
S	LNB/13	106738	LNB/23	106776	14 - 15 cm
M	LNB/15	106745	LNB/25	106783	16.5 - 18 cm
L	LNB/17	106752	LNB/27	106790	20 - 21.5 cm
XL	LNB/18	106769	LNB/28	106806	23 - 24 cm





Air X Wrist Thumb

Product Code: **AWT**

A two-piece wrist thumb brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Aluminium palmar bar and thumb stay provides excellent support to wrist and first CMC & MCP joints. The thumb stay can be easily adjusted for the required degree of abduction and the thumb closure is easily adjusted for a comfortable fit

- Breathable
- Vapour transfer for reduced moisture
- Skin Friendly
- Comfort, compression & Conforms

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Air X Thumb Restriction Splint


Product Code: **ATRS**

A low profile AirX wrap that allows uninterrupted finger function and provides adjustable support to the thumb CMC joint. Excellent in the treatment of osteoarthritis and medial subluxation. The adjustable strap can be wrapped around the first CMC and MCP joints to provide targeted joint support.

Size	Left	Right	Wrist Circumference
Infant	*/11	*/21	11 - 14 cm
Child	*/12	*/22	14 - 17 cm
S	*/13	*/23	15 - 18 cm
S+	*/13+	*/23+	17 - 18 cm
M	*/15	*/25	18 - 20 cm
M+	*/15+	*/25+	20 - 21 cm
L	*/17	*/27	20 - 23 cm
L+	*/17+	*/27+	21 - 24 cm
XL	*/18	*/28	23 - 25 cm

* Insert product code

- Breathable
- Comfortable
- Skin Friendly



Moisture-transporting material

The unique AirX™ fabric has the ability to transport moisture efficiently away from the body. The layering of fibre structures promotes a moisture wicking process, which results in a dry comfortable feeling to the surface of the skin so that body always feels fresh and dry. AirX™ gives the body optimal freedom to breathe – reduces the risk of skin maceration which results in improved patient compliance.

A pleasant surprise

It is hard not to be impressed. AirX™ is made of a soft, lightweight material that perfectly follows the contours of the body. This means that AirX™ adapts itself perfectly to the body so that it feels extremely comfortable to wear – even for long periods at a time – a feature, which is often a problem with orthopaedic and sport medical braces and supports. In this way, AirX™ is a comfortable material that helps to increased patient compliance without you even knowing!

AirX™
Let your body breathe.



Air X One Piece Wrist

Product Code: **AWB**

A one piece wrist brace utilising the excellent properties of AirX to provide a comfortable well fitting brace. The textile is very durable and, together with the reinforced bar casing, makes this brace an excellent choice for patients using braces longterm.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Air X Two Piece Wrist

Product Code: **AFS**

This two piece brace is shaped to achieve excellent fit, compression and stability around the hand, wrist and forearm. The 3D structure of the AirX textile gives a high degree of cushioning and increases comfort and compression. It also allows air transfer making the brace comfortable to wear for extended periods

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Air X Thumb Spica

Product Code: **ATS**

Adjustable thumb Spica offering support around the wrist, MCP and CMC joints. AirX offers improved properties

- Breathable
- Vapour transfer for reduced moisture
- Skin Friendly
- Comfort, compression & Conforms

Size	Left	Right	Wrist Circumference
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	15 - 16.5 cm
L	*/17	*/27	16.5 - 18 cm
XL	*/18	*/28	18 - 19 cm

* Insert product code



Scotland
NDC Codes

Jura Wrist Brace

Product Code: **JWB**

Dual compression elastic for improved support and fit. A two-piece brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Hard wearing webbing bar casing with aluminium palmar bar. Hook/loop fastening

Size	Left	NDC Code	Right	NDC Code	Wrist Circumference
XS	*/12	182794	*/22	182800	11.5 - 13 cm
S	*/13	182848	*/23	182855	14 - 15 cm
M	*/15	182824	*/25	182831	16.5 - 18 cm
L	*/17	166084	*/27	182817	20 - 21.5 cm
XL	*/18	182770	*/28	182787	23 - 24 cm

* Insert product code

Jura Wrist Thumb

Product Code: **JWT**

Dual compression elastic for improved support and fit. A two-piece wrist thumb brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Aluminium palmar bar and thumb stay provides excellent support to wrist and first CMC & MCP joints. The thumb stay can be easily adjusted for the required degree of abduction and the thumb closure is easily adjusted for a comfortable fit.

Size	Left	NDC Code	Right	NDC Code	Wrist Circumference
XS	*/12	-	*/22	-	11.5 - 13 cm
S	*/13	165766	*/23	165810	14 - 15 cm
M	*/15	165773	*/25	165827	16.5 - 18 cm
L	*/17	165780	*/27	165834	20 - 21.5 cm
XL	*/18	165797	*/28	165841	23 - 24 cm

* Insert product code



Scotland
NDC Codes

Jura Wrist Thumb (Long)

Product Code: **JWTL**

Similar design and function as the standard Jura Wrist Thumb Brace. Increasing length to 23cm offers improved wrist stability and immobilisation.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code





Thumbs Up

Product Code: **TUJ**

Lightweight structure and comfortable fitting thumb spica. Innovative design with breathable holes. Protection glove and contact closure make this thumb spica an easy to use, well fitting solution for thumb care. Anatomically shaped, provides the right immobilisation to relieve pain during everyday activities.

Contains Latex

Size	Left	Right	Wrist Circumference
S	*/13	*/23	13 - 17 cm
M	*/15	*/25	17 - 19 cm
L	*/17	*/27	19 - 23 cm
XL	*/18	*/28	23 - 28 cm

* Insert product code



Jura Palm Protector

Product Code: **JPC**

Fleece lined lightweight palm cover. Designed to allow full finger movement whilst reducing finger contractures and skin damage/breakdown in the palmar region. The Jura Palm Cover is easily laundered and is easily applied using hook/loop fastening. One size fits all.

Code (Left)	Code (Right)
JPC/L	JPC/R



Intrinsic Wrist Brace

Product Code: **IWB**

Padded elastic wrist brace with a full-length adjustable aluminium resting splint. The splint is wide enough to contain all four fingers and can be easily adjusted to hold the wrist and fingers in the desired position. There is a wide, padded adjustable dorsal closure strap that can be used to hold the fingers in position.



Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Wrightington Neobrace

Product Code: **WNB**

An 18cm wrist brace with easy-open overstitched hook/loop tabs designed for the patient with limited dexterity.

Made in Neoprene which retains heat to soothe painful joints, it conforms well and can be cut to fit swollen arthritic joints without fraying.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code (Supplied in boxes of 5)



Ventoprene Neobrace

Product Code: **VNB**

An 18cm wrist brace made in ventilated Neoprene with a towelling lining to reduce perspiration build-up. Incorporating a malleable aluminium palmar bar which can be shaped to fit the wrist.

Size	Left (Beige)	Right (Beige)	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code (Supplied in boxes of 5)

Size	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12/BK	*/22/BK	11.5 - 13 cm
S	*/13/BK	*/23/BK	14 - 15 cm
M	*/15/BK	*/25/BK	16.5 - 18 cm
L	*/17/BK	*/27/BK	20 - 21.5 cm
XL	*/18/BK	*/28/BK	23 - 24 cm

* Insert product code



Ventoprene Wrist Brace

Product Code: **VWNB (Beige) R/VWNB (Blue)**

An 18cm wrist brace with Wrightington easy-open, overstitched hook/loop tabs designed for the patient with limited dexterity.

Made in beige or blue ventilated Neoprene with a towelling lining to reduce heat build-up and perspiration. Incorporating a malleable aluminium palmar bar which can be shaped to fit the wrist.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code (Supplied in boxes of 5)

Wrist/Ulnar Deviation Support

Product Code: **UDSA**

A unique functional splint to correct early MCP Ulnar Deviation (U.D.) deformity and provide wrist support for patients with Rheumatoid arthritis. Made in ventilated Neoprene with a towelling to reduce heat build up and reduce perspiration. Incorporating a malleable aluminium palmar bar which can be shaped to fit the wrist.



Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code

Wristlet

Product Code: **WL**

A 15cm wrist support made in Neoprene which retains heat to soothe painful joints. It conforms well giving light compression and can be cut without fraying, to fit swollen joints.



Size	Left	Right	Wrist Circumference
S-M	*/14	*/24	14 - 18 cm
M-L	*/16	*/26	17 - 22 cm
L-XL	*/LXL/L	*/LXL/R	23 - 24 cm

* Insert product code

Elastic Wrist Thumb Brace

Product Code: **WTE**

An elastic wrist brace incorporating a padded malleable aluminium thumb stay.

Holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint.

Contains Latex



Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Neoprene Wrist Thumb Brace

Product Code: **WTO**

A malleable aluminium thumb stay holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint.

Size	Left	NDC Code	Right	NDC Code	Wrist Circumference
XS	*/12	027231	*/22	082186	11.5 - 13 cm
S	*/13	027279	*/23	027262	14 - 15 cm
M	*/15	027255	*/25	027248	16.5 - 18 cm
L	*/17	027293	*/27	027286	20 - 21.5 cm
XL	*/18	027316	*/28	027309	23 - 24 cm
XXL	*/19	-	*/29	-	24 cm +

* Insert product code

Wrist Thumb Brace

Product Code: **WTN**

A Neoprene wrist brace incorporating a padded, malleable aluminium thumb stay. Holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint. The thumb section opens fully for easy application and to accommodate swelling.



Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code

Bedford Wrist Thumb Brace

Product Code: **BWTB**

A unique wrist immobilisation brace which can be used to support undisplaced fractures of the scaphoid, trapezium and thumb metacarpal. It features a mouldable bar which holds the thumb firmly in abduction and has been improved to reduce pressure points. It also now features an open thumb seam for ease of application and comfort.



Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code

Abducted Thumb Wrap

Product Code: **ATW**

A Neoprene wrist/thumb wrap with a hook/loop closure, incorporating a flexible stay to hold the thumb firmly in abduction.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 22 cm

* Insert product code



Wrist Thumb Wrap

Product Code: **WTW**

Holds the thumb lightly in abduction and stabilises the wrist with a double wrist wrap in Neoprene.

Universal design fits left and right.

Size	Code	Wrist Circumference
S	*/03	14 - 15 cm
M	*/05	16.5 - 18 cm
L	*/07	20 - 22 cm
XL	*/08	22 - 24 cm

* Insert product code



Wrist Wrap

Product Code: **WW**

The double wrist wrap in Neoprene provides warmth and compression to soothe arthritic joints.

Universal design fits left and right.

Size	Code	Wrist Circumference
One Size	WW/10	13 - 21.5 cm





Elastic Thumb Spica

Product Codes: **ETS** (Beige) **ETSG** (Blue with green trim)

Dual compression beige elastic supports the wrist and thumb, leaving the palm free. The 11cm removable flexible thumb stay restricts movement at the CMC and MCP joints, leaving the IP joint free.

Size	Left	Right	Wrist Circumference
XS	ETS/12	ETS/22	11 - 13 cm
S	ETS/13	ETS/23	14 - 15 cm
M	ETS/15	ETS/25	15 - 16.5 cm
L	ETS/17	ETS/27	16.5 - 18 cm
XL	ETS/18	ETS/28	18 - 19 cm

Size	Left	Right	Wrist Circumference
XS	ETSG/12	ETSG/22	11 - 13 cm
S	ETSG/13	ETSG/23	14 - 15 cm
M	ETSG/15	ETSG/25	15 - 16.5 cm
L	ETSG/17	ETSG/27	16.5 - 18 cm
XL	ETSG/18	ETSG/28	18 - 19 cm

Contains Latex



Elastic Thumb Spica - Modified

Product Code: **ETSM**

One-way stretch elastic supports the wrist and thumb, leaving the palm completely free. An 10cm removable flexible thumb stay restricts movement at the CMC and MCP joints, leaving the IP joint free. Angled hook/loop attachment at wrist and thumb ensures perfect fit.

Contains Latex

Size	Left	NDC Code	Right	NDC Code	Wrist Circumference
XS	*/12	-	*/22	-	11 - 13 cm
S	*/13	052714	*/23	052721	14 - 15 cm
M	*/15	052660	*/25	052677	16 - 17 cm
L	*/17	052684	*/27	052691	17 - 18 cm
XL	*/18	126682	*/28	126699	18 - 19 cm

* Insert product code



Promedics Thumb Shield

Product Code: **TS**

A Neoprene wrist/thumb wrap with adjustable cinch strap and removable thermoplastic insert. The cinch strap allows greater compression at the site of the pain, without complete immobilisation.

The thermoplastic insert, made from Orfit NS 2.0mm, increases the support at the MCP joint and can be removed and remoulded as required.

Size	Code	Wrist Circumference
XS	*/02	11 - 13 cm
S	*/03	14 - 16 cm
M	*/05	17 - 19 cm
L	*/07	20 - 22 cm
XL	*/08	23 - 25 cm

* Insert product code



Pollo Form Spica - Standard

Product Code: **SPF**

The unique two-way stretch material gives the Pollo Form excellent properties for supporting the wrist and thumb. Outer material is durable and supports an inner pocket with a pre-cut low temperature thermoplastic insert. The insert is activated by pouring warm water into the pocket. After removing the water the brace is then fitted to the patient whilst forming the insert to hold the thumb in the required position.

The insert can be reactivated to allow serial adjustments or removed later in the treatment programme. The Long version gives greater control around the wrist.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	Up to 14 cm
S	*/13	*/23	14 - 16 cm
M	*/15	*/25	16 - 19 cm
L	*/17	*/27	19 - 21 cm
XL	*/18	*/28	21 - 24 cm

* Insert product code



Pollo Form Spica - Long

Product Code: **LPF**

Details above.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	Up to 14 cm
S	*/13	*/23	14 - 16 cm
M	*/15	*/25	16 - 19 cm
L	*/17	*/27	19 - 21 cm
XL	*/18	*/28	21 - 24 cm

* Insert product code



Pro-rheuma Thumb Spica

Product Code: **PRHS**

Made from a breathable high compression material with microfleece lining. Removable semi-flexible thumb stay helps stabilise CMC & MCP joints. Simple hook and loop closures.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	13 - 14 cm
S	*/13	*/23	14 - 16 cm
M	*/15	*/25	16 - 18 cm
L	*/17	*/27	18 - 20 cm
XL	*/18	*/28	20 - 22 cm

* Insert product code



Neoprene Thumb Spica

Product Code: **NTS**

A Neoprene thumb spica featuring a removable flexible stay which restricts movement at the CMC and MCP joints, leaving the IP joint free.

Size	Left	Right	Wrist Circumference
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	15 - 16.5 cm
L	*/17	*/27	16.5 - 18 cm
XL	*/18	*/28	18 - 19 cm

* Insert product code



Jura Ulnar Deviation Support

Product Code: **JUDS**

Helps support the MCP joints and reduce ulnar drift whilst allowing functional use of the fingers. Manufactured from neoprene to provide warmth and support

Size	Left	Right	Width at MCP Joints
S - M	*/14	*/24	6 - 9 cm
M - L	*/16	*/26	9 - 11 cm

* Insert product code



Juraprene Wrist Thumb Wrap

Product Code: **JTWWS**

A latex free neoprene with terry plush lining makes this a very comfortable support. Full finger flexion whilst giving good circumferential compression and support to the wrist and thumb.

The wrap is easily fastened using hook and loop closures.

Fits either left or right hand and can be trimmed in areas for a custom fit. Available in Beige, Black and Blue.

Size	Colour	Code	Wrist Circumference
S	Beige	*/03/BE	Up to 17 cm
M	Beige	*/05/BE	17 - 19 cm
L	Beige	*/07/BE	19 cm +
S	Blue	*/03/BL	Up to 17 cm
M	Blue	*/05/BL	17 - 19 cm
L	Blue	*/07/BL	19 cm +
S	Black	*/03/BK	Up to 17 cm
M	Black	*/05/BK	17 - 19 cm
L	Black	*/07/BK	19 cm +

* Insert product code



Fits either left or right hand and can be trimmed in areas for a custom fit. Available in Beige, Black and Blue.

Juraprene Wrist Thumb Wrap - Long

Product Code: **JTWWL**

A latex free neoprene with terry plush lining makes this a very comfortable support. Full finger flexion whilst giving good circumferential compression and support to the wrist and thumb. This longer wrist wrap adds to the compression and firmness around the wrist and is easily fastened using hook and loop closures.

Size	Colour	Code	Wrist Circumference
S	Beige	*/03/BE	Up to 17 cm
M	Beige	*/05/BE	17 - 19 cm
L	Beige	*/07/BE	19 cm +
S	Blue	*/03/BL	Up to 17 cm
M	Blue	*/05/BL	17 - 19 cm
L	Blue	*/07/BL	19 cm +
S	Black	*/03/BK	Up to 17 cm
M	Black	*/05/BK	17 - 19 cm
L	Black	*/07/BK	19 cm +

* Insert product code



Juraprene Short Wrist Wrap - Universal

Product Code: **JSWWU**

A latex free neoprene with terry plush lining makes this a very comfortable universal support. Full finger and thumb flexion whilst giving good circumferential compression and support to the wrist.

Size	Colour	Code
Universal	Beige	*/10/BE
Universal	Blue	*/10/BL
Universal	Black	*/10/BK

* Insert product code



Juraprene Long Wrist Wrap - Sizes


Product Code: **JLWW**

A latex free neoprene with terry plush lining makes this a very comfortable support. Full finger and thumb flexion whilst wider reinforced hook and loop closure gives improved circumferential compression and support to the wrist.

Available in Beige, Black and Blue.

Size	Colour	Code	Wrist Circumference
S	Beige	*/03/BE	Up to 17 cm
M	Beige	*/05/BE	17 - 19 cm
L	Beige	*/07/BE	19 cm +
S	Blue	*/03/BL	Up to 17 cm
M	Blue	*/05/BL	17 - 19 cm
L	Blue	*/07/BL	19 cm +
S	Black	*/03/BK	Up to 17 cm
M	Black	*/05/BK	17 - 19 cm
L	Black	*/07/BK	19 cm +

* Insert product code



Procool

Procool is a range of lightweight and breathable splints producing excellent compression and support for a wide range of conditions and injuries.

Made with tiny perforations in the Neoprene for comfortable wear and a terry cloth lining that feels good against the skin.

Procool offers you just the fit you want - excellent stability without excessive bulk.



Procool

This brace can be trimmed if required. Available in black and beige.

Indications Include:

- CMC Joint Arthritis
- CMC Instability
- Tendinitis, Pain associated with repetitive motion.

Procool Thumb Restriction Splint

Product Codes: PTRS (Black) PTRSB (Beige)

Whilst providing full finger function, this very low profile wrap around support provides adjustable support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support.

An extensive size range covers paediatric to adult sizes and plus sizes are a full size bigger in thumb circumference to fit oedematous thumb conditions. Function and fit without bulk.

Size	Left	Right	Wrist Circumference
Infant	*/11	*/21	11 - 14 cm
Child	*/12	*/22	14 - 17 cm
S	*/13	*/23	15 - 18 cm
S+	*/13+	*/23+	17 - 18 cm
M	*/15	*/25	18 - 20 cm
M+	*/15+	*/25+	20 - 21 cm
L	*/17	*/27	20 - 23 cm
L+	*/17+	*/27+	21 - 24 cm
XL	*/18	*/28	23 - 25 cm

Infant & child sizes not available in beige

* Insert product code



Procool

Indications Include:

- Arthritic conditions
- Post reconstructive surgery
- Sprains / strains
- Tendinitis
- Post wrist fracture

Procool Wrist Support

Product Code: PWS

This support provides effective circumferential compression and support to the wrist whilst allowing full thumb and finger motion.

Useful as a step-down splint in later stage rehab or to encourage function with patients who struggle with transition from immobilisation to movement due to pain or anxiety.

The brace can be trimmed if required.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



 Procool

This brace can be trimmed if required.

Indications Include:

- Sprains
- Arthritis
- Carpal Tunnel Syndrome
- Strains

Procool D-ring Wrist Splint

Product Codes: **PDR** (Standard) **PDRF** (Firm)

Low profile support with circumferential slip-on design and D-Ring fastening to allow easier one handed application. Aluminium palmar bar can be removed for further adjustment. Dorsal stockinette lining facilitates application and cushions strap closures. Design of the Standard and Firm braces allow full finger and thumb function whilst limiting wrist flexion.

The Firm splint has the addition of adjustable dorsal stays to reduce unwanted flexion and extension and improve overall immobilisation of acute soft tissue injuries.

Size	Left (Standard)	Right (Standard)	Left (Firm)	Right (Firm)	Wrist Circumference
XS	*/12	*/22	*/12	*/22	> 15 cm
S	*/13	*/23	*/13	*/23	15 - 17 cm
M	*/15	*/25	*/15	*/25	17 - 19 cm
L	*/17	*/27	*/17	*/27	19 - 22 cm
XL	*/18	*/28	*/18	*/28	22 cm +

* Insert product code



 Procool

Indications Include:

- Post scaphoid fracture
- De Quervian's
- Post surgical repair
- Carpal Tunnel Syndrome
- Skier's Thumb
- Tendinitis

Procool Deluxe Wrist Thumb Brace

Product Code: **PDWT**

A uniquely shaped wrist thumb brace offering excellent support and compression to the wrist, CMC and MCP joints.

Malleable palmar bar positions wrist in functional extension and flexible stay provides abduction and limited motion to support and protect CMC & MCP joints. Adjustable thumb closure allows easy application of the brace.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	12 - 14 cm
S	*/13	*/23	14 - 16 cm
M	*/15	*/25	16 - 18 cm
L	*/17	*/27	18 - 20 cm
XL	*/18	*/28	20 - 22 cm

* Insert product code



 Procool

Indications Include:

- Post scaphoid fracture
- De Quervian's
- Post surgical repair
- Carpal Tunnel Syndrome
- Skier's Thumb
- Tendinitis

Procool Wrist Thumb Restriction

Product Code: **PWTRS**

This extended version of the thumb restriction splint provides additional wrist support for a combination of wrist and thumb conditions. A wide 'pull through' wrist strap allows easier brace application, adjustment and fastening, offering extra compression and support at the wrist. Function and fit without bulk.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	12 - 14 cm
S	*/13	*/23	15 - 18 cm
M	*/15	*/25	18 - 20 cm
L	*/17	*/27	20 - 24 cm
XL	*/18	*/28	23 - 25 cm

* Insert product code

Black Wrist Brace

Product Code: **JRWB**

The Jura Black Wrist Brace immobilises and protects the wrist following tendonitis and other soft tissue injuries such as sprains/strains, and post-trauma healing. Lightweight, breathable and contoured for a comfortable fit. Slide-on design with D-ring closures allow easy application.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	> 15 cm
S	*/13	*/23	15 - 17 cm
M	*/15	*/25	17 - 19 cm
L	*/17	*/27	19 - 22 cm
XL	*/18	*/28	23 - 25 cm

* Insert product code



NEW



Adjustable Palmar Bar

Jura Black Wrist / Thumb Brace

Product Code: **JRWTB**

The Jura Black Wrist/Thumb Brace immobilises and protects the wrist and thumb following tendonitis and other soft tissue injuries such as sprains/strains, and post-trauma healing. Lightweight, breathable and contoured for a comfortable fit. Slide-on design with D-ring closures allow easy application.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	> 15 cm
S	*/13	*/23	15 - 17 cm
M	*/15	*/25	17 - 19 cm
L	*/17	*/27	19 - 22 cm
XL	*/18	*/28	23 - 25 cm

* Insert product code



Breathable Material





NEW

Ventilated Thumb Restriction Splint

Product Codes: **TRSV**

A ventilated, low profile wrap around support providing adjustable support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support.

An extensive size range covers paediatric to adult sizes and plus sizes are a full size bigger in thumb circumference to fit oedematous thumb conditions. Function and fit without bulk.

Size	Left	Right	Wrist Circumference
XXS	*/11	*/21	11 - 14 cm
XS	*/12	*/22	14 - 17 cm
S	*/13	*/23	15 - 18 cm
S+	*/13+	*/23+	17 - 18 cm
M	*/15	*/25	18 - 20 cm
M+	*/15+	*/25+	20 - 21 cm
L	*/17	*/27	20 - 23 cm
L+	*/17+	*/27+	21 - 24 cm
XL	*/18	*/28	23 - 25 cm

* Insert product code



NEW

Reversible Wrist

Product Codes: **RWB**

A universal left or right design. Easily converted by moving the aluminium palmar bar. The unique elastic design offers excellent compression and conforms well around the hand and wrist.

The additional elastic wrist wrap offers improved compression and wrist stabilisation.

Contains Latex

Size	Code	Wrist Circumference
XS	*/02	11.5 - 13 cm
S	*/03	14 - 15 cm
M	*/05	16.5 - 18 cm
L	*/07	20 - 21.5 cm
XL	*/08	23 - 24 cm

* Insert product code





NEW

Jura Universal Wrist Brace

Product Codes: **JUWB**

A universal support that is easily adjusted for a good fit. Made from a comfortable ventilated material that offers compression and support.

- Sprains / strains
- Carpal Tunnel
- Post Trauma
- Arthritic pain

Size	Code
Left	JUWB/L
Right	JUWB/R



NEW

Mini Thumb Wrap

Product Code: **MTW**

A one piece low profile wrap offering mild compression and support to the thumb CMC and MCP.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11 -13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16 - 17 cm
L	*/17	*/27	17 - 18 cm
XL	*/18	*/28	18 - 19 cm

* Insert product code



NEW

Two Piece Neoprene Spica

Product Code: **TPNS**

A unique two piece construction offers excellent support and fit around the wrist and thumb.

Mouldable stay provides support to the thumb CMC and MCP.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11 -13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16 - 17 cm
L	*/17	*/27	17 - 18 cm
XL	*/18	*/28	18 - 19 cm

* Insert product code



NEW

Neobrace PLUS Wrist Brace

Product Code: **NNWB**

This enhanced version of the Standard Neobrace offers easy grip-tabs to the hook closures. Bound edges and a hard wearing palmar bar casing ensures long term wear. Adjustable aluminium palmar bar.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16 - 17 cm
L	*/17	*/27	17 - 18 cm
XL	*/18	*/28	18 - 19 cm

* Insert product code



NEW

Air X Low Profile Spica

Product Code: **LPS**


A low profile thumb spica offering targeted support to the thumb MCP and CMC joints.

The unique construction of the closure area offers excellent adjustment and stability around the joint.

The 3D structure of the AirX textile gives a high degree of cushioning and increases comfort and compression.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16 - 17 cm
L	*/17	*/27	17 - 18 cm
XL	*/18	*/28	18 - 19 cm

* Insert product code



Low profile hook



Hook attractive reinforced section for greater adjustment



 Procool

Procool Wrist Thumb Support

Product Code: **PWTS**

An adjustable wrist thumb wrap offering good support around the wrist, MCP and CMC joints. Adjustable thumb closure with excellent shaping around the palm.

Indications Include:

- Basal thumb arthritis
- Thumb CMC & MCP arthritis
- Post UCL injury
- Degenerative joints
- Rheumatoid arthritis
- Osteo-arthritis

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	12 - 14 cm
S	*/13	*/23	14 - 16 cm
M	*/15	*/25	16 - 18 cm
L	*/17	*/27	18 - 20 cm
XL	*/18	*/28	20 - 22 cm

* Insert product code



CARBON SX

Carbon SX Wrist Brace

Product Code: **CSX**

A newly developed ultra thin material that allows the CARBON SX brace to conform exactly to the contours of the wrist and hand.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code

This unique material has high compression features and provides excellent support to the injured wrist and allows full function of the fingers and thumb at the same time as protecting the wrist.



BLACK SEAL

Black Seal SX Wrist Brace

Product Code: **SSX**

A four way stretch high compression Lycra is combined with a fully breathable microfleece to provide a material that offers great wrist support and compression. Lycrafleece is a comfortable and durable material and has unique properties that keep the patient's wrist and hand at the optimal temperature. The inner microfleece also helps wick away moisture from the skin to allow greater patient comfort.

Use CSX size chart above



air

Air SX Wrist Brace

Product Codes: **ASX**

Using a 3D four way stretch fabric the Air brace is designed to offer compression, support and comfort. The Air brace offers excellent breathability and temperature control due to the open 3D construction of the textile. It is well suited to those patient's that prefer the security of a thicker more cushioned material.

The Air brace utilises the 3D material to offer a comfortable, well fitting brace with excellent breathability and overall temperature control.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Functional Resting Splint

Product Codes: **JRS** (Padded Cover), **JRSE** (Easy Liner)

A functional resting splint that is available in two options:

- **Full Padded Cover** - Preformed plastic shell with fully removable padded cover and comfortable wide closure straps. This is excellent for patients using the splint on a longer term basis and the fully padded cover offers excellent comfort and compliance to those patients that are very sensitive.

- **Easy Liner** – Preformed plastic shell with a towel liner that is attached using hook/loop tabs. The shell is supplied with two towel lined padded liners that can be easily removed and replaced. Removable straps can be positioned as required. This option is more suited to conditions that benefit from easily changeable liners such as burns or trauma and is also a lower cost option

Size	Left (Beige)	Right (Beige)	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Push CMC

The Push CMC consists of a cylindrical section around the ball of the thumb and a fixed connection running along the palm to the outside of the hand. This can be closed at the back of the hand with two non-elastic bands. The cylindrical section consists of two components: an aluminium strip and a flexible, synthetic part. The aluminium strip can be shaped and individually adapted for maximum stabilisation of the base of the thumb.

The selected materials and rounded edges guarantee maximum comfort. The simple closing system enables patients to fit the brace themselves without any effort.

The Push CMC is a slim and compact product. The brace can be used in the water. Hygiene is important around the hand. Therefore the synthetic material is antibacterial and the brace can be washed in a washing machine at 40°.



Judy C. Colditz

OT/L, CHT, FAOTA International hand therapy expert, President IFSHT

The Push CMC is an ingenious design that solves the problem of pain for those with CMC osteoarthritis, including me! I am eager to wear it because of its comfort, lack of functional impediment, and durability. In my experience it is the only commercially produced thumb brace that is superior to what I could custom mold for a patient.

Size	Code (Left)	Circumference
1	310111	16 - 19.5 cm
2	310112	19.5 - 22.5 cm
3	310113	22.5 - 26 cm

Size	Code (Right)	Circumference
1	310121	16 - 19.5 cm
2	310122	19.5 - 22.5 cm
3	310123	22.5 - 26 cm



Circumference





Single Armchair

Product Code: **CAP/2**

A dynamic splint with coil design for PIP and DIP joint flexion contractures and Boutonniere's deformity. Maintains assisted extension at the PIP joint while allowing active flexion and can be used for soft tissue injury affecting the PIP joint with a possible flexion deformity.

Size	Code	Circumference at PIP
XS	*XS	4 - 5 cm
S	*S	5.5 - 6 cm
M	*M	6 - 6.5 cm
L	*L	6.5 - 7 cm
XL	*XL	7 - 7.5 cm

* Insert product code



Double Armchair

Product Code: **CAP/3**

A double armchair design for index/middle or middle/ring fingers which maintains assisted extension at the PIP joint while allowing active flexion. A dynamic splint with coil design for PIP and DIP joint flexion contractures and Boutonniere's deformity. Can be used for soft tissue injury affecting the PIP joint with a possible flexion deformity.

Size	Code	Circumference at PIP
S	*S	5.5 - 6 cm
M	*M	6 - 6.5 cm
L	*L	6.5 - 7 cm
XL	*XL	7 - 7.5 cm

* Insert product code



Triple Armchair

Product Code: **CAP/8**

A dynamic splint with coil design for PIP and DIP joint flexion contractures which allows DIP joint to be free for active range of movement. The shorter, stronger spring supports PIP, leaving the DIP joint free to flex. Can be used as a progressive splint after intensive use of the Single Armchair Splint CAP/2.

Size	Code	Circumference at PIP
XS	*XS	5 - 5.5 cm
S	*S	5.5 - 6 cm
M	*M	6 - 6.5 cm
L	*L	6.5 - 7 cm
XL	*XL	7 - 7.5 cm
XXL	*XXL	7.5 - 8 cm

* Insert product code



Cambridge Finger Splint

Product Code: **CAP/5**

A dynamic splint with coil design for PIP joint flexion contractures, which allows DIP joint to be free for active range of movement. Offers the same level of support as the Simpson Armchair Splint CAP/4, providing additional support to the middle phalanx by the application of a padded strap at MP.

Size	Code	Circumference at PIP
S	*/S	6 cm
M	*/M	7 cm
L	*/L	8 cm

* Insert product code



Simpson Modified Armchair

Product Code: **CAP/4**

A dynamic splint with coil design for PIP and DIP joint flexion contractures which allows DIP joint to be free for active range of movement. The shorter, stronger spring supports PIP, leaving the DIP joint free to flex. Can be used as a progressive splint after intensive use of the Single Armchair Splint CAP/2.

Size	Code	Circumference at PIP
XS	*XS	5 - 5.5 cm
S	*S	5.5 - 6 cm
M	*M	6 - 6.5 cm
L	*L	6.5 - 7 cm
XL	*XL	7 - 7.5 cm
XXL	*XXL	7.5 - 8 cm

* Insert product code



Plastic Mallet Splint

Product Codes: **MFS**

Supports the DIP joint in neutral, to help correct mallet deformity, while protecting the distal phalanx and extending the PIP joint. Perforated design prevents skin maceration.

Size	Code	NDC Code	Circumference at PIP
S	*/1	052745	4.5 - 5 cm
S-M	*/2	052752	5 - 5.5 cm
M	*/3	052769	6 - 6.5 cm
M-L	*/4	052776	6.5 - 7 cm
L	*/5	052783	7 - 7.5 cm
L-XL	*/5.5	081608	7.5 - 8 cm
XL	*/6	090402	8 - 8.5 cm

* Insert product code

(Supplied in packs of 10)



Double Bedford Finger Stall

Product Code: **DBF**

Cotton/Elastic fabric conforms to the size of the finger. The stalls support both the PIP and DIP joints.

Contains Latex

Size	Code	Width	Length
Paediatric	*/01	4.5 cm	5 cm
Child	*/02	4.5 cm	7 cm
Small	*/03	6.25 cm	7 cm
Medium	*/05	6.25 cm	9 cm
Large	*/07	6.25 cm	10 cm

* Insert product code

(Supplied in packs of 10)



Whitechapel Finger Splints

Product Code: **WFS**

To increase range of movement by securing the injured finger next to a stronger finger ('buddy strap' effect). Use for soft tissue injury to PIP joint. Made in 12.5cm elasticated loop and plain hook.

Contains Latex

Size	Code	Length
S	*/15	15 cm
M	*/30	30 cm
L	*/50	50 cm

* Insert product code

(Supplied in packs of 50)



Swan Neck



Boutonniere



The economical Oval-8 Kit includes 42 splints in sizes 3-14, a carrying case, a ring gauge, a pen light and instructions. Packaged with most splints in the more commonly used mid-range sizes.

The Oval-8 Sizing Set includes one each of splint size 3-14 and the ring gauge.

Oval - 8 Finger Splint


Complete Kit Code: **P1008-K** Sizing Set Code: **P1008-S**

Made of 1.6mm high-temperature polypropylene, Oval-8 splints can be adjusted with a heat gun for a small angle alteration. With a simple turn of the splint, each Oval-8 Finger Splint can be used to correct Swan Neck, Boutonniere and Mallet deformities, reduce Lateral Deviation and rest Trigger Finger.

Oval-8 splints feature a wide, flat band for comfortable pressure distribution and control. Available in 12 sizes, each splint fits both a full and half ring size. Sizes correspond to jeweller's ring sizes, though correlation may vary due to swelling or joint changes.

Refills (Pack of 5)	Size	Individual	Size
P1008-2-5	2	P1008-2	2
P1008-3-5	3	P1008-3	3
P1008-4-5	4	P1008-4	4
P1008-5-5	5	P1008-5	5
P1008-6-5	6	P1008-6	6
P1008-7-5	7	P1008-7	7
P1008-8-5	8	P1008-8	8
P1008-9-5	9	P1008-9	9
P1008-10-5	10	P1008-10	10
P1008-11-5	11	P1008-11	11
P1008-12-5	12	P1008-12	12
P1008-13-5	13	P1008-13	13
P1008-14-5	14	P1008-14	14

Cervical & Upper Limb





Universal Two Piece Collar

Product Code: **UTP/10**

A universal two piece height adjustable collar that offers good support for the upper cervical spine.

- Preformed polyethylene shells with comfort liners
- Multiple height adjustment
- Large anterior aperture to allow airway management
- Lateral hook/loop closures with slim posterior panel for ease of fitting

Size
Universal



AC Collar

Product Code: **AC**

Made of 3mm polythene plastic with foam-filled, padded vinyl edges and ventilation holes. The two-part design allows variable height adjustment from 7.5cm to 11cm, providing a strong, rigid support for the unstable neck.

Size	Code	Collar Length
S	*/1	40.5 cm
M	*/2	46 cm
L	*/3	51 cm

* Insert product code



Budget Collar

Product Code: **BC**



Made in fully washable, medium density closed cell foam, which will not collapse with use. Stabilises the neck after trauma, or for long term use as a night resting collar.

Size	Code	NDC Code	Depth x Length
S	*/1	019847	7.5 x 43 cm
M	*/2	019854	9 x 47 cm
L	*/3	019861	9 x 51 cm
L/Deep	*/4	027026	10 x 51 cm

* Insert product code

Philly Collar

Product Code: **PCC**



A lightweight, two-piece collar in Plastazote, with rigid occipital and mandibular posts. Strong comfortable support for the head and neck, incorporating a tracheal access hole. A Philly Collar Stabiliser (Product code: PCS/10), one size, is available on request; its anterior and posterior attachments further restrict cervical movement.

Code	NDC Code	Neck Circumference	Collar Depth
*/33	106868	25.5 x 30.5 cm	8.5 cm
*/53	106882	33 x 38 cm	8.5 cm
*/73	-	41 x 48 cm	8.5 cm
*/34	-	25.5 x 30.5 cm	11 cm
*/54	106899	33 x 38 cm	11 cm
*/74	106905	41 x 48 cm	11 cm

* Insert product code

(5.5cm & 13.5cm depths available on request)

Halton Collar

Product Code: **FC**

Available in 8 sizes to fit most patients. Made in fully washable medium density closed cell foam, which will not collapse with use.

Size	Code	Depth x Length
Large	FC/1	9 x 51 cm
Medium	FC/2	9 x 46 cm
Small	FC/3	9 x 42 cm
Long/Medium	FC/4	7.5 x 51 cm
Medium/Medium	FC/5	7.5 x 46 cm
Long/Shallow	FC/6	7 x 48.5 cm
Medium/Shallow	FC/7	7 x 44.5cm
Extra Small	FC/8	7 x 39.5 cm

Universal Lancaster Shoulder Sling

Product Code: **LSS/10**

- Padded, Non-stretch straps
- Fast release buckles
- 'Y'-Tab Adjustments
- Fold lines



Size	Code	NDC Code	Dimensions
One Size	LSS/10	165605	13 x 41 cm



Adjust sling length
using fold lines and circular hook pads



Cervical & Upper Limb



Fast release buckles for easy patient removal & re-application



Easily adjust body and shoulder strap length



Removable Hook/Loop
'Y' - Tabs



Scan QR Code
to watch video

www.youtube.com/promedicsorthopaedic





Deluxe Shoulder Immobiliser

Product Code: *DSI*

A fully-adjustable body binder and humeral strap attachment holds the arm firmly to the body. The wrist sling is made in padded foam for comfort, and is attached to a strong webbing strap to position the wrist and stabilise the arm and shoulder.

Universal design fits left and right.

Size	Code		Dimensions
One Size	DSI/10	Body Binder	15 x 119 cm
		Wrist Cuff	15 x 18 cm
		Humeral Strap	8 x 51 cm



Bradford Arm Elevator Sling

Product Code: *AESB*

Made in hook/loop attractive foam, this sling supports the forearm and upper arm close to the body, using the shoulder strap and body strap.

The strap extending from the wrist section includes a D-ring which can be attached to a drip stand to allow arm elevation.

Universal design fits left and right arm.

Size	Code	NDC Code	Dimensions
One Size	AESB/10	120802	28 x 43 cm



Economy Shoulder Immobiliser

Product Code: *ESI*

An arm sling made in coutil cotton with an adjustable shoulder strap and body strap, to stabilise the arm close to the body.

Universal design fits right and left arm.

Size	Code	Dimensions
S	*/03	14 x 38 cm
M	*/05	17 x 43 cm
L	*/07	19 x 48 cm

* Insert product code





Sling and Swathe


Product Code: **SSW**

A simple padded sling supports the arm, and webbing shoulder straps cross at the back to hold the arm in position. An adjustable body binder holds the upper arm close to the body for increased shoulder stability.

Universal design fits right and left arm.

Size	Code	Body Binder	Padded Sling
S	*/03	13 x 41 cm	20 x 28 cm
M	*/05	13 x 94 cm	23 x 31 cm
L	*/07	15 x 119 cm	25 x 33 cm

* Insert product code



Universal Pro Sling

Product Code: **UPS**

A universal collar and cuff type sling that is fully adjustable to fit all sizes from child to adult. Hook/loop fastening with adjustable loop buckle to secure sling around the wrist.

Universal design fits right and left arm.

Size	Code	NDC Code	Dimensions
Standard	UPS/10	058723	5 x 112 cm

(Supplied in packs of 10)



Scotland
NDC Codes



Boscombe Sling

Product Code: **BOS**

A development of the Pro Sling design, incorporating an adjustable double sling which supports shoulder, elbow, forearm and wrist.

Size	Code	Dimensions
S	*/03	5 x 102 cm
M	*/05	5 x 117 cm
L	*/07	5 x 132 cm

* Insert product code



Elevation strap for positioning.


Trauma Sling

Product Code: **TAU/10**

The Trauma Sling is a high quality universal sling that has multiple uses within the healthcare field. The Trauma Sling can be used as a traditional sling or for arm elevation following trauma, burns or post-operative management of the shoulder or upper arm.

The design and quality of materials used offer protection, comfort and ease of access to the damaged limb.

The Trauma Sling is easy to apply and is supplied vacuum packed to reduce storage volume.



Adjustment at elbow allows for varying degree of elbow flexion/extension.

Size
Universal



Paediatric Trauma Sling

Product Code: **TAU/P**

Same excellent features as Trauma Sling, but with closed seam behind elbow.

The design and quality of materials used offer protection, comfort and ease of access to the damaged limb. The Trauma Sling is easy to apply and is supplied vacuum packed to reduce storage volume.



Size
Universal

- Comfortable
- Protective
- Universal



Jura Abduction Sling PLUS

Product Code: **JASSP/10**

The Jura Plus has a three point offloading harness that offers excellent support and comfort. The contoured abduction pad is profiled for improved patient compliance and holds the arm in 45° of abduction. A body swathe combined with upper/lower arm wraps offer good stability and shoulder fixation. Hand resting splint supplied.

- Rotator cuff repairs
- Immobilisation following joint replacement
- Glenohumeral dislocations / subluxation
- Posterior/ Anterior capsule repairs
- Total shoulder reconstructions
- Soft tissue repairs / strains

Size
Universal



 Scotland
NDC Code

Jura Shoulder Immobiliser

Product Code: **JSI**

The Jura Shoulder Immobiliser can be used as an arm sling or as a shoulder immobiliser and it is designed to be used on both left and right arms.

Strong cotton sling with fully adjustable padded shoulder and body straps.

Size	NDC Code	Forearm Length
XS	-	< 33 cm
S	-	33 - 38 cm
M	-	37 - 43 cm
L	205691	42 - 48 cm
XL	-	48 cm >

* Insert product code



Jura Shoulder Lock

Product Code: **JSL/10**

A universal fast-fit low profile shoulder immobiliser. The Jura Shoulder Lock comes complete with a body wrap to which the upper arm wrap is secured to effectively immobilise the shoulder. The wrist wrap fixes to the anterior of the body wrap and this in turn is secured to the shoulder strap using a fast-fit snap buckle.

The shoulder strap can be adjusted posteriorly to achieve a comfortable fit for the patient and also to alter the brace to fit left or right. There is a large comfort pad to distribute pressure over the shoulder /clavicle.

This is a shoulder immobiliser that is easy to fit and the modular design greatly improves patient compliance.

Size
Universal



SlingShot® 3 Shoulder Brace

The SlingShot® 3 is a clinician-driven innovation in shoulder bracing. Its design provides superior comfort for the patient recovering from surgery while providing multiple options of post-operative support.

- Innovative Offloading Shoulder Harness for superior patient comfort
- Convertible Abduction Pillow for Post-Operative options (15° Abduction, 90° Neutral and 45° Neutral)
- Comprised entirely of Airmesh® for enhanced breathability
- Lined with a moisture wicking blue fabric for maximum comfort
- Quick release shoulder and waist straps
- Four sizes (S-XL)
- Universal left or right includes exercise ball

Code	Description
0004X	SlingShot 3, S-XL

X = Insert size # from table below

Size #	Size	Length (Olecranon to Knuckles)
X = 2	S	29 - 33 cm
X = 3	M	34 - 37 cm
X = 4	L	38 - 41 cm
X = 5	XL	42 - 44 cm



Atlas Universal Shoulder Brace

The Atlas Universal Sling is designed to provide superior comfort for the patient recovering from shoulder surgery. It features the offloading shoulder harness (patents pending) to promote comfort and all day wear. It also features a comfortable, breathable Airmesh sling that is universal in size, to fit most patients. The product contains a 15° abduction pillow and quick release shoulder and waist strap buckles for easy application.

- Innovative Offloading Shoulder Harness for superior patient comfort
- Universally sized shoulder sling
- 15° Abduction pillow (with detachable waist strap)
- Comprised entirely of Airmesh®, a moisture wicking fabric for enhanced breathability
- Quick release shoulder and waist straps
- Universal left or right

Code
00052



Universal Combo Sling

Product Code: **UCS**

A comfortable offloading sling designed for patients recovering from shoulder surgery. Two wedges are included which can be applied in a choice of 3 different combinations: 45° Neutral, 90° Neutral or 15° Abduction. Lightweight and breathable 3D spacer fabric sling.

- Exercise ball
- Detachable waist strap
- Opening to relieve elbow pressure



Fold to fit length

Size
Universal



Evans Abduction Sling

Product Code: **EAS**

This is a low cost universal option that utilizes the popular Lancaster Shoulder Immobiliser. In addition to the comfortable sling with shoulder and body straps – the Evans Abduction Sling incorporates a 15° abduction pad that is easily secured to the inside of the sling.

The inner profile of the abduction pad is shaped to rest comfortably against the patient. The Evans Abduction sling is useful in the treatment of various shoulder conditions

- Rotator cuff repairs
- Glenohumeral dislocations / subluxation
- Posterior/ Anterior capsule repairs
- Soft tissue repairs / strains

Size
Universal



Jura Abduction Sling

Product Code: **JASS**

The Jura Abduction has a three-point offloading harness that offers excellent support and comfort. The sling is made from a lightweight breathable 3D spacer fabric and attaches to the abduction pads. The sling is supplied with a contoured 15° pad that can be used on its own.

In addition to this there is an additional pad that increases abduction to 45° and positions the shoulder in a neutral position. A detachable exercise ball is also provided. The shoulder harness is fitted with easy release buckles to facilitate easy application and removal of the sling.

Size	Code	Forearm Length
S	*/03	38 cm
M	*/05	40 cm
L	*/07	43 cm
XL	*/08	46 cm

- Rotator cuff repairs
- Glenohumeral dislocations/ subluxation
- Posterior/ Anterior capsule repairs
- Total shoulder reconstructions
- Soft tissue repairs / strains

* Insert product code



Jura Clavicle Splint

Product Code: **JCS/10**

Soft padded shoulder straps with return loops and padded back pad. Comfort pads can be positioned anywhere along the length of the straps – for example in the axilla area.

Straps can be easily trimmed to length and refastened using the hook Y-tab closures.

Size
Universal




Subluxation Cuff

Product Code: **SCF**

A humeral cuff with soft cotton lining supports the injured arm, while the padded webbing strap fits under the opposite axilla. Universal design fits right and left arm.

Size	Code	Mid Upper Arm Circumference
S	*/03	24 - 29 cm
M	*/05	29 - 34 cm
L	*/07	34 - 40 cm

* Insert product code



Clavicle Splint

Product Code: **CLS**

Two padded webbing straps adjust in a figure of eight design fixed to a pad over the scapulas. Designed to assist shoulder retraction.

The length of the shoulder strap is measured from top of back pad, under axilla, to bottom of back pad.

Size	Code	Length of Shoulder Straps
XS	*/02	48 cm
S	*/03	56 cm
M	*/05	68 cm
L	*/07	87 cm
XL	*/08	92 cm

* Insert product code



Neoprene Shoulder Brace

Product Code: **NSB**

Two adjustable upper arm cuffs fix with hook/loop at the back of the neck, to prevent shoulder subluxation and hold the head of humerus in the glenoid cavity.

Size	Code	Mid Upper Arm Circumference
S	*/03	23 - 28 cm
M	*/05	28 - 33 cm
L	*/07	33 - 38 cm
XL	*/08	38 - 43 cm

* Insert product code



BREG™

T Scope® Elbow Premier

Product Code: **0725x**

The T Scope Elbow Premier is designed for fixed or controlled range of motion in the treatment of ligamentous injuries and stable fractures of the elbow or upper arm. This brace incorporates the patented T Scope Premier hinge and telescoping uprights.

Designed for exceptional patient fit and enhanced compliance, the T Scope Elbow Premier is one of the most lightweight, low profile post-op elbow braces on the market and features a quick lock to prevent flexion / extension at any angle. Malleable aluminum bars can be manipulated to accommodate patient swelling and post-operative bandages.

Code	Description
07254	T Scope Elbow Premier, Left Reg.
07255	T Scope Elbow Premier, Right Reg.
07256	T Scope Elbow Premier, Left Long
07257	T Scope Elbow Premier, Right Long



Cubital Tunnel Splint

Product Code: **CTS/10**

A universal easy to fit splint that effectively restricts elbow flexion to help reduce the pain associated with Cubital Tunnel Syndrome. A soft padded wrap with an adjustable aluminium stay. The splint is easy to apply, lightweight and breathable.

Size
Universal



Hinged Elbow Wrap

Product Code: **HEW**

A front-opening Neoprene sleeve with hinges, which limit medial/lateral movement and elbow flexion. Hook/loop straps above and below the elbow ensure a good fit. Specifically designed for the patient with limited dexterity.

Size	Code	Circumference 13cm Below Elbow Crease	Circumference 13cm Above Elbow Crease
S	*/03	13 - 18 cm	23 - 28 cm
M	*/05	18 - 23 cm	28 - 33 cm
L	*/07	23 - 28 cm	33 - 38 cm
XL	*/08	28 - 33 cm	38 - 43 cm

* Insert product code

Orthoprene Elbow

Product Code: **OES**

A simple pull-on Neoprene sleeve to give light support and compression around the elbow area and ease the pain of inflammation.



Size	Code	Elbow Circumference
S	*/03	18 - 20 cm
M	*/05	23 - 25 cm
L	*/07	28 - 30 cm

* Insert product code

Tennis Elbow Strap

Product Code: **TEP**

A D-ring on strong cotton webbing adjusts tension to ease the pain of medial or lateral epicondylitis.



Size	Code	Strap Length
One Size	TEP/10	5 x 46 cm

(Supplied in packs of 6)



Jura Clasp

Product Code: **JUC**

Designed to reduce the pain and discomfort associated with elbow tendonopathies, conditions more commonly referred to as epicondylitis and tennis or golfer's elbow.

Contains Latex

Size	Code	Arm Circumference
XS	*/02	22 - 23 cm
S	*/03	24 - 25 cm
M	*/05	26 - 27 cm
L	*/07	28 - 30 cm
XL	*/08	31 - 34 cm

* Insert product code



 Scotland
NDC Codes

Vari Pad Tennis Elbow Support

Product Code: **VPS**

Made in Neoprene with a D-ring hook/loop strap and adjustable pressure pad, which is easily positioned on the insertion point to ease the pain of medial or lateral epicondylitis.

This support has a high level of patient compliance, as the stretch Neoprene accommodates the changing muscle bulk of the moving forearm.

Size	Code	NDC Code	Circumference 2.5cm Below Elbow Crease
XS	*/02	-	18 - 20 cm
S	*/03	052547	23 - 25 cm
M	*/05	052554	28 - 30 cm
L	*/07	052561	33 - 35 cm
XL	*/08	052738	37 - 40 cm



Replay Tennis Elbow Support


Product Code: **RVP**

A tennis elbow support in 'Ventoprene', a ventilated Neoprene providing compression and heat retention. The adjustable pressure pad can be positioned to ease the pain of medial or lateral epicondylitis.

Stretch 'Ventoprene' accommodates the changing muscle bulk of the moving forearm when playing racket sports.

Size	Code	Circumference 2.5cm Below Elbow Crease
S	*/03	23 - 25 cm
M	*/05	28 - 30 cm
L	*/07	33 - 35 cm
XL	*/08	37 - 40 cm

* Insert product code



Lower Limb



T Scope® Premier Post-op Knee Brace

The Breg T Scope Premier Post-Op knee brace is designed to provide protected, controlled range-of-motion for patients recovering from knee surgery or those who have knee injuries or instabilities. The brace is comfortable, lightweight, simple to apply, and easily adjustable for a personal fit.

It features a patented, user-friendly hinge to control the range of motion of the patient's knee. It is universal, and fits most patients so there are less SKUs to stock in inventory. The innovative BridgeTech Incision Pad (an optional accessory) provides added comfort around the incision site.

Code	Description	NDC Code
08814	T Scope Premier 17"-27" (43.2 cm - 68.6 cm) fits up to a 30.5" (77.5cm) thigh, universal Right or Left	204588
08815	T Scope Premier, XL 17"-27" (43.2 cm - 68.6 cm) fits up to a 35.5" (90.2cm) thigh, universal Right or Left	-
08816	T Scope Premier, Full Foam, XL 17"-27" (43.2 cm - 68.6 cm) fits up to a 35.5" (90.2cm) thigh, universal Right or Left	-
73621	BridgeTech Incision Pad (accessory), universal	-

Common Examples of Use

- ACL, PCL, MCL, LCL repairs/injuries
- Tibial plateau fractures
- Osteochondral repairs
- Meniscal repairs
- Patella tendon repairs
- Condylar fractures
- Sprains/strains of the knee
- HTOs (High Tibial Osteotomy)



BridgeTech Incision Pad

Performance

- Easy to use ROM hinge with quick-adjusting flexion and extension stops to control and limit knee flexion and extension
- Telescoping calf and thigh sleeves for sizing a wide range of patients from 5' to 6'4" tall
- Easy to use extension drop lock allows the brace to be locked out in 5 positions (in -10°, 0°, 10°, 20°, and 30°) of knee extension with the push of a button
- Stiff telescoping bars and overall brace design limits knee range of motion
- Accommodates both right and left leg procedures, reducing the need to stock additional products

Easy to Use

- Extremely quick application time and simple to adjust to enhance patient compliance and comfort
- Easily and quickly customized for fit across a wide range of patient sizes
- Trimmable straps for size customization
- Strap lock clips keep straps positioned on bars and the hinges in the correct placement around the knee
- Quick clip buckles allow for easy application and removal of the brace
- QR Code video helps patients to fully understand the use and care of the brace -



Scan QR Code
to watch video

Comfort

- Comfortable padding throughout the brace, designed with patient compliance in mind
- Brace is lightweight, weighing approx 990g (standard version)
- BridgeTech Incision pad relieves pressure and provides added comfort around the incision site (optional accessory)
- Excellent overall suspension due to the design of the pads, strapping and overall brace
- More flexible cuffs allow for better comfort and suspension of the brace
- Near universal sizing, with the hinge bars extending from 17" to 27" (43.2 cm - 68.6 cm) to accommodate a wide range of patient anatomies
- Standard version has reduced foam vs. other post-op braces for greater breathability

The closer you look, the more you'll like the NEW T-Scope Premier...



Axiom-D Elite Ligament Knee Brace



Dynamic Bracing, a new vision for conservative treatment

- For ACL, PCL or combined instability injuries
- Aluminium / Magnesium lightweight frame
- Provides full-bodied support and protection
- Multiple configurations for customisation
- Low profile brace contours to the leg
- Over centre buckles & pivoting straps



PCL
Counterforce
Strap



Dynamic Technology
Enhanced control of tibial translation that provides stability at the moment of clinical need for ACL, PCL or combined ligament injuries

The Axiom-D Elite delivers innovation for ligament bracing. The unique dynamic hinge delivers a progressive counterforce to the lower leg to reduce tibial translation in patients with ligament deficiencies or recovering from ligament reconstruction. This maintains the proper anatomic relationship between the femur and tibia as it goes through a full range of motion.

The Axiom-D Elite is rigid, durable and comfortable providing full-bodied support and protection. It is available in a dual upright 43 cm frame for robust and extended support. The redesigned low profile brace contours to the leg and is offered in aluminum or magnesium, with multiple configurations for customization. Standard Fit, Athletic Cut or Custom Manufacture, the Axiom-D Elite is a knee brace as demanding as the patient requires it to be. Whether you're competing in contact sports, or playing as a weekend warrior your patients are supported for a full range of activities.

Code	Description
PK428YXX	Axiom-D Elite, Aluminum
PK528YXX	Axiom-D Elite, Combined Instability, Aluminum
PK628YXX	Axiom-D Elite, Sport
PK423YXX	Axiom-D Elite, Athletic
PK523YXX	Axiom-D Elite, Athletic, Combined Instability
PK623YXX	Axiom-D Elite, Sport, Athletic
MG428YXX	Axiom-D Elite, Magnesium
MG528YXX	Axiom-D Elite, Combined Instability, Magnesium
MG423YXX	Axiom-D Elite, Athletic, Magnesium
MG523YXX	Axiom-D Elite, Athletic, Combined Instability
MG428000	Axiom-D Elite Custom, Magnesium
PK428000	Axiom-D Elite Custom, Aluminum
PK628000	Axiom-D Elite Custom, Sport, Aluminum

X = Insert size # from table below



Code	Size	Circumference Thigh	Circumference Calf
XX = 01	XS	34 - 41 cm	32 - 36 cm
XX = 03	S	41 - 48 cm	36 - 39 cm
XX = 05	M	48 - 55 cm	39 - 43 cm
XX = 07	L	55 - 62 cm	43 - 47 cm
XX = 09	XL	62 - 69 cm	47 - 51 cm
XX = 11	XXL	69 - 75 cm	51 - 53 cm
XX = 13	3XL	75 - 89 cm	55 - 58 cm

Thigh circumference measure 15 cm above mid-patella.
Calf circumference measure 15 cm below mid-patella.



NEW



Freerunner Knee Brace

The FreeRunner knee brace is the premier soft good for patellofemoral stabilization. It is designed for patients who have patellofemoral discomfort resulting from maltracking, subluxation or other patellofemoral injuries. It combines a patent pending offset CAM hinge for dynamic support and introduces innovative Smart-Zone™ compression fabrics, for a low profile and ultra lightweight silhouette. The difference in patellofemoral bracing is dynamic; now athletes and patients can focus on the game and not the brace.

- Patella dislocation/subluxation
- Chronic patellar maltracking
- Patella tendonitis
- Mild chondromalacia
- Runners Knee
- Patellofemoral pain syndrome
- Post lateral release

Code	Description
100379-1XX	FreeRunner, Left
100379-2XX	FreeRunner, Right

X = Insert size # from table below

Code	Size	Circumference Thigh
XX = 05	XXS	29 - 34 cm
XX = 10	XS	34 - 39 cm
XX = 20	S	39 - 44 cm
XX = 30	M	44 - 49 cm
XX = 40	L	49 - 55 cm
XX = 50	XL	55 - 62 cm
XX = 60	XXL	62 - 70 cm

Circumference taken 10 cm above knee center



Freestyle OA Knee Brace

Promedics is excited to announce the latest addition to the osteoarthritis product offering. The Freestyle OA provides offloading of the medial or lateral compartment for mild to moderate osteoarthritis. Specifically designed with the patient's needs in mind.

The simple wrap around design allows patients to apply and remove the brace with ease. Comfortable, breathable material provides focused compression around the knee joint to address swelling and promote blood flow. The single-sided low profile hinge encourages patient compliance, while the innovative thumbwheel dial allows precision unloading adjustments without the use of a tool.

Code	Description
1171X	Freestyle OA Medial, Left
1172X	Freestyle OA Medial, Right
1176X	Freestyle OA Lateral, Left
1177X	Freestyle OA Lateral, Right

X = Insert size # from table below

Code	Size	Circumference 15cm Above Mid Patella	Circumference 15cm Below Mid Patella
X = 2	S	38 - 46 cm	30 - 36 cm
X = 3	M	46 - 53 cm	36 - 38 cm
X = 4	L	53 - 61 cm	41 - 46 cm
X = 5	XL	61 - 69 cm	46 - 51 cm
X = 6	XXL	69 - 76 cm	51 - 57 cm

Fusion® OA Plus

The Fusion® OA Plus knee brace provides patients with medial compartment off-loading for unicompartmental osteoarthritis. It features Breg's exclusive adjustable hinge technology, which provides an effective valgus load to the knee. The hinge utilizes Breg's thumbwheel dial design that allows patients to easily adjust the level of offloading without using a tool. Breg's ProForm technology provides an intimate, contoured fit, that keeps the brace in place all day long. Fusion® OA Plus provides medial / lateral, anterior / posterior stabilization and support.



Code	Description
01221	Fusion® OA Plus Custom
01222	Fusion® OA Plus Custom Color / Pattern
130XX	Fusion® OA Plus, Prefabricated, Left
131XX	Fusion® OA Plus, Prefabricated, Right
787XX	Fusion® OA Plus, Refurbish, Left
788XX	Fusion® OA Plus, Refurbish, Right

X = Insert size # from table below

Code	Size	Circumference 15cm Above Mid Patella	Circumference 15cm Below Mid Patella
XX = 20	S	38 - 46 cm	30 - 36 cm
XX = 30	M	46 - 50 cm	36 - 38 cm
XX = 35	M+	50 - 53 cm	38 - 41 cm
XX = 40	L	53 - 61 cm	41 - 46 cm
XX = 50	XL	61 - 69 cm	46 - 51 cm
XX = 60	XXL	69 - 76 cm	51 - 58 cm

LPR Knee Brace

Breg makes another bold statement as a leader in innovative excellence – with the LPR ligament knee brace. Crafted of high grade, aerospace aluminum and shaped to offer the stiffness and support your patients need for ligament stabilization, the LPR is one of the lightest knee braces available.

Game changing, diagonal strapping opens up the popliteal space and offers full range of motion for all levels of activity. Combining ligamentous support and patient comfort and compliance – the LPR – Ligament Protection Redefined. Off-the-shelf frame color is dark gray, custom color shown in red.



Code	Description
022XX	LPR OTS RT
056XX	LPR OTS LT
00028	LPR Custom
00029	LPR Colour Custom

X = Insert size # from table below

Code	Size	Circumference Knee Joint	Circumference Thigh	Circumference Calf
XX = 10	XS	28 - 32 cm	33 - 39 cm	25 - 30 cm
XX = 20	S	32 - 37 cm	39 - 46 cm	30 - 36 cm
XX = 30	M	37 - 39 cm	46 - 50 cm	36 - 38 cm
XX = 35	M+	39 - 41 cm	50 - 53 cm	38 - 41 cm
XX = 40	L	41 - 46 cm	53 - 60 cm	41 - 46 cm
XX = 50	XL	46 - 50 cm	60 - 70 cm	46 - 51 cm
XX = 60	XXL	50 - 56 cm	70 - 79 cm	51 - 58 cm



Dynamic Unloading Osteoarthritis (DUO) Knee Brace

A unique design, our DUO (Dynamic Unloading Osteoarthritis) brace is the only dual-upright dynamic OA brace on the market. This brace creates a load across the joint when the knee is straight and reduces the load as the knee flexes. Thus, when the patient is bearing weight, the brace is actively working. DUO reduces force when sitting, making it ideal for patients who need to wear a brace for daily activities.

Common Examples of Use

- Unicompartamental OA
- Meniscal instabilities with ligamentous injuries

Product Features

- Low-profile frame
- Dynamic loading hinge
- Incremental arm adjustments
- Up to 16 degrees of correction can be applied
- Unique strapping design minimizes migration



Code	Description
ED312000	DUO Custom Short, Medial
ED314000	DUO Custom Short, Lateral
ED112000	DUO Custom, Medial
ED114000	DUO Custom, Lateral
ED114YXX-B	DUO, Lateral, XS-XXL
ED112YXX-B	DUO, Medial, XS-XXL
ED113YXX-B	DUO, Lateral Athletic, XS-XXL
ED111YXX-B	DUO, Medial Athletic, XS-XXL
ED314YXX-B	DUO, Lateral, Short, XS-XXL
ED312YXX-B	DUO, Medial, Short, XS-XXL
ED313YXX-B	DUO, Lateral Athletic, Short, XS-XXL
ED311YXX-B	DUO, Medial Athletic, Short, XS-XXL

X = Insert size # from table opposite

Code	Size	Circumference Thigh	Circumference Calf
XX = 01	XS	34 - 41 cm	32 - 36 cm
XX = 03	S	41 - 48 cm	36 - 39 cm
XX = 05	M	48 - 55 cm	39 - 43 cm
XX = 07	L	55 - 62 cm	43 - 47 cm
XX = 09	XL	62 - 69 cm	47 - 51 cm
XX = 11	XXL	69 - 75 cm	51 - 53 cm

Thigh circumference measure 15 cm above mid-patella.
Calf circumference measure 15 cm below mid-patella.



Thruster RLF

This single-upright OA brace provides the greatest amount of OA correction on the market. The unique design of the hinge allows for 36° of brace angulation. This allows the brace to unload knee joints that have excessive tissue. It can be made to unload the medial or lateral compartment of the knee.

Code	Description
AG061012	Custom
AG061016	Custom, Short
TA008YXX-B	Lateral, XS-3XL
TA009YXX-B	Medial, XS-3XL
TA018YXX-B	Lateral Short, XS-3XL
TA019YXX-B	Medial Short, XS-3XL
TA028YXX-B	Lateral Athletic, XS-3XL
TA029YXX-B	Medial Athletic, XS-3XL

X = Insert size #
from table below

Code	Size	Circumference Thigh	Circumference Calf
XX = 01	XS	34 - 41 cm	32 - 36 cm
XX = 03	S	41 - 48 cm	36 - 39 cm
XX = 05	M	48 - 55 cm	39 - 43 cm
XX = 07	L	55 - 62 cm	43 - 47 cm
XX = 09	XL	62 - 69 cm	47 - 51 cm
XX = 11	XXL	69 - 75 cm	51 - 53 cm
XX = 13	3XL	75 - 89 cm	55 - 58 cm

Thigh circumference measure 15 cm above mid-patella.
Calf circumference measure 15 cm below mid-patella.



Fusion Knee Brace

Ideal for normal daily activities as well as athletic endeavors, this brace is designed using ProForm technology, which creates a comfortable, precise fit that won't compromise mobility. Available in custom and off the shelf (OTS) versions. Off-the-shelf frame color is black, custom color shown in blue.

- AirTech™ Frame pads – innovative windows and channels in the pad permit cooling air flow while allowing moisture and heat to move away from the skin
- Pivot Point Strap tabs – pivoting strap tabs allow straps to contour to the leg for superior fit, comfort and suspension
- ProForm Technology – medial side structure, now redesigned for superior contoured fit and protection

Code	Description
01209	Fusion Custom Color/Pattern
01207	Fusion Custom
005XX	Fusion OTS, Left
007XX	Fusion OTS, Right
976XX	Refurbish Kit (Thigh), Left*
979XX	Refurbish Kit (Thigh), Right*
984XX	Refurbish Kit (Calf), Left*
985XX	Refurbish Kit (Calf), Right*
00522	Condyle pad (set of 2)

Common Examples of Use

- ACL
- PCL
- Collateral Ligament protection

X = Insert size # from table below

* Kit contains frame pads, straps and instructions

Size #	Size	Circumference 15cm Above Mid Patella	Circumference 15cm Below Mid Patella
XX = 10	XS	33 - 38 cm	25 - 30 cm
XX = 20	S	38 - 46 cm	30 - 36 cm
XX = 30	M	46 - 50 cm	36 - 38 cm
XX = 35	M+	50 - 53 cm	38 - 41 cm
XX = 40	L	53 - 61 cm	41 - 46 cm
XX = 50	XL	61 - 69 cm	46 - 51 cm
XX = 60	XXL	69 - 76 cm	51 - 58 cm





GENURANGE



Lower Limb



Scotland
NDC Codes

Genurange Universal

Product Code: **GN/0510**

To reduce stock-holding of different circumferences, the one size Universal wraps can be cut, without fraying, to fit most sizes. Open-cell padding tibial and femoral liners are attached by four or six straps for an excellent fit.

Hinges, applied to the liners, align with the tibial condyles and the liners can be removed for replacement or laundering.

Available in four lengths. Replacement wraps and straps available, see below.

Size	Code	NDC Code	Length
S	*/18	080106	45 cm
M	*/20	080205	50 cm
L	*/22	-	55 cm
XL	*/24	-	60 cm

* Insert product code



Genuranger

Product Code: **GN/04**

Towelling-lined tibial and femoral liners are attached by four straps; six straps on 24", for an excellent fit. The hinges are applied to the liners to align with the tibial condyles, and the liners can be removed for replacement or laundering.

Available in five circumferences and five lengths. Replacement wraps and straps available.

Replacement Liners and Straps

To order replacement liners and straps for different sizes of Genurange Universal and Genuranger, add 'W' after 'GN' in the code, e.g. GNW/0510/18 or GNW/0405/18

Size	40 cm	45 cm	50 cm	55cm	60 cm	Knee Circumference
S	*03/16	*03/18	*03/20	*03/22	*03/24	35 - 40 cm
M	*05/16	*05/18	*05/20	*05/22	*05/24	40 - 45 cm
L	*07/16	*07/18	*07/20	*07/22	*07/24	45 - 50 cm
XL	*08/16	*08/18	*08/20	*08/22	*08/24	55 - 60 cm
XXL	*09/16	*09/18	*09/20	*09/22	*09/24	60 - 65 cm

* Insert product code



Genurehab

Product Code: **GN/0310**

The 40cm Genurehab is a short version of the Genurange Universal, designed for long term support and motion control during the more active phase of recovery. The one-size Genurehab wraps can be cut without fraying, to fit most sizes.

To order replacement liners and straps, add ' **W** ' after ' **GN** ' in the code, e.g. **GNW/0310**

Size
Universal




Genurange Lite

Product Code: **GN/08**

A 'lite' wrap-around knee brace with ROM hinges. 3D spacer fabric with open mesh posterior panel makes this a supportive, comfortable well fitting brace. Suitable for the treatment various knee conditions including mild/moderate collateral ligament injuries and cruciate mild instability of the cruciate ligaments.

Size	Code	Circumference 15cm above mid-patella
XS	*02	32 - 38 cm
S	*03	38 - 44 cm
M	*05	44 - 51 cm
L	*07	51 - 58 cm
XL	*08	58 - 66 cm
XXL	*09	66 - 75 cm

* Insert product code



GENURANGE

- A geared, polycentric hinge which mimics the action of the knee.
- Four small pins which are inserted into apertures to set the hinge at the required angles of flexion and extension, between 10° and 110°.
- A plastic hinge cover, with screw closure, to prevent patients tampering with the hinge mechanism.
- Features plastic tags to prevent loss of pins.



Jura Knee Immobiliser

Product Code: **BTPK**

A universal three panel knee immobiliser manufactured from a quality foam laminate with terry towel lining. This brace is designed to offer a lower cost alternative to hospital departments with high volume usage.

Code	Splint Length
BTPK/12	30 cm
BTPK/14	36 cm
BTPK/16	41 cm
BTPK/18	46 cm
BTPK/20	51 cm
BTPK/22	56 cm
BTPK/24	60 cm
BTPK/26	66 cm
BTPK/28	71 cm



Post-Op Knee Splint

Product Code: **PKS/20**

A three panel splint made in foam laminate with two detachable side panels which can be positioned to fit most leg circumferences. The four steels are removable for the patient to launder the splint.

The length of the splint is 51cm.

Size
Universal



Three Panel Knee Immobiliser

Product Code: **TPK**

A three panel splint made in high quality foam laminate with bound edges, designed for longer term use. The two detachable side panels can be positioned to fit most leg circumferences.

The four steels are removable for the patient to launder the splint.

Code	NDC Code	Splint Length	Code	NDC Code	Splint Length
*/12	-	30 cm	*/22	-	56 cm
*/14	027323	36 cm	*/24	027347	61 cm
*/16	-	41 cm	*/26	-	66 cm
*/18	-	46 cm	*/28	-	71 cm
*/20	027330	51 cm			

* Insert product code





Jura One Piece Immobiliser

Product Code: **OPK**

A well fitting comfortable knee immobiliser designed to maintain knee extension following surgery, trauma and soft tissue injuries.

- One piece foam laminate wrap with towel lining
- Detachable hook / loop straps for easy positioning

Size	Code	Length	Circumference (Mid Thigh)
Paed	*/12	30 cm	Up to 32 cm
S	*/16	41 cm	Up to 41 cm
S	*/20SM	51 cm	Up to 52 cm
M-L	*/20ML	51 cm	47 - 66 cm
S	*/24SM	56 cm	Up to 52 cm
M-L	*/24ML	56 cm	47 - 66 cm
S	*/28SM	71 cm	Up to 52 cm
M-L	*/28ML	71 cm	47 - 66 cm

* Insert product code



Canvas Knee Immobiliser

Product Code: **CNI**

A one-piece knee immobiliser made in soft padded canvas with four removable steels which can be removed for the patient to launder the splint.

The length of the splint is 48cm.

Size	Code	NDC Code	Knee Circumference
S	*/03	094882	28 - 34 cm
M	*/05	094769	34 - 40 cm
L	*/07	094738	40 - 46 cm
XL	*/08	094806	46 - 56 cm

* Insert product code



Neosplint

Product Code: **NS**

A one-piece knee immobiliser in Neoprene, which conforms very well to the leg. Designed for long term wear, the four steels can be removed to launder the splint.

The length of the splint is 48cm.

Size	Code	Knee Circumference
S	*/03	28 - 34 cm
M	*/05	34 - 40 cm
L	*/07	40 - 46 cm
XL	*/08	46 - 56 cm

* Insert product code



Air X Polycentric Knee Brace

Product Code: AMP

A 30cm long wrap-around knee brace offering good medial and lateral support using quality polycentric hyperstop hinges. Hook/loop return closures make this brace easy to apply and the posterior straps can be adjusted quickly to allow accurate midline positioning of the polycentric joints.

AirX gives great compression and support and the properties of the textile offer a high degree of breathability and reduction in moisture. This makes the brace comfortable to wear for extended periods of time and improves patient compliance.

Size	Code	Knee Circumference
S	*/03	28 - 34 cm
M	*/05	34 - 40 cm
L	*/07	40 - 46 cm
XL	*/08	46 - 56 cm
XXL	*/09	56 - 61 cm

* Insert product code



Air X SK Knee Brace

Product Code: ASK

The SK Knee Brace has been the industry leading wrap-around knee brace for many years and thousands of patients have benefited from this easy to use comfortable brace in the treatment of mediolateral instability.

The SK Knee Brace is now available using the excellent properties of the AirX 3D spacer fabric. AirX gives great compression and support and the properties of the textile offer a high degree of breathability and reduction in moisture. This makes the brace even more comfortable to wear for extended periods of time and improves patient compliance.

Size	Code	Knee Circumference
XS	*/02	25 - 28 cm
S	*/03	28 - 34 cm
M	*/05	34 - 40 cm
L	*/07	40 - 46 cm
XL	*/08	46 - 56 cm
XXL	*/09	56 - 61 cm

* Insert product code



Air X Long ROM Knee Brace

Product Code: ARK

A 50cm long wrap-around knee brace with Range of Motion (ROM) Joints. AirX 3D spacer fabric offers excellent comfort, compression support and fit. The brace has a range of rigid and elastic straps that can be adjusted anteriorly and posteriorly to allow a perfect fit and joint alignment. The ROM joints have a wide range of flexion and extension settings and the brace is suitable for the treatment of various knee conditions including mild/moderate collateral ligament injuries and mild instability of cruciate ligaments.

Size	Code	Knee Circumference
XS	*/02	25 - 28 cm
S	*/03	28 - 34 cm
M	*/05	34 - 40 cm
L	*/07	40 - 46 cm
XL	*/08	46 - 56 cm
XXL	*/09	56 - 61 cm

* Insert product code



Air X VK Long Knee Brace

Product Code: **AVK**


The VK Knee Brace has been one of the industry leading long wrap-around knee braces for many years and thousands of patients have benefited from this easy to use comfortable brace in the treatment of mediolateral instability. The longer length (41cm) provides increased medial and lateral support for the unstable valgus or varus knee. The front opening design makes it easy to apply – even for those patients with limited dexterity

The VK Knee Brace is now available using the excellent properties of the AirX 3D spacer fabric.

AirX gives great compression and support and the properties of the textile offer a high degree of breathability and reduction in moisture. This makes the brace even more comfortable to wear for extended periods of time and improves patient compliance.

Size	Code (Left)	Code (Right)	Knee Circumference
XS	*/12	*/22	25 - 28 cm
S	*/13	*/23	28 - 34 cm
M	*/15	*/25	34 - 40 cm
L	*/17	*/27	40 - 46 cm
XL	*/18	*/28	46 - 56 cm
XXL	*/19	*/29	56 - 61 cm

* Insert product code



- **Breathable**
- **Comfortable**
- **Skin Friendly**

Moisture-transporting material

The unique AirX™ fabric has the ability to transport moisture efficiently away from the body. The layering of fibre structures promotes a moisture wicking process, which results in a dry comfortable feeling to the surface of the skin so that body always feels fresh and dry. AirX™ gives the body optimal freedom to breathe – reduces the risk of skin maceration which results in improved patient compliance.

A pleasant surprise

It is hard not to be impressed. AirX™ is made of a soft, lightweight material that perfectly follows the contours of the body. This means that AirX™ adapts itself perfectly to the body so that it feels extremely comfortable to wear – even for long periods at a time – a feature, which is often a problem with orthopaedic and sport medical braces and supports. In this way, AirX™ is a comfortable material that helps to to increased patient compliance without you even knowing!



Ventilated SK Knee Plus

Product Code: **VSKP**

The SK Knee Brace has been the industry leading wrap-around knee brace for many years and thousands of patients have benefited from this easy to use comfortable brace in the treatment of mediolateral instability.

The Ventilated SK Plus takes the brace to the next level- using multi-perforated towel-lined latex free neoprene - the SK Plus provides excellent medial and lateral support and great patient comfort. The ventilated material reduces moisture build up whilst the towel lining improves comfort

The brace and hinge casings are fully bound to enhance the overall look of the brace whilst the hook/loop straps have grip-taps to assist patients with limited dexterity.

Size	Code	Knee Circumference
XS	*/02	25 - 28 cm
S	*/03	28 - 34 cm
M	*/05	34 - 40 cm
L	*/07	40 - 46 cm
XL	*/08	46 - 56 cm
XXL	*/09	56 - 61 cm

* Insert product code



SK Knee Brace

Product Code: **SK**



The front opening design and easy-open tabs on the hook/loop straps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity. Hyperstop hinges provide medial and lateral support for the unstable knee, while allowing flexion.

Size	Code	NDC Code	Knee Circumference
XS	*/02	-	25 - 28 cm
S	*/03	052479	28 - 34 cm
M	*/05	052455	34 - 40 cm
L	*/07	052516	40 - 46 cm
XL	*/08	052509	46 - 56 cm
XXL	*/09	-	56 - 61 cm

* Insert product code



Nuffield SK Knee Brace

Product Code: **NSK**

The front opening design and easy-open tabs on the hook/loop straps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity.

The brace incorporates a polycentric hinge which follows the true motion of the knee. The top opening design allows easy removal of the hinges for laundering.

Size	Code	Knee Circumference
XS	*/02	25 - 28 cm
S	*/03	28 - 34 cm
M	*/05	34 - 40 cm
L	*/07	40 - 46 cm
XL	*/08	46 - 56 cm

* Insert product code



VK Knee Brace

Product Code: VK

The longer 41cm metal hyperstop hinges provide increased stability and medial and lateral support for the unstable valgus or varus knee. The front opening design and easy-open tabs on the hook/loopstraps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity.

The VK Knee brace is sided left and right for improved fit.

Size	Code (Left)	Code (Right)	Knee Circumference
XS	*/12	*/22	25 - 28 cm
S	*/13	*/23	28 - 34 cm
M	*/15	*/25	34 - 40 cm
L	*/17	*/27	40 - 46 cm
XL	*/18	*/28	46 - 56 cm
XXL	*/19	*/29	56 - 61 cm

* Insert product code



Autolock II

Product Code: AUT

- A front opening design with patella flap to prevent the varus/valgus knee protruding.
- Distal hinge insertion for ease of hinge removal.
- Four circumferential straps ensure a good fit and Easy Open straps on the hook/loop straps make this brace suitable for the patient with limited dexterity.
- Available in beige Neoprene.

Size	Code (Left)	Code (Right)	Knee Circumference
S	*/13	*/23	28 - 34 cm
M	*/15	*/25	34 - 40 cm
L	*/17	*/27	40 - 46 cm
XL	*/18	*/28	46 - 56 cm

* Insert product code



Lycrafleece Patella Lux

Product Code: BLPL

Provides medial or lateral stabilisation of the knee. The crescent shaped buttress pad can be placed on the medial or lateral side of the patella. A non slip rubber lining is used to ensure the support does not move after fitting.

The Patella Lux is ideal for lateral or medial subluxation of the patella.

Universal design fits left and right.

Size	Code	Knee Circumference
XS	*/02	25 - 28 cm
S	*/03	28 - 36 cm
M	*/05	33 - 41 cm
L	*/07	38 - 46 cm
XL	*/08	45 - 55 cm

* Insert product code



Economy Knee Sleeve

Product Code: **EN**

A Neoprene sleeve with open patella providing warmth and light compressional support.

Size	Code	NDC Code	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/02	-	30 cm	30 - 33 cm	38 cm
S	*/03	052523	33 cm	33 - 36 cm	46 cm
M	*/05	052493	38 cm	36 - 38 cm	50 cm
L	*/07	052462	41 cm	38 - 43 cm	53 cm
XL	*/08	052486	43 cm	43 - 48 cm	58 cm
XXL	*/09	-	52 cm	48 - 56 cm	61 cm

* Insert product code

Standard Knee Sleeve - Open Patella

Product Code: **SNO**

A Neoprene sleeve with a reinforced open patella section providing warmth and light compressional support.

Size	Code	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/02	30 cm	30 - 33 cm	38 cm
S	*/03	33 cm	33 - 36 cm	46 cm
M	*/05	38 cm	36 - 38 cm	50 cm
L	*/07	41 cm	38 - 43 cm	53 cm
XL	*/08	43 cm	43 - 48 cm	58 cm
XXL	*/09	52 cm	48 - 56 cm	61 cm

* Insert product code

Spiral Knee Sleeve

Product Code: **SP**

A Neoprene sleeve with two flexible spirals, providing light compression and medial/lateral support.

Size	Code	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/02	30 cm	30 - 33 cm	38 cm
S	*/03	33 cm	33 - 36 cm	46 cm
M	*/05	38 cm	36 - 38 cm	50 cm
L	*/07	41 cm	38 - 43 cm	53 cm
XL	*/08	43 cm	43 - 48 cm	58 cm
XXL	*/09	52 cm	48 - 56 cm	61 cm

* Insert product code



Buttress Knee Sleeve

Product Code: BN

A standard Neoprene knee sleeve with open patella, incorporating a buttress pad. The Buttress Knee Sleeve can be selected with a buttress pad to oppose medial/lateral patella migration; caused by ligament deficiency.

The buttress pad positions are defined as viewed by the patient, for example a left buttress pad would give lateral support to the left leg and medial support to the right leg. A full buttress is also available.

Size	Code (Left)	Code (Right)	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/12	*/22	30 cm	30 - 33 cm	38 cm
S	*/13	*/23	33 cm	33 - 36 cm	46 cm
M	*/15	*/25	38 cm	36 - 38 cm	50 cm
L	*/17	*/27	41 cm	38 - 43 cm	53 cm
XL	*/18	*/28	43 cm	43 - 48 cm	58 cm
XXL	*/19	*/29	52 cm	48 - 56 cm	61 cm

* Insert product code



Stabilised Knee Sleeve

Product Code: ST

A pull-on Neoprene sleeve incorporating a buttress pad, spiral stays for medial and lateral support, and hook/loop straps above and below the knee for improved fit. The Stabilised Knee Sleeve can be selected with a buttress pad to oppose medial/lateral patella migration; caused by ligament deficiency.

The buttress pad positions are defined as viewed by the patient, for example a left buttress pad would give lateral support to the left leg and medial support to the right leg. A full buttress is also available.

Size	Code (Left)	Code (Right)	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/12	*/22	30 cm	30 - 33 cm	38 cm
S	*/13	*/23	33 cm	33 - 36 cm	46 cm
M	*/15	*/25	38 cm	36 - 38 cm	50 cm
L	*/17	*/27	41 cm	38 - 43 cm	53 cm
XL	*/18	*/28	43 cm	43 - 48 cm	58 cm
XXL	*/19	*/29	52 cm	48 - 56 cm	61 cm

* Insert product code



Hinged Knee Sleeve

Product Code: HN

A pull-on Neoprene sleeve, with strong metal hyperstop hinges to control medial and lateral instability.

Hook/loop straps above and below knee give improved fit.

Size	Code	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/02	30 cm	30 - 33 cm	38 cm
S	*/03	33 cm	33 - 36 cm	46 cm
M	*/05	38 cm	36 - 38 cm	50 cm
L	*/07	41 cm	38 - 43 cm	53 cm
XL	*/08	43 cm	43 - 48 cm	58 cm
XXL	*/09	52 cm	48 - 56 cm	61 cm

* Insert product code



T Scope Hip Brace

The Breg T Scope Post-Op Hip brace delivers a patient-centric design to provide unprecedented comfort, simplicity and support during post-operative hip rehabilitation. Developed in conjunction with industry-leading physicians specializing in injuries, instabilities, and arthroscopic surgery of the hip, the T Scope Post-Op Hip brace incorporates Breg's premium and time-tested technologies.

It features Breg's telescoping design for a fast, universal fit; the patented T Scope Range-of-Motion hinge with quick-adjusting stops; a Quick Lock to lock out brace positions with the push of a button; and a unique waist compression lacing system that gives patients an easy one-handed pull to provide an intimate fit and offer additional back support. Padding throughout the brace, and thoughtful details, like cushioned hand protection when using crutches and easy to use features, help ensure patient compliance.

Size	Code
08910	T Scope Hip, Left, 31"– 60" Waist
08911	T Scope Hip, Right, 31"– 60" Waist
08912	T Scope Hip Small, Left, 24"– 31" Waist
08913	T Scope Hip Small, Right, 24"– 31" Waist
08914	T Scope Hip Full Kit, Left, 31"– 60" Waist
08915	T Scope Hip Full Kit, Right, 31"– 60" Waist
08916	T Scope Hip Full Kit Small, Left, 24"– 31" Waist
08917	T Scope Hip Full Kit Small, Right, 24"– 31" Waist



Hipo

This modular hip orthosis system offers functional and effective support in the treatment of various hip pathologies.

- Preformed modular plastic waist and thigh components
- Combination Kinesio Joint, controls flexion / extension, abduction / adduction
- Hip and waist components supplied with comfort liners
- Component sizes can be mixed for optimum fit.
- Kinesio Offset hinge with abduction/adduction settings (7.5° increments) and flexion/extension settings (from 0° to 120° in 15° increments)

Hip / Waist Component

Size	S	M	L	XL
Circ (Waist)	60 - 82 cm	80 - 96 cm	94 - 115 cm	99 - 120 cm
Circ (Hip)	68 - 90 cm	90 - 110 cm	106 - 126 cm	111 - 131 cm
Code (Left)	M.2122	M.2123	M.2124	M.2125
Code (Right)	M.2112	M.2113	M.2114	M.2115

Kinesio Hinge

Size	Height	Code
Reg	< 180 cm	M.2303
Long	> 181 cm	M.2304

Thigh Component

Size	S	M	L	XL
Circ (Waist)	44 - 53 cm	48 - 60 cm	57 - 68 cm	62 - 76 cm
Circ (Hip)	34 - 41 cm	36 - 45 cm	42 - 49 cm	46 - 56 cm
Code (Left)	M.2222	M.2223	M.2224	M.2225
Code (Right)	M.2212	M.2213	M.2214	M.2215



Promedics Hip Joint

Product Code: **PHJ/10**

The Promedics Hip Joint is 40cm in length and allows control of abduction (between -5 to +25°) with variable flexion and extension stops. The abduction axis may be set at various positions while the brace is on the patient, by simply turning the central adjustment screw with the allen key provided.

Universal Hip Abduction Waist Belt

Product Code: **JBB/10**

Universal Hip Abduction Belt with Shoulder Strap

Product Code: **JBBS/10**



Ankle & Footwear



NEW



airstep™

Promedics Airstep™ Walker

Product Codes: **ASW** (Tall) / **ASSW** (Short)

Airstep™ is an innovative walker with a fully integrated inflation system that allows independent adjustment.

- Super lightweight - ventilated shell system
- Breathable fitted contoured liner
- Low profile sole unit - minimal height discrepancy
- Easy three strap closure
- Fast and easy to fit
- No separate pump required
- Excellent stability
- Achilles wedges available

Size	Code (Tall)	Shoe Size
XS	ASW/02	1.5 - 3
S	ASW/03	3 - 5.5
M	ASW/05	6 - 9
L	ASW/07	9.5 - 11
XL	ASW/08	11 +

Size	Code (Short)	Shoe Size
XS	ASSW/02	1.5 - 3
S	ASSW/03	3 - 5.5
M	ASSW/05	6 - 9
L	ASSW/07	9.5 - 11
XL	ASSW/08	11 +



Dual integrated pumps
for fast, easy inflation



Short version



Scan QR Code to watch the
Airstep Walker Instructional video
on our Youtube channel.



Low profile, hard wearing sole unit
Improved gait function



Adjustable / removable toe guard

Genesis Full Shell Walker

Product Codes: **BL5150xx** (3-Strap) / **BL5140xx** (4-Strap)

The Genesis Walker is the lightest full shell pneumatic walker on the market. It has been computer engineered to provide strength where required while eliminating extra weight. Similar to popular athletic shoes, the Genesis Walker uses a dynamic air chamber in the heel of the tread which actively absorbs the shock of heel strike and allows for a more comfortable walking experience.

Common Examples of Use

- Ankle sprains and strains
- Forefoot sprains
- Stable foot fractures
- Distal tibial fractures

Product Features

- Built-in air pump
- Durable hard shell design
- Trim-able shell
- Narrow heel width
- Dynamic air stride technology
- Dual air chambers
- Scalloped no-slip tread

Size	Code	Shoe Size
XS	*01	1 - 2
S	*03	3 - 5
M	*05	6 - 9
L	*07	10 - 12
XL	*09	12 +

* Insert product code

Genesis Full Shell Walker (Short)

Product Code: **BL52500xx**

The mid-calf version of the cutting-edge Genesis Full Shell Walker that includes all of the same great features. Designed to be even lighter and more maneuverable, the Genesis Mid-Calf is ideal for those patients who require comfortable foot support.



Ankle & Footwear



Dual air chambers



Durable hard shell design



Built in pump



Cut line trimable shell



DUAL AIRBAGS
- INDEPENDENTLY ADJUSTABLE



JURA Medical



Jura Air Shield Walker

Product Code: JWA

Available with charcoal or red liners, this lightweight low profile walker allows good gait pattern and has independently adjustable dual airbags for improved stability and support.

Features a comfortable breathable liner with toe cover. Anterior & posterior shells offer improved stability and immobilisation with velour strap fastening. The shells can be removed for staged rehabilitation.

Replacement liners available in charcoal.

Jura Achilles Wedges

Product Code: JWW

Extension straps included.

Left	Right
JWW/RL	JWW/RR



Jura Air Shield Walker - Charcoal

Size	Code	Replacement Liners Code	Shoe Size
S	*/03/C	*/03/L/C	3 - 5.5
M	*/05/C	*/05/L/C	6 - 9
L	*/07/C	*/07/L/C	9.5 - 11
XL	*/08/C	*/08/L/C	11 +

* Insert product code

Jura Air Shield Walker - Red

Size	Code	Shoe Size
S	*/03/R	3 - 5.5
M	*/05/R	6 - 9
L	*/07/R	9.5 - 11
XL	*/08/R	11 +

* Insert product code

Jura Air Shield Walker - Short

Product Code: JWAS

Only available with charcoal liners.

Size	Code	Replacement Liners Code	Shoe Size
S	*/03/C	*/03/L/C	3 - 5.5
M	*/05/C	*/05/L/C	6 - 9
L	*/07/C	*/07/L/C	9.5 - 11
XL	*/08/C	*/08/L/C	11 +

* Insert product code

Jura Air Shield Walker Extension straps

Supplied as standard with Jura Achilles Wedges.

Code
JWA-ES



JURA Medical





Jura Fixed Walker

Product Code: JWF

One piece injection moulded walker with 'paddle' uprights to reduce pressure on lower limb and provide improved stability and immobilisation. Foot section can accommodate full dressings and padding.

New lightweight design sole unit to greatly improve shock absorption, gait pattern to reduce the overall weight of the walker to an absolute minimum.

Reduced strap fastening to simplify removal / application and improve patient acceptance. Comfortable detachable liners with replacements available.

Available with standard or Air Liner.

Fixed Walker

Size	Code	Replacement Liners Code	Shoe Size
XS	*/02	*/02/L	1.5 - 3
S	*/03	*/03/L	3 - 5.5
M	*/05	*/05/L	6 - 9
L	*/07	*/07/L	9.5 - 11
XL	*/08	*/08/L	11 +

* Insert product code

Fixed Walker with Air Liner

Size	Code	Replacement Liners Code	Shoe Size
XS	*/02/A	*/02/L/A	1.5 - 3
S	*/03/A	*/03/L/A	3 - 5.5
M	*/05/A	*/05/L/A	6 - 9
L	*/07/A	*/07/L/A	9.5 - 11
XL	*/08/A	*/08/L/A	11 +

* Insert product code



Jura Fixed Walker - Short

Product Code: JWFS

Short version of the Jura Fixed Walker.

Size	Code	Replacement Liners Code	Shoe Size
XS	*/03	*/03/L	3 - 5.5
S	*/03	*/03/L	3 - 5.5
M	*/05	*/05/L	6 - 9
L	*/07	*/07/L	9.5 - 11
XL	*/08	*/08/L	11 +

* Insert product code



ROM Walker

Product Code: **JWR**

Fully moulded plastic foot section with rubber sole and padded insole.

Easily adjustable 'Push-Dial-Set' range of motion settings. Chosen settings can be 'locked' in by addition of tie wrap through flexion/extension controls.

Plantarflexion Settings: 0, 10, 30, 45 degrees.

Dorsiflexion Settings: 0, 10, 30 degrees

The lightweight ROM walker features detachable foam liners and return-loop hook / loop Fastening. Replacement liners available.

Size	Code	Replacement Liners Code	Shoe Size
XS	*/02	*/02/L	1 - 3
S	*/03	*/03/L	3 - 5.5
M	*/05	*/05/L	6 - 9
L	*/07	*/07/L	9.5 - 11
XL	*/08	*/08/L	11 +

* Insert product code



Dorsal Night Splint

Product Code: **1144**

Designed to help provide support while sleeping for patients with plantar fasciitis and Achilles tendonitis.

Common Examples of Use

- Plantar fasciitis
- Plantar flexion contracture
- Achilles tendonitis

Product Features

- Lightweight, low profile design increases comfort
- Helps patient sleep more comfortably therefore increasing compliance
- Slips onto the foot and is easily adjusted with hook closure
- Fits either left or right foot

Size	Code	Shoe Size
S/M	*2	5 - 8
L/XL	*4	9 - 14

* Insert product code





Podus Boot

Product Code: **PB**

Complete pressure free heel suspension support, lined with Anti-Decubitus Kodel type liner to prevent pressure sores. The heel features a rotator bar which, when correctly positioned, corrects internal / external rotation of the hip, foot and ankle. The liner is removable for washing and can be replaced. Measure length of foot from back of calcaneus to the tip of the hallux.

Code	Item / Size
*/STDS	Podus Boot Complete Standard, (With Straps)
*/XL	Podus Boot Complete Ex Large, fits 25 - 33cm foot
*/XLS	Podus Boot Complete Ex Large, (With Straps)
*/STD-L	Replacement Liner for Standard Podus Boot
*XL-L	Replacement Liner for Ex Large Podus Boot

* Insert product code



Boot Foot Support

Product Code: **BFS**

A simple foot support to prevent drop foot. Fleece lined and snug fit for added comfort. The lining can be removed for washing or replaced. One size fits all. Only available in Grey.

Code	Item / Size
*G/10	Boot Foot Support Grey
*G/10-L	Replacement Liner Grey

* Insert product code



Heel Relief Orthosis

Product Code: **ANTFORTH**

Features an adjustable bar to control external/internal rotation.

A fully aerated heel and adjustable sole plate make this a very effective and comfortable brace.

Universal left and right. Replacement liners available.

Size	Code	Replacement Liner Code	Shoe Size
S	*/S	*/S-L	Up to 5
M	*/M	*/M-L	6 - 9
L	*/L	*/L-L	9.5 +

* Insert product code



Universal Plantar Fasciitis Orthosis

Product Code: UFO

Holds foot in dorsiflexion and allows load application to soft tissues of plantar fascia and achilles tendon. Dorsi flexion easy to adjust using double hook/loop straps additional wedges (3 supplied)

Fits left and right. Replacement liners available.

Size	Code	Replacement Liners Code	Shoe Size
S	*/S	*/S-L	3 - 5.5
M	*/M	*/M-L	6 - 9
L	*/L	*/L-L	9.5 - 11

* Insert product code



Ankle Contracture AFO

Product Code: LEB

A quality AFO constructed from lightweight polycarbonate that provides heel suspension to assist in the treatment and prevention of pressure sores. The clear polycarbonate frame allows the heel to be monitored without the need to remove the orthosis. The AFO is also effective in the management of foot drop and non-fixed plantar flexion contractures.

The liner is treated to wick away moisture and provide a drier cleaner environment thus reducing the risk of bacterial build up and improve patient compliance. There is a padded anti-rotation bar for positioning and prevention of lateral/medial hip rotation. The wide design of the frame also provides lateral ankle protection if patients roll onto their side during sleep. Adjustable toe rest and removable anti-slip sole unit

Size	Code	Size (Liner only)	Code
Small	LEB/S	Small	LEB/S/L
Regular	LEB/R	Regular	LEB/R/L



Multi Use Boot

Product Code: LEM

An adjustable night resting splint providing passive stretch for the treatment of plantar fasciitis and Achilles tendonitis. Adjustable closure straps and dorsiflexion- assist straps. Padded liner for patient comfort. Non-slip coating on sole unit to assist in patient transfers.

Size
Universal

Heal-up® Protection Boot

Product Code: **HUP/10**

- Offloads heel to prevent and manage ulcers
- Simple and effective closure straps
- Floats the heel
- Open heel provides visibility with additional opening for greater inspection visibility
- Universal size, left or right.

Size	Code	Code
Universal	1	HUP/10
Universal	Box of 10	NH/HUP/10

* Insert product code



Treated with **Agion®** and **Silvadur™** antimicrobial to inhibit the growth of bacteria.



Numbered tabs for easy 1 - 2 - 3 application




heal-up.co.uk



Anatomical strapping with configuration to inhibit "drop-foot"

Ankle & Footwear

Friction free textile



Additional opening for improved inspection

Floated heel

HEALUP®



Swedish Ankle Foot Orthosis

Product Code: **M.9050**

Provides static dorsiflexion assistance and lateral stability for the ankle/foot area. Made in polythene, it is thicker on the vertical aspect for rigidity and thinner on the footplate to allow for trimming.

The low arch and open heel give a neat shape which fits into a shoe easily and the calf section is padded for comfort.

Can be trimmed to size with scissors and shaped with a heat gun if necessary.

Size	Code (Left)	Code (Right)	Shoe Size
S	*/SL	*/SR	3 - 6
M	*/ML	*/MR	7 - 8.5
L	*/LL	*/LR	7.5 - 10

* Insert product code



Fixed Ankle

Product Code: **821**

A lightweight, polypropylene drop foot splint which provides dorsiflexion assistance and prevents plantarflexion.

The splint is heat mouldable.

Size	Code (Left)	Code (Right)	Shoe Size
S	*0/0	*5/0	4 - 6
M	*1/0	*6/0	5 - 7
L	*2/0	*7/0	6 - 8
XL	*3/0	*8/0	7 - 9

* Insert product code



Leaf Spring Orthosis

Product Code: **825**

A lightweight, polypropylene splint which provides static dorsiflexion assistance and allows limited plantarflexion.

The low arch and open heel give a neat shape which fits into a shoe easily, the splint is also heat mouldable.

Size	Code (Left)	Code (Right)	Shoe Size
S	*0/0	*5/0	4 - 6
M	*1/0	*6/0	5 - 7
L	*2/0	*7/0	6 - 8
XL	*3/0	*8/0	7 - 9

* Insert product code



Air/Gel Ankle

Product Code: **AGA/10**

Two plastic medial and lateral ankle supports, shaped to fit comfortably over the malleoli, are lined with closed-cell foam. The two detachable air/gel bags line the splint and can be used for hot or cold therapy.

The plastic supports are attached to a 'stirrup' of closed cell foam, fixed with hook/loop and worn under the heel.

One size; universal design fits left and right.

Size
Universal



Barnsley Ankle Brace

Product Code: **EWB**

The stitched forefoot section allows this brace to slip on like a sock, and the hook/loop closure applies variable tension for comfort.

Made in 10cm quality elastic, giving light ankle support.

Contains Latex

Size	Left	Right	Shoe Size
S	*/13	*/23	3 - 5
M	*/15	*/25	5 - 7
L	*/17	*/27	7 - 9

* Insert product code



Plantar Fasciitis Wrap

Product Code: **1136x**

The Plantar Fasciitis Wrap provides support for patients with plantar fasciitis, is low profile and fits inside most shoes for daily wear.

Contains Latex

Size	Code	Ankle Circumference
S	*2	28 cm
M	*3	30 cm
L	*4	33 cm

* Insert product code





Figure of 8 Ankle Brace

Product Code: **F8LA**

A lightweight ankle support which laces up and features an elastic figure of 8 strap for added support and compression. Fits either ankle.

Contains Latex

Size	Code	Shoe Size
XS	*/02	1 - 3
S	*/03	3 - 5
M	*/05	6 - 8
L	*/07	9 - 12
XL	*/08	12 +

* Insert product code



Malleostromg

Product Code: **M.A1028**

Ankle brace for limiting pro- and supination. Thin breath-o-prene structure with padded medial and lateral shells. Additional rear padding for improved comfort. Circular strap for limiting talar shift (can be shortened to fit)

- Functional therapy of 2nd degree distortions
- Functional recovery following surgical/conservative treatment of ligament injuries (3rd degree distortions)
- Functional recovery of consequences of malleolar fractures (surgically/conservatively treated)

Contains Latex

Size	Code (Left)	Code (Right)	Shoe Size
XS	*/XSL	*/XSR	1 - 3
S	*/SL	*/SR	3 - 5
M	*/ML	*/MR	6 - 8
L	*/LL	*/LR	9 - 12
XL	*/XLL	*/XLR	12 +

* Insert product code



Hallux Valgus Night Splint

Product Code: **HVN**

A plastic semi flexible splint with padded closed cell foam toe section. Features adjustable forefoot pad and adjustable hook/loop fastening.

Size	Code (Left)	Code (Right)
XS	*/12	*/22
S	*/13	*/23
M	*/15	*/25
L	*/17	*/27

* Insert product code



Universal Soft Hallux Valgus

Product Code: **JSHV**

A soft padded universal size splint that helps maintain the large toe in a neutral position. Easy to apply with adjustable hook/loop fastening, adjustable aluminum stay provides stability and assists in the positioning of the alignment of the large toe. The splint can be worn at night or during the day in footwear

Size	Code (Left)	Code (Right)
Universal	JSHV/L	JSHV/R



Anklet

Product Code: **ANK**

A pull-on Neoprene support for light compression and ankle support.

Size	Code	Shoe Size
S	*/03	3 - 5
M	*/05	5 - 7
L	*/07	7 - 9
XL	*/08	9 - 11

* Insert product code



Lincoln Ankle Brace

Product Code: **LAB**

A zipped ankle support in Ventoprene, which is ventilated, towelling-lined Neoprene. The zip is worn on the lateral aspect of the ankle, and a flexible spiral stay gives good ankle support on the medial aspect of the ankle.

The 'figure of eight' strapping is designed to limit plantar flexion.

Size	Code (Left)	Code (Right)	Shoe Size
S	*/13	*/23	3 - 5
M	*/15	*/25	5 - 7
L	*/17	*/27	7 - 9

* Insert product code



McDavid Ankle

Product Code: **MCD**

The McDavid Ankle Guard offers outstanding support, protection, and comfort in a laced ankle brace with several features including:

- Durable vinyl, nylon fabric
- New ventilated tongue
- Spring steel stays for support
- Sewn-in arch support
- Removable medial & lateral plastic inserts for added rigidity

Size	Code	Shoe Size
XS	*/02	3 - 4
S	*/03	5 - 6
M	*/05	7 - 9
L	*/07	9 - 11
XL	*/08	11 +

* Insert product code



The coloured graphic depicts the pressure distribution on the sole of the foot.



Jura Post-Op Shoe

Product Code: **P-OP**

- Square toe configuration offers protection and improved fit/ accommodation.
- Rocker-sole profile aids gait and reduces forces on foot.
- Outsole is non-slip and material composition provides reduced movement under the metatarsal heads.
- Simple closure system allows for additional dressings.
- Insole can be removed for alteration, improved accommodation or fitting 'Off-Loading Insole'

Mens Size	Code	Shoe Size	Ladies Size	Code	Shoe Size
S	*/MS	5 - 7	S	*/WS	3 - 5
M	*/MM	7.5 - 9	M	*/WM	5.5 - 7
L	*/ML	9.5 - 11	L	*/WL	7.5 - 9
XL	*/XL	11.5 - 13			

* Insert product code




Jura High Post-Op Shoe

Product Code: **HP-OP**

- Square toe configuration offers protection and improved fit/ accommodation.
- Rocker-sole profile aids gait and reduces forces on foot.
- Outsole is non-slip and material composition provides reduced movement under the metatarsal heads.
- Simple closure system allows for additional dressings.
- Insole can be removed for alteration, improved accommodation or fitting 'Off-Loading Insole'
- Adjustable instep strap for improved seating of heel.
- Ankle high upper for improved fit and control for patients with compression dressings and total contact casts.

Mens Size	Code	Shoe Size	Ladies Size	Code	Shoe Size
S	*/MS	5 - 7	S	*/WS	3 - 5
M	*/MM	7.5 - 9	M	*/WM	5.5 - 7
L	*/ML	9.5 - 11	L	*/WL	7.5 - 9
XL	*/XL	11.5 - 13			

* Insert product code



The coloured graphic depicts the pressure distribution on the sole of the foot.




The coloured graphic depicts the pressure distribution on the sole of the foot.

Off-loader Heel Shoe

Product Code: **OWH**

- Hard wearing padded nylon upper.
- Adjustable hook / loop fastening.
- Hard-wearing sole unit shaped to prevent force on forefoot.
- Extendable forefoot closure panel to accommodate dressings padding.
- Lower profile wedge.
- Insole can be removed for alteration, improved accommodation or fitting 'Off-Loading Insole'

Size	Code	Shoe Size
XS	OWH/XS	2 - 4
S	OWH/S	5 - 6
M	OWH/M	7 - 8
L	OWH/L	9 - 11
XL	OWH/XL	12 - 14



Off-loading Insoles

- Easily remove sections of insole for custom off-loading in treatment of various foot pathologies.
- Greatly reduces pressure on plantar surface when patient is walking.
- Self-adhesive stabilising base board.
- Available for all Jura post op shoes.



Off-loading insoles for Post-op & High Post-op Shoes

Mens Size	Code	Shoe Size	Womans Size	Code	Shoe Size
S	P-OP/MSI	5 - 7	S	P-OP/WSI	3 - 5
M	P-OP/MMI	7.5 - 9	M	P-OP/WMI	5.5 - 7
L	P-OP/MLI	9.5 - 11	L	P-OP/WLI	7.5 - 9
XL	P-OP/XLI	11.5 - 13			

Off-loading insoles for Off-loader Heel Shoe

Size	Code	Shoe Size
XS	OWH/XSINS	2 - 4
S	OWH/SINS	5 - 6
M	OWH/MINS	7 - 8
L	OWH/LINS	9 - 11
XL	OWH/XLINS	12 - 14



NEW

Heel Relief Shoe

Product Code: **HRS**




- Alleviates pressure on the heel
- Comfortably padded
- Insole can be removed for alteration or to accommodate a custom orthotic
- Ankle strap reduces heel slippage
- Removable / adjustable closure to accommodate dressings etc

Common Examples of Use

- Trauma to the rear of the foot
- Wounds or ulcerations on the heel area
- Post surgical healing
- Pain or bruising to the heel

Size	Code	Shoe Size (Women's)	Shoe Size (Men's)
XS	*/XS	4 - 7.5	-
S	*/S	8 - 10	6 - 7.5
M	*/M	10.5 - 13	8 - 10
L	*/L	-	10.5 - 12
XL	*/XL	-	12.5 +

* Insert product code



NEW

Comfortable padded heel counter

Return Loop buckles

SuperStep™



Durable hard-wearing sole unit



Repositionable hook tabs

Superstep Post Op Shoe

Product Code: **RPO**

A padded durable nylon upper with return loop buckles.

The upper will accommodate swelling and dressings.

Padded insole with fully rigid, lightweight hard wearing outer sole unit.

The sole unit is designed for improved stability to aid rehabilitation following foot surgery, trauma or any condition that requires immobilisation to aid healing.

Size	Code	Sole Length	Sole Width
Paediatric	*/PAED	20 cm	8 cm
S (Ladies)	*/WS	24 cm	8 cm
M (Ladies)	*/WM	25 cm	9 cm
L (Ladies)	*/WL	27 cm	9 cm

Size	Code	Sole Length	Sole Width
S (Gents)	*/MS	28 cm	10 cm
M (Gents)	*/MM	29 cm	10 cm
L (Gents)	*/ML	30 cm	10 cm
XL	*/MXL	32 cm	11 cm

* Insert product code



Hi Sided Cast Sandal

Product Code: **HSCS**

Designed with hard wearing, non-stretch nylon faced padded upper for long term wear. Loop-lock D-ring fastening for secure and easy fastening. Non-slip rocker sole.

Also available in paediatric sizes, see pg. 116.

Size	Code	Length	Width
XS	*/AXS	23 cm	8 cm
S	*/AS	24.5 cm	9 cm
M	*/AM	26.5 cm	10 cm
L	*/AL	29 cm	10.5 cm
XL	*/AXL	31 cm	11 cm

* Insert product code



Open Toe Cast Sandal with Buckles

Product Code: **CS/30**

Designed for long term wear, a strong moulded sole is stitched to a hard-wearing canvas upper. The return loop buckles ensure a good fit on the patient.

For use with patients in a weight bearing cast. Colours may vary.

Size	Code	Sole Length	Sole Width
XS	*XS	18 cm	9 cm
S	*S	23 cm	10 cm
M	*M	27 cm	11.5 cm
L	*L	31 cm	13 cm
XL	*XL	36 cm	14 cm

* Insert product code



Medi Shoe

Product Code: **811**

A quality padded upper with return loop buckles accommodates the oedematous or bandaged foot very comfortably. The ridged non-slip sole is semi-flexible to allow normal gait, and lined with closed cell foam for comfort.

Size	Code	Sole Length	Sole Width
XS	*42	20 cm	8 cm
S	*43	23 cm	8 cm
M	*47	27 cm	9 cm
L	*35	29 cm	10 cm
XL	*37	31 cm	12 cm

* Insert product code

Basic Post-Op Shoe

Product Code: **BPO**

A padded nylon upper with return loop buckles accommodates the oedematous or bandaged foot very comfortably. The rigid non-slip sole is semi flexible to allow normal gait.



Size	Code	Sole Length	Sole Width	NDC Code
XS	*XS	20 cm	8 cm	171965
S	*S	23 cm	8 cm	165599
M	*M	27 cm	9 cm	165582
L	*L	29 cm	10 cm	165575
XL	*XL	31 cm	12 cm	165568

* Insert product code



Open Toe Cast Sandal

Product Code: **CS/08**

A hard-wearing canvas upper with hook/loop strapping is stitched to a rigid non-slip rocker sole. For use with patients in a weight-bearing cast.



Size	Code	NDC Code	Sole Length	Sole Width
XS	*XS	020744	18 cm	9 cm
S	*S	010653	23 cm	10 cm
M	*M	010646	27 cm	11.5 cm
L	*L	020751	31 cm	13 cm
XL	*XL	-	36 cm	14 cm

* Insert product code



Open Toe Cast Sandal with Buckles

Product Code: **CS/48**

Return loop buckles ensure a good fit on the patient. A hard-wearing canvas upper with hook/loop strapping is stitched to a rigid non-slip rocker sole. For use with patients in a weight bearing cast.



Size	Code	NDC Code	Sole Length	Sole Width
XS	*/XS	-	18 cm	9 cm
S	*/S	027101	23 cm	10 cm
M	*/M	027217	27 cm	11.5 cm
L	*/L	027224	31 cm	13 cm
XL	*/XL	-	36 cm	14 cm

* Insert product code





Size	Code	NDC Code	Sole Length	Sole Width
XS	*/XS	-	22 cm	7 cm
S (Woman)	*/WS	091102	24 cm	8 cm
M (Woman)	*/WM	090938	25 cm	9 cm
L (Woman)	*/WL	091157	27 cm	9 cm

Supershoe

Product Code: **442**

A padded nylon upper with return loop buckles accommodates the oedematous or bandaged foot comfortably and the removable tongue can be positioned anywhere on the inner for comfort.

A moulded heel counter provides a smooth stable contour around the heel.

The semi-rigid non-slip sole is strong and will not flex when walking; useful when toes need to be immobilised eg. after Keller's surgery.

A pressure sensitive longitudinal arch support is included for comfort.

Size	Code	NDC Code	Sole Length	Sole Width
S (Man)	*/MS	091140	28 cm	10 cm
M (Man)	*/MM	091027	29 cm	10 cm
L (Man)	*/ML	091041	30 cm	10 cm
XL	*/XL	-	32 cm	11 cm



Rigid Post Op Shoe

Product Code: **445**

A padded durable nylon upper with return loop buckles. The upper will accommodate swelling and dressings.

Padded insole with fully rigid lightweight composite outer sole unit. The sole unit is rigid to prevent any flexion following foot surgery that requires immobilisation to aid healing and rehabilitation.

Size	Code	Sole Length	Sole Width
S (Ladies)	*/WS	24 cm	8 cm
M (Ladies)	*/WM	25 cm	9 cm
L (Ladies)	*/WL	27 cm	9 cm

Size	Code	Sole Length	Sole Width
S (Gents)	*/MS	28 cm	10 cm
M (Gents)	*/MM	29 cm	10 cm
L (Gents)	*/ML	30 cm	10 cm
XL	*/XL	32 cm	11 cm



Dress Shoe Comfort Inserts

- Visco elastic polymer providing excellent absorption of shock and shear force
- Fits easier into thinner footwear, eg dress shoes
- Unique shaping allows multiple functions
- Scallop pad cushion MT heads
- Heel rim band absorbs and dissipates shock/shear
- Heel cup area aids painful heel and heel spur area

Size	Shoe Size (Ladies)	Code	Shoe Size (Gents)	Code
S	3 - 5	L-DSF1	5 - 7	L-DSM1
M	6 - 7	L-DSF2	8 - 10	L-DSM2
L	8 - 10	L-DSF3	11 - 13	L-DSM3

PQ - 6mm Heel Lift



- Visco elastic polymer
- Shock absorbing properties
- 6mm raise for shortage accommodation
- sold in single units

Shoe Size	Code
3 - 5	L-PQ61
6 - 8	L-PQ62
9 - 11	L-PQ63

PQ - 9mm Heel Lift



- Visco elastic polymer
- Shock absorbing properties
- 9mm raise for shortage accommodation
- sold in single units

Shoe Size	Code
3 - 5	L-PQ91
6 - 8	L-PQ92
9 - 11	L-PQ93



Star Heel Pad

- 6mm thick Visco elastic heel pad
- Unique star raise pattern improves dissipation of forces
- Well tapered towards pad edge

Size	Shoe Size	Code
S	3 - 5	L-SHPS
M	6 - 8	L-SHPM
L	9 - 11	L-SHPL



Star Spur Pad

- 6mm thick Visco elastic heel pad
- Unique star raise pattern improves dissipation of forces
- Well tapered towards pad edge
- Heel relief reduces pain associated with heel spurs and plantar fasciitis

Size	Shoe Size	Code
S	3 - 5	L-SSPS
M	6 - 8	L-SSPM
L	9 - 11	L-SSPL



Epoda Dual Heel Cups

Product Code: **EPCB**

Dual density design with heel cradle for extra comfort and to aid foot positioning.

Size	Shoe Size	Code
S	4 - 5	*/03
M	6 - 7	*/05
L	8 - 9	*/07
XL	10 - 11	*/08

* Insert product code



Epoda Dual Heel Pads

Product Code: **EPB**

Dual density wedge shaped design with integral extra soft blue area.

Size	Shoe Size	Code (5mm Thick)	Code (10mm Thick)
S	4 - 5	*/53	*/103
M	6 - 7	*/55	*/105
L	8 - 9	*/57	*/107
XL	10 - 11	*/58	*/108

* Insert product code

JURAGel

Jura Gel is a therapeutic range of pure silicone gel containing mineral oil. The products are washable and reusable, manufactured to a high standard to give effective therapy, together with continuous lubrication of the tissue.

Each product is extremely soft and elastic and is comfortable to wear. All products are designed to be worn against the skin and fit easily in the shoe.



Silicone Heel Cup with Soft Heel Spot

- Silicone heel support with stabilising walls
- High viscosity/elasticity cushions loading points when walking/standing
- New design of silicone distributes load more evenly
- Softer silicone area in centre of heel
- Hard-wearing, unbreakable under normal foot loads
- Sold in pairs


Size	Code
S	L-405
M	L-406
L	L-407
XL	L-408



Silicone Heel Cup

- Absorbs shocks and vibrations, distributing the weight evenly over the whole foot
- Silicone heel support with stabilising walls
- High viscosity/elasticity cushions loading points when walking/standing
- New design of silicone distributes load more evenly
- Hard-wearing, unbreakable under normal foot loads
- Sold in pairs

Size	Code
S	L-401
M	L-402
L	L-403
XL	L-404




Silicone Gel Metatarsal Pad with Elastic Band

- Comfortable elastic band with large metatarsal silicone gel pad
- Guaranteed shock absorbing action on the entire plantar areas of the forefoot
- Useful for all types of metatarsal pain
- Sold in single units

Size	Code
S - M	L-473
L - XL	L-474

Contains Latex



Cap

- Elasticated protective toecap lined with silicone gel
- Used in the treatment of toe or nail lesions
- Sold singly

Shoe Size	Code
3 - 5	L-455
6 - 8	L-456
9 - 12	L-457




Elastic Sheath

- Seamless elasticated tubular protection lined with silicone gel on one side
- Can be cut to size
- Useful in the treatment of callouses, blisters, soft corns and hammer toes
- Sold singly

Shoe Size	Code
3 - 5	L-461
6 - 8	L-462
9 - 12	L-463

Contains Latex



Tubular Protection A

- Seamless elasticated tubular protection lined with silicone gel on one side
- Can be cut to size
- Useful in the treatment of callouses, blisters, soft corns and hammer toes
- Sold singly

Shoe Size	Code
3 - 5	L-449
6 - 8	L-450
9 - 12	L-451

Contains Latex




Tubular Protection B

- Seamless elasticate tubular protection
- Cut to length lined with silicone gel
- Useful in the treatment of callouses, blisters, soft corns and hammer toes
- Sold singly

Shoe Size	Code
3 - 5	L-443
6 - 8	L-444
9 - 12	L-445


Contains Latex



Toe Retractor

- Silicone gel pad avoids overlapping of toes
- Reduces pain caused by callouses between toes
- Sold singly

Shoe Size	Code
3 - 5	L-430
6 - 8	L-431
9 - 12	L-432



Silicone Gel Discs

- Fitted on adhesive fabric for multiple use
- Protects against rubbing or direct pressure on particular areas of the foot
- Sold in pairs

Shoe Size	Code
Universal	L-495



Silicone Gel Retractor / Protector

- Helps to realign the hallux valgus, preventing the first toe from overlapping the second
- Prevents rubbing between toes and reduces pain caused by callouses
- Sold singly

Shoe Size	Code
3 - 5	L-412
6 - 8	L-413
9 - 12	L-414



Silicone Gel Hallux Valgus Protector

- Pad design helps distribute load about joint
- Reduces and removes pain associated with bursitis of first MP joint
- Sold in single units

Size	Code
S - M	L-439
M - L	L-440



Anti-slip Pad

- Self-adhesive silicone gel pad
- Easily attached to inner heel of shoe
- Prevents the shoe from slipping off and protects the heel
- Sold in pairs

Size	Code
Universal	L-478

Spinal Support





MP45 Plus 3-Point Hyperextension Brace

Product Code: **ORT-**

Lightweight aluminium construction with unique spring-loaded 'floating' pelvic band. Thorocolumbar pad can be positioned horizontally or vertically. Unique 'Plus-lock' closure ensures controlled posterior loading.

Sternal component can be adjusted without bending, using serrated disc for angular adjustments in 5° settings.

Size	S	M	L	XL
Short	-	*/PMS	*/PLS	*/Pxls
Standard	*/PS	*/PM	*/PL	*/PXL
Long	-	*/PML	-	-

* Insert product code



Nhyco Cruciform Hyperextension Orthosis

Product Code: **NHYCO**

A strong, lightweight aluminium alloy frame brace, weighing less than 700g. The replaceable EVA padding is non toxic, waterproof and latex-free.

Features new 'Plus' safety lock closure system, easy for patient to apply/remove brace and ensures the same lumbar pressure each time.

Size	Code
Short	TO410CS
Regular	TO410CR
Long	TO410CL

* Insert product code



Chairback Brace

Product Code: **CBB**

A lightweight, spinal orthosis to maintain the spine in neutral and provide sacro-lumbar stability. Made in heat-mouldable Kydex for a custom fit and lined with closed cell foam, this brace is lighter and thinner than conventional metal and leather braces, and ideal for elderly and osteoporotic patients who cannot tolerate heavier braces.

Size	Code	Hip Circumference	Brace Height (Center back)
S	*/03	74 - 81 cm	29 cm
M	*/05	84 - 91 cm	32 cm
M-L	*/06	94 - 102 cm	33 cm
L	*/07	104 - 112 cm	34 cm

* Insert product code



NEW

Ortholux + Extender

Product Code: **ORT-TLS-**

Bledsoe's Ortholux spine orthoses feature the advanced engineering health care providers and patients alike have come to expect from Bledsoe brace products. With a sleek, breathable, low profile design Ortholux delivers superior patient comfort and better overall results. Straps with D-rings for easy donning, Extender Chest piece attached, and adjustable anterior height maximises brace utility.

Features include:

- Plush, breathable, moisture-wicking, tri-laminate foam
- Corrugated and padded posterior kydex inserts
- Straps with D-rings for easy donning
- Extender chest piece attached

Indications:

Chronic Back Pain, Compression Fracture, Degenerative Disc Disease, Herniated or Bulging Disc, Kyphosis, Ligament Strain/Sprain, Post-Laminectomy, Spinal Stenosis, Spondylolisthesis, Spondylolysis, Osteoporosis

Size	Code	Waist Circ (cm)	Hip Circ (cm)	Front Height (cm)	Back Height (cm)
S	*SM-EXT	70 - 85	70 - 90	24	45
M	*MD-EXT	75 - 90	80 - 95	26	46
L	*LG-EXT	85 - 100	90 - 105	28	47
XL	*XL-EXT	85 - 110	100 - 115	28	47
XXL	*XXL-EXT	105 - 125	110 - 130	28	52

* Insert product code



Corrugated and padded posterior kydex inserts



Spinal Support



Pro Tem Belt

Product Code: **TB**

Simple but effective, the high quality elastic and two flexible posterior stays make this our most popular belt.

The edges have quality Lycra binding and a back depth of 20cm cut away to 18cm at the front for increased comfort when sitting.

Size	Code	NDC Code	Waist Measurement
S	*/03	-	66 - 74 cm
M	*/05	-	76 - 84 cm
L	*/07	052424	86 - 94 cm
XL	*/08	052431	96 - 104 cm
XXL	*/09	-	107 - 114 cm

* Insert product code





Proform Belt

Product Code: **PFM**

An elastic front with hook/loop closure and a cotton back section, with detachable comfort pad and two malleable steels which can be shaped to fit the spinal curve.

A webbing strap which can be positioned in the lumbar or sacral area, for extra compression on the spinal curve. Back depth is 28cm to support the thoracic area reducing to 20cm at the front for comfort.

Contains Latex

Size	Code	NDC Code	Waist Measurement
S	*/83	-	58 - 66 cm
M	*/85	-	68 - 76 cm
L	*/87	-	79 - 86 cm
XL	*/88	020522	89 - 97 cm
XXL	*/89	-	99 - 107 cm

* Insert product code



Crossform Belt

Product Code: **CFB**

Two layers of bias-cut elastic give excellent shaping around the waist and hips, also features with a padded cotton front, 17cm deep with hook/loop closure to give mild abdominal pressure.

The 20cm deep back pocket encloses a thermoplastic panel which is moulded in hot water at 65° and shaped to the spine.

Contains Latex

Size	Code	Waist Measurement
S	*/83	64 - 71 cm
M	*/85	74 - 81 cm
L	*/87	84 - 91 cm
XL	*/88	94 - 102 cm

* Insert product code



Sacro Lumbar Support

Product Code: **SLS**

To accommodate the waist/hip differential, two layers of elastic are sewn at an angle giving this belt increased strength and good conformity. The back depth is 23cm, front depth is 18cm.

The posterior pocket is lined with Neoprene for extra comfort and warmth and a flexible panel (FPP/8) or a thermoplastic panel (TPP/8) can be supplied separately if required.

Contains Latex

Size	Code	Waist Measurement
S	*/03	64 - 74 cm
M	*/05	76 - 84 cm
L	*/07	86 - 94 cm
XL	*/08	96 - 104 cm

* Insert product code



TLSO Back Brace

The TLSO Back Brace provides full-circumferential rigid support from the T2/T3 scapular thoracic region to the sacrococcygeal junction. It features rigid anterior and posterior panels, a telescoping sternal aluminum Y-Bar and dual lacer compression pads.

- Easy application with dual side pulls
- Even, progressive compression from proprietary teardrop Lacer Compression System
- Easily adjustable for all-day wear
- Comprised of a breathable Airmesh liner for maximum comfort
- Midline opening in posterior panel eases pressure over the incision site, allowing for wound inspection without removing the orthosis
- Heat moldable panels

Code	Size	Circumference Waist
60002	S	25" - 30"
60003	M	30" - 35"
60004	L	35" - 42"
60005	XL	42" - 48"
60006	XXL	48" - 57"



BOA with Chairback Back Brace

The BOA Back Brace with Chairback provides exceptional support from T9 - S1. The modular construction of this brace is designed to support a three-stage treatment protocol.

- Easy application with single-pull design
- Even, progressive compression from proprietary teardrop Lacer Compression System
- Easily adjustable for all day wear
- Breathable Airmesh fabric for maximum comfort
- Midline opening eases pressure over the incision site, allowing for wound inspection without removing the orthosis
- Heat moldable chairback panel for a customized fit
- Lightweight

Code	Description
6002X	BOA w/ Chairback 8"
6001X	BOA w/ Chairback 10"
6013X	BOA Chairback Liner Available in size S-XL only
60128	Anterior Extension Panel (16" L x 7" H)
6008X	Chairback Accessory, Rigid Frame, S - XL (13" Tall)
6009X	Chairback Accessory, Short Rigid Frame, M - XL (13" Tall)

X = Insert size # from table below

Size #	Size	Circumference Waist
X = 2	S	25" - 30"
X = 3	M	30" - 35"
X = 4	L	35" - 40"
X = 5	XL	40" - 45"
X = 6	XXL	45" - 50"
X = 7	XXXL	50" - 55"





The Proflex Belt

Features Include:

- Hook/loop fastening elastic belt, with an elastic back pocket to accommodate a panel.
- Available with a choice of back panels.
- A cotton fulcrum strap which can be positioned in the lumbar or sacral area from extra compression on the spine.
- The front depth is 20cm.



Proflex Belt (No panel)

Product Code: **PFN/8** (20cm back depth)

Size	Code	Waist Measurement
S	*3	64 - 76 cm
M	*5	79 - 89 cm
L	*7	89 - 99 cm
XL	*8	99 - 109 cm
XXL	*9	109 - 119 cm

* Insert product code

Contains Latex

Proflex Belt (Flexible Panel)

Product Code: **PF/8** (20cm back depth)

The back pocket encloses a panel of polythene lined with closed cell foam, which becomes flexible with body heat.

Size	Code	Waist Measurement
S	*3	64 - 76 cm
M	*5	79 - 89 cm
L	*7	89 - 99 cm
XL	*8	99 - 109 cm
XXL	*9	109 - 119 cm

* Insert product code

Contains Latex



The Lumbarform Belt

Features Include:

- Triple elastic fulcrum straps at the back which attach to one simple strap at the front, giving extra compression on the spinal curve.
- Two flexible stays which keep the belts smooth at the front for comfort when sitting.
- Ventilated elastic which has Lycra binding to prevent chafing.
- The Lumbarform belt range is available with a choice of back panels and stays.

Standard Lumbarform Belt (Flexible Panel)

Product Code: **LBF/9** (23cm back depth)

The double-layer back pocket encloses a panel of polythene lined with closed cell foam, which becomes flexible with body heat.

Size	Code	Waist Measurement
S	*3	69 - 76 cm
M	*5	79 - 86 cm
L	*7	89 - 97 cm
XL	*8	99 - 107 cm

* Insert product code

Contains Latex

Lumbarform Belt with Stays

Product Code: **LBS/9** (23cm back depth)

Two flexible posterior stays and a double-layer ventilated elastic back pocket stabilise the spine.

Size	Code	Waist Measurement
S	*3	69 - 76 cm
M	*5	79 - 86 cm
L	*7	89 - 97 cm
XL	*8	99 - 107 cm

* Insert product code

Contains Latex

Deep Lumbarform Belt (Flexible Panel)

Product Code: **LBF/12** (31cm back depth)

The double-layer back pocket encloses a panel of polythene lined with closed cell foam, which becomes flexible with body heat.

Size	Code	Waist Measurement
S	*3	69 - 76 cm
M	*5	79 - 86 cm
L	*7	89 - 97 cm
XL	*8	99 - 107 cm

* Insert product code

Contains Latex

A close-up, profile view of a pregnant woman with long dark hair, wearing a blue long-sleeved sweater. She is gently holding her pregnant belly with both hands. Her eyes are closed, and she has a soft smile. The background is a soft, out-of-focus white.

Maternity Support



- Adaptable to a standing or sitting position and can be worn over clothing. Provides strong circumferential support around the pelvis holding the two sacroiliac joints and the symphysis pubis joint firmly.

- Also available in black



 Scotland
NDC Codes



 Scotland
NDC Codes

Pelvic Instability Belt

Product Code: **SPD**

- Designed by obstetric physiotherapists and trialed by expectant mothers, this 'many-tailed bandage' pattern allows variable tension on each set of straps, with strong tension over the pubic area and lighter tension under the abdomen.

Contains Latex

Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*/1	*/1/BK	90 - 96 cm
*/2	*/2/BK	96 - 102 cm
*/3	*/3/BK	102 - 108 cm
*/4	*/4/BK	108 - 116 cm
*/5	*/5/BK	116 - 126 cm
*/6	*/6/BK	126 - 136 cm
*/7	*/7/BK	136 - 146 cm
*/8	*/8/BK	146 - 156 cm

* Insert product code

Wiltshire Maternity Belt

Product Code: **WMB**

- A pelvic support belt with an elastic back section and two padded side straps, plus two variable tension straps.
- The slimline design of this belt makes it comfortable to wear when sitting.
- Provides strong circumferential support around the pelvis holding the two sacroiliac joints and the symphysis pubis joint firmly.
- Back depth 10cm, front depth 8cm.
- Also available in black.

Contains Latex

Code (White)	NDC Code (White)	Code (Black)	Circumference
*/1	205707	*/1/BK	96 - 116 cm
*/2	205714	*/2/BK	116 - 130 cm
*/3	-	*/3/BK	130 - 146 cm
*/4	-	*/4/BK	146 - 160 cm

* Insert product code

Trochanteric Belt

Product Code: **TROCH**

- A strong webbing strap with foam padding gives firm circumferential support at trochanteric level to ease the pain of separated symphysis.
- Can also be used to treat sport-induced symphysis pubis pain.
- Also available in black.

Size	Code (Navy Blue)	NDC Code (Navy Blue)	Code (Black)
One size	*/10	052448	*/10/BK

* Insert product code



Maternity Belt Style A

Product Code: **MBAS**

- Made in quality elastic with oblique-angled hook/loop closure for circumferential support around the pelvis. Depth 20cm back and front. Also available in black.

Contains Latex

Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*/1	*/1/BK	90 - 96 cm
*/2	*/2/BK	96 - 102 cm
*/3	*/3/BK	102 - 108 cm
*/4	*/4/BK	108 - 116 cm
*/5	*/5/BK	116 - 126 cm
*/6	*/6/BK	126 - 136 cm
*/7	*/7/BK	136 - 146 cm
*/8	*/8/BK	146 - 156 cm
*/9	*/9/BK	156 - 165 cm
*/10	*/10/BK	165 - 178 cm

* Insert product code



- Made in stronger lenoweave elastic/cotton mix for extra support.
- Back depth 20cm, front depth 17cm.
- Also available in black.

Contains Latex

Maternity Belt Style B

Product Code: **MBB**

- Cut away design at the front gives comfort when sitting and prevents belt from rolling up.
- Angled posterior seams give improved shaping around the sacrum.

Code (White)	NDC Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*/1	-	*/1/BK	90 - 96 cm
*/2	-	*/2/BK	96 - 102 cm
*/3	002214	*/3/BK	102 - 108 cm
*/4	002221	*/4/BK	108 - 116 cm
*/5	002238	*/5/BK	116 - 126 cm
*/6	002245	*/6/BK	126 - 136 cm
*/7	-	*/7/BK	136 - 146 cm
*/8	-	*/8/BK	146 - 156 cm
*/9	-	*/9/BK	156 - 165 cm
*/10	-	*/10/BK	165 - 178 cm
*/11	-	*/11/BK	178 - 185 cm
*/12	-	*/12/BK	185 - 196 cm

* Insert product code



- Made in stronger lenoweave elastic/cotton mix for extra support.
- Back depth 16cm, front depth 11cm.
- Also available in black.

Contains Latex

Maternity Belt Style C

Product Code: **MBC**

- The low-cut design of this belt makes it comfortable to wear when sitting.
- Angled posterior seams give improved shaping around the sacrum.

Code (White)	NDC Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*/1	-	*/1/BK	90 - 96 cm
*/2	076796	*/2/BK	96 - 102 cm
*/3	076802	*/3/BK	102 - 108 cm
*/4	076819	*/4/BK	108 - 116 cm
*/5	076826	*/5/BK	116 - 126 cm
*/6	076833	*/6/BK	126 - 136 cm
*/7	076840	*/7/BK	136 - 146 cm
*/8	-	*/8/BK	146 - 156 cm
*/9	-	*/9/BK	156 - 165 cm

* Insert product code



Contains Latex

Deluxe Maternity Belt

Product Code: **DMB**

- Anterior foam laminate panels - low profile shape for improved fit, especially when sitting.
- Hook/loop closure.
- Adjustable stretch fulcrums give additional support.
- Helps relieve lower back pain and associated symptoms.

Code	Circumference (At widest part of abdomen)
*/1	90 - 96 cm
*/2	96 - 102 cm
*/3	102 - 108 cm
*/4	108 - 116 cm
*/5	116 - 126 cm
*/6	126 - 136 cm
*/7	136 - 146 cm
*/8	146 - 156 cm
*/9	156 - 165 cm
*/10	165 - 178 cm

* Insert product code



Duofit Maternity Belt

Product Code: **DUB**

- 20cm deep circumferential elastic wrap with anterior cut-away to improve comfort when sitting.
- Angled hook/loop closure to allow accurate positioning under tummy.

Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*/1	*/1/BK	90 - 96 cm
*/2	*/2/BK	96 - 102 cm
*/3	*/3/BK	102 - 108 cm
*/4	*/4/BK	108 - 116 cm
*/5	*/5/BK	116 - 126 cm
*/6	*/6/BK	126 - 136 cm
*/7	*/7/BK	136 - 146 cm
*/8	*/8/BK	146 - 156 cm
*/9	*/9/BK	156 - 165 cm
*/10	*/10/BK	165 - 178 cm

* Insert product code

- Additional narrow pelvic wrap provides extra support around the pelvis for improved stability at SI joints.
- Also available in black.



Lo-pro Maternity Belt

Product Code: **LPM**

- A low profile belt manufactured using a newly designed high strength, super thin, breathable elastic.
- Anterior shaping makes this belt unobtrusive and extremely comfortable to wear whilst standing or sitting.

Code	Circumference (At widest part of abdomen)
*/1/BK	90 - 96 cm
*/2/BK	96 - 102 cm
*/3/BK	102 - 108 cm
*/4/BK	108 - 116 cm
*/5/BK	116 - 126 cm
*/6/BK	126 - 136 cm
*/7/BK	136 - 146 cm
*/8/BK	146 - 156 cm
*/9/BK	156 - 165 cm

* Insert product code

- The Lo-pro Maternity Belt provides excellent support, with improved patient compliance.
- Only available in black.

Contains Latex



Fracture Bracing & Splinting



Standard Humeral Fracture Brace

Product Codes: **SHC** (Full Kit), **SHS** (Sleeve Only)

A polypropylene humeral sleeve with closed cell foam lining, and a polythene insert for comfort and improved fit. The sleeve can be cut with bandage shears and is heat mouldable. Kit includes two humeral cast socks and a fully adjustable sling.

Size	Code (Full Kit)	Code (Sleeve Only)	Length	Circumference of upper arm	
				Distal	Proximal
S	SHC/03	SHS/03	15 cm	18 - 23 cm	23 - 28 cm
M	SHC/05	SHS/05	17 cm	23 - 28 cm	28 - 33 cm
L	SHC/07	SHS/07	18 cm	28 - 33 cm	33 - 38 cm
XL	SHC/08	SHS/08	23 cm	33 - 38 cm	38 - 43 cm



Extended Humeral Fracture Brace

Product Codes: **EHC** (Full Kit), **EHS** (Sleeve Only)

An extended polypropylene humeral sleeve lined with closed cell foam, and a polythene insert to treat non-union of high humeral fractures. The sleeve can be cut with bandage shears and is heat mouldable.

Kit includes two humeral cast socks and a fully adjustable sling.

Size	Code (Full Kit)	Code (Sleeve Only)	Length	Circumference of upper arm	
				Distal	Proximal
S	EHC/03	EHS/03	15 cm	18 - 23 cm	23 - 28 cm
M	EHC/05	EHS/05	17 cm	23 - 28 cm	28 - 33 cm
L	EHC/07	EHS/07	18 cm	28 - 33 cm	33 - 38 cm
XL	EHC/08	EHS/08	23 cm	33 - 38 cm	38 - 43 cm

Sling strap

Product Code: **SS**

Sling Strap is a non-stretch padded 50mm wide material that is combined with Hook Y-Tabs to provide an adjustable and versatile sling system.

Sling Strap can be cut to the required length to produce a simple collar n' cuff or similar type of sling. This reduces waste and the adjustability of the material ensures that the sling will fit and can be easily adjusted. Sling Strap is padded for comfort around the neck and the modern black colour improves patient compliance and looks better for longer. Sling Strap is supplied in a pack that contains: 2 x 12m rolls & 20 Y-tab closures.



Length	Width
2 x 12 Metres	5 cm



Humeral Stabilising Orthosis

Product Codes: *HSO/10/K (Full Kit)* , *HSO/10 (Sleeve Only)*

A unique lightweight, comfortable humeral sleeve with a ventilated 'skin friendly' protective lining, which means there is no humeral sock required.

The universal humeral is fully customisable, the thermo formable external sleeve can be altered using standard scissors (A sizing guide is included with sleeve).

An optional kit is available which ensures immobilisation of the arm, eliminating swinging of the arm. Kit includes a Deluxe Shoulder Immobiliser, product code: DSI (See pg. 33).



Code	Description
HSO/10/K	Full Kit
HSO/10	Sleeve Only
HSO/10/E	Sleeve Only (Extended)



Soft Humeral

Product Code: *PSH*

A lightweight preformed brace made from a tri-lamination with outer fabric facing, semi-flexible foam and inner smart-textile lining.

- Inner lining has Coolmax treatment to reduce perspiration and Ion/silver treatment to improve temperature control within the brace.
- The kit is supplied with lightweight humeral sock and fully adjustable hook/loop strap closures.
- If required, the brace can be trimmed.

Size	Code	Length	Proximal Circ	Distal Circ
XS	*/02	28 cm	20 - 23 cm	15 - 19 cm
S	*/03	30 cm	23 - 28 cm	19 - 23 cm
M	*/05	32 cm	28 - 33 cm	23 - 28 cm
L	*/07	34 cm	33 - 38 cm	28 - 33 cm
XL	*/08	35 cm	38 - 43 cm	33 - 38 cm

* Insert product code



Promedics Hip Joint

Product Code: *PHJ/10*


- The Promedics hip hinge allows control of abduction (between -5 to +25°) with variable flexion and extension stops.
- The abduction axis may be set at various positions while the brace is on the patient by simply turning the central adjustment screw with the allen key provided.
- The flexion/extension control can be altered by positioning the screws in the relevant marked stops to give limited range of motion.
- A pre-attached plastic headplate is fitted to the hip hinge for use in conjunction with casting materials and thermoplastics.
- Length 40cm.
- Waist and shoulder straps to be ordered separately.

Universal Hip Abduction Waist Belt

Product Code: *JBB/10*

Universal Hip Abduction Waist Belt with Shoulder Strap

Product Code: *JBBS/10*



- A geared, polycentric hinge which can be set to limit the range of flexion and extension at the knee.
- Available in offsets of 12mm, 18mm and 32mm, with strong aluminium side steels which can be shaped with a bending iron.
- Detachable headplates allow easy removal of the hinges from the cast brace.
- Order hinges in pairs.



Alignment Jig for Genurange Hinges

Product Code: **AJ/10**

To assist with correct positioning of the Genurange hinges for cast bracing.



Genurange Fracture Bracing Hinges

Product Code: **GNF**

Order Genurange Range of Motion Hinges in pairs.

Size	Medial Offset	Lateral Offset
*/12/R	12 mm	12 mm
*/12/L	12 mm	12 mm
Size	Medial Offset	Lateral Offset
*/18/R	18 mm	12 mm
*/18/L	12 mm	18 mm
Size	Medial Offset	Lateral Offset
*/32/R	32 mm	12 mm
*/32/L	12 mm	32 mm

* Insert product code



Replacement Headplates

Product Code: **GNP/10**

Replacement headplates for Genurange hinges. Detachable for easy removal of the cast brace; supplied in sets of 4.

Replacement Screws

Product Code: **GNS/10**

Replacement screws to fix headplates on the Genurange hinge; supplied in sets of 4.

Replacement Hinge Cover Screws

Product Code: **GEN/CS**

Replacement screws for Genurange hinge cover.



Derby Heel Cup

Product Codes: **DHC**

Made in polypropylene plastic with straight or butterfly stems, to fix the heel cup in plaster or synthetic bandage and prevent cast migration.

The hinged ankle gives ankle stability, allowing plantar/dorsiflexion but limiting medial/lateral movement.

Universal design fits left and right.

Size	Code	To fit shoe size
Universal/Standard (Butterfly Stems)	*/1	3 - 6
Universal/Large (Butterfly Stems)	*/2	6 - 9
Universal/Standard (Straight Stems)	*/3	3 - 6
Universal/Large (Straight Stems)	*/4	6 - 9

* Insert product code



Femoral Sock

Product Code: **FEM**

- Full length femoral sock
- Acrylic/lycra fabric
- Range of sizes for individual fitting
- Available singly
- Wicking properties help reduce perspiration

Size	Code	Dimensions
1	*/1	91 x 15 cm
2	*/2	89 x 23 cm
3	*/3	99 x 23 cm
4	*/4	107 x 23 cm
5	*/5	117 x 23 cm

* Insert product code

Humeral Sock

Product Code: **HS/10**

- Acrylic/lycra fabric
- Universal fitting
- Available in pairs only
- Wicking properties help reduce perspiration

Size
Universal

Tibial Sock

Product Code: **BK/10**

- Below knee tibial sock
- Acrylic/lycra fabric
- Universal fitting
- Available singly
- Wicking properties help reduce perspiration

Size
Universal

neofrakt®

The rigid foam splint system

The principle of the NEOFRAKT® system is the combination of polyurethane rigid foam as a splint, and the elastic stocking as bandage surface. The polyurethane rigid foam NEOFRAKT® employs was developed especially for medical purposes. The NEOFRAKT® bandage is applied in soft condition. It hardens on the body and is easy to modify.

NEOFRAKT® can not only be used for simple splints and cylinder bandages. It is also very suitable for complicated temporary orthopaedic applications.

The polyurethane hard and soft foams were especially developed for medical applications. These polyurethanes harden by mixing 2 components which are packed in 2 separate cans. These cans are connected to each other by a plastic ring and they are not sold separately.



- Removable through zip(s) and/or hook / loop fastener(s)
- Self moulding - perfect fit
- Easy to apply
- X-ray transparent
- Low weight
- Water resistant
- Durable & Comfortable, 100%



Removable full leg walker



Removable corset



Lower arm with thumb hole



Thorax corset with forehead support indicators

*For further details on any Neofrakt product,
please contact Customer Services:*

Call: 01475 746 400 or E-mail: enquiries@promedics.co.uk



Loop Strap

Product Code: **LS**

Loop strap is a hardwearing strap that will accommodate multiple hook closures. Loop Strap is hook receptive on both sides and can be cut to the required length. Loop Strap really is a hardwearing long lasting strap solution that will have multiple applications in the splinting workshop. Loop Strap will last and last!

Width	Length	Code
20 cm	10 Meter roll	LS/20
25 cm	10 Meter roll	LS/25
30 cm	10 Meter roll	LS/30
40 cm	10 Meter roll	LS/40
50 cm	10 Meter roll	LS/50



Terry Strap

Product Code: **TS**

Terry Strap is a padded, adhesive backed, terry towel lined material that is available in three widths. 3mm thick Terry Strap is a versatile material that is effective when used to line splints or to pad the edges of casts or splints. The structure of the material allows for a limited amount of stretching and shaping to help when forming it to the edges or inner surface of a splint or cast. This is a really useful material to have in the plaster room and the splinting workshop and it will greatly improve the comfort and finish of any splint or cast

Width	Length	Code
3 cm	15 Metre roll	TS/30
5 cm	15 Metre roll	TS/50
8 cm	15 Metre roll	TS/80



Goniometer

Product Code: **GON/10**

Perfect for precision finger measurements and great for measuring difficult finger angles.



Alphastrap

Product Code: **APS**

- A hook / loop attractive laminate material with limited stretch, which gives soft cushioned support and holds splints in place
- Edges self-seal when cut

Size	Size	Colour
APS/1	1" x 60 ft	Beige
APS/2	2" x 30 ft	Beige
APS/3	1.5" x 40 ft	Beige
APS/4	1" x 60 ft	Blue



Proform Plus

Resistance to Stretch

- Minimum to moderate resistance to stretch allows the splintmaker to control the material throughout the splinting process.
- Proform Plus tolerates minimum to moderate handling during splint moulding.
- Handle on a horizontal plane and form splint in a gravity-assisted position.
- Gentle, smooth motions yield superior results.
- Proform Plus has moderate drape and offers conformability with only light to moderate pressure. Light weight, 3/32" Proform Plus provides more drape than 1/8" Proform Plus.

Memory/Recovery

- With moderate memory and recovery, splints can be repeatedly remoulded without thinning.
- Reheated splints will not thin or stretch, making it easy to adjust for bandage or edema changes.
- Splints made from Proform Plus have a moderate degree of rigidity.
- A circumferential splint made of Proform Plus is flexible enough for easy donning and doffing, yet strong enough for immobilisation.
- Splints will give slightly against strong force and will resist breaking.

Bonding

- The unique coating gives Proform Plus a light tack to hold the material in place, allowing the splintmaker to achieve contours and positioning without the material pulling or falling away.
- A wet bond will hold firmly but can be pulled apart in warm water if necessary.
- Dry heating two surfaces will form a strong bond. For a permanent bond, solvent may be required.

Surface Finish

- Low-profile, beige colour has moderate resistance to fingerprints and resists showing stains.
- Smooth, self sealing edges do not require additional finishing, resulting in a professional finish, even after reheating.
- Working time is sufficient enough to mould the splint until it holds its shape, remove it to do necessary trimming and then reposition it to complete the moulding process.

Applications

- 3.2mm Proform Plus can be used for almost any splinting application. Use for dynamic and static splinting; for adult and paediatric positioning; and for post-burn and post-surgical care.
- 2.4mm Proform Plus is thin and lightweight, yet strong, and can be used for forearm and hand-based splints. Ideal for arthritis splints, post-operative splints, and for paediatric splints for children who are able to co-operate.

Code	Description	Size
PPF-1S2	3.2 mm Proform Plus Smooth	46 x 61 cm
PPF-1S1	3.2 mm Proform Plus Smooth	61 x 91 cm
PPF-1P2	3.2 mm Proform Plus Perforated	46 x 61 cm
PPF-1P1	3.2 mm Proform Plus Perforated	61 x 91 cm



Leader in
thermoplastic
innovations

Based on innovative technologies, Orfit develops and produces the most precise and reliable thermoplastic materials for medical devices that improve patient treatment around the world.

Orthotic fabrication materials from Orfit offer the highest possible performance to both therapists and patients in physical rehabilitation.

Our wide product range consists of engineered thermoplastics that will help you solving your most challenging splinting requirements. Through workshops and courses, our experts provide lots of tips and tricks to make your orthotic fabrications functional and successful.

All of our thermoplastic materials are developed and produced to make sure that they are perfectly suited to their purpose and your needs.



Orflight (NS)

- Immobilization systems for cancer patients in Radiation Oncology
- Orthotic fabrication materials for patients in Physical Rehabilitation
- Prosthetic socket materials for amputee patients

www.orfit.com



Orficast



Orfit Classic



Aquafit NS



Orfit Colours NS



Orfit Eco (NS)



Dynamics - Dynamic Splinting

For further details on any Orfit product,
please contact Customer Services:

Call: **01475 746 400** or E-mail: enquiries@promedics.co.uk



Standard Hook / Loop

Product Code: **HAL**

- 25m reels of plain nylon hook and loop materials
- Available in White, Black, Navy Blue, Royal Blue, Beige, Red, Dayglo Pink, Dayglo Yellow, Olive Green
- A stock matrix is available upon request
- When ordering please specify colour and width, hook or loop

Code (Hook)	Code (Loop)	Width
*/15/H	*/15/L	15 mm
*/20/H	*/20/L	20 mm
*/25/H	*/25/L	25 mm
*/30/H	*/30/L	30 mm
*/50/H	*/50/L	50 mm

* Insert product code



Self Adhesive Hook / Loop

Product Code: **SA**

- 25m reels of adhesive backed nylon hook and loop materials
- Available in White, Black, Navy Blue, Royal Blue, Red, Dayglo Pink, Dayglo Yellow
- A stock matrix is available upon request
- When ordering please specify colour and width, hook or loop

Code (Hook)	Code (Loop)	Width
*H/15	*L/15	15 mm
*H/20	*L/20	20 mm
*H/25	*L/25	25 mm
*H/30	*L/30	30 mm
*H/50	*L/50	50 mm

* Insert product code

Hook / Loop
Colours



Black



White



Navy Blue



Beige



Red



Olive Green



Royal Blue



Dayglo Yellow



Dayglo Pink



Juraprene Splinting Material

Juraprene is a welcome addition to Promedics splinting range. A latex free neoprene material which can be used to create custom splints. Juraprene offers creative splinting options while providing warmth, compression, comfort and support.

Features of Juraprene include:

- 3mm thickness to provide light, restrictive support and compression as required.
- Plush lined neoprene for improved comfort and to reduce pressure points.
- Latex Free
- Machine Washable
- Sheet size 45 x 61 cm

Code	Description
JPR/BEI	Beige
JPR/PUR	Purple
JPR/RBL	Royal Blue
JPR/RED	Red
JPR/BLK	Black



Cold Therapy





Polar Care Kodiak® Cold Therapy

Product Code: **100291**

The Polar Care Kodiak Cold Therapy System is the most convenient and versatile offering in Breg's Polar Care line. Its easy to use, compact design makes it great for clinic, hospital, and home use. With the addition of a little ice and water, patients will enjoy 6-8 hours of effortless cold therapy.

Offering the only battery powered option on the market, each battery pack comes with four replaceable AA batteries ready to power 10-14 hours of motorized cold therapy. Proper use requires an insulation barrier between the Intelli-Flo® pad and the skin.

Common Examples of Use

- Post Operative
- Arthroscopic Procedures
- Reconstructive Procedures
- Plastic Surgery
- General Surgery
- Post-Trauma
- Chronic Pain
- Physical Therapy

Product Features

- Proven reliability and effectiveness
- Easy to use, intuitive design
- Customized Intelli-Flo pads provide temperature control specific to each treatment area resulting in superior coverage and comfort
- Provides motorized cold treatment for 6-8 hours



Battery pack

Code	Accessory
100316-000	6V Power Supply
97050	Battery Pack

Code	Description
100291-000	Polar Care Kodiak® Cold Therapy Unit

Code	IntelliFlo Pads
100291-000	Shoulder
100293-000	Shoulder, X-Large
100294-000	Knee
100295-000	Knee, Compression
100296-000	Multi-use
100297-000	Multi-use, Long stem
100298-000	Ankle
100392-000	Ankle, Compression
100300-000	Back
100301-000	3 x 5 (8 cm x 13 cm)
100302-000	Hip
100303-000	Hand/Wrist



Universal Cold Therapy Wraps

Product Codes: **CTKW/10** (Knee) **CTSW/10** (Shoulder)

The Promedics Cold Therapy Wraps feature two reusable 12 x 24 cm Hot/Cold packs that slip into the internal pouches. Hook and loop secures hot/cold packs to inside of wrap. Latex Free, Washable, and Reusable. Universal size, fits left or right.

Code	Description	Size
CTKW/10	Knee Wrap	Universal
CTSW/10	Shoulder Wrap	Universal





Paediatric



Lycrafleece™

Paediatric Thumb Spica

Product Code: **PLTS**

- Manufactured in Lycrafleece material
- Colourful hook / loop design
- Flexible splint to support thumb and provide greater comfort with thumb movement

Size	Code (Left)	Code (Right)	Wrist Circumference
XS	*/12	*/22	12 - 13 cm
S	*/13	*/23	13 - 14 cm
M	*/15	*/25	14 - 15 cm
L	*/17	*/27	15 - 17 cm
XL	*/18	*/28	17 - 18 cm

* Insert product code



Lycrafleece™

Paediatric Wrist Splint

Product Code: **PLWB**

- Manufactured in Lycrafleece material
- Colourful hook / loop design
- 15 cm depth
- Sold in single units

Size	Code (Left)	Code (Right)	Wrist Circumference
XS	*/12	*/22	12 - 13 cm
S	*/13	*/23	13 - 14 cm
M	*/15	*/25	14 - 15 cm
L	*/17	*/27	15 - 17 cm
XL	*/18	*/28	17 - 18 cm

* Insert product code



Lycrafleece™

Paediatric Wrist/Thumb Splint

Product Code: **PLWTB**

- Manufactured in Lycrafleece material
- Colourful hook / loop design
- 15 cm depth
- Sold in single units

Size	Code (Left)	Code (Right)	Wrist Circumference
XS	*/12	*/22	12 - 13 cm
S	*/13	*/23	13 - 14 cm
M	*/15	*/25	14 - 15 cm
L	*/17	*/27	15 - 17 cm
XL	*/18	*/28	17 - 18 cm

* Insert product code



Contains Latex

Paediatric Elastic Wrist Brace

Product Code: **PEW**

Manufactured in single compression elastic. Contrasting colourful bar casing and edge binding. 15cm length. Sold in single units.

Size	Code (Left)	Code (Right)	Wrist Circumference
1	*/13	*/23	11.5 cm
2	*/15	*/25	14.5 cm
3	*/17	*/27	16 cm

* Insert product code



Contains Latex

Paediatric Elastic Wrist/Thumb Brace

Product Code: **PEWT**

Manufactured in single compression blue elastic. Contrasting colourful bar casing and edge binding. 15cm length. Sold in single units.

Size	Code (Left)	Code (Right)	Wrist Circumference
1	*/13	*/23	11.5 cm
2	*/15	*/25	14.5 cm
3	*/17	*/27	16 cm

* Insert product code



Paediatric Neoprene Wrist Brace

Product Code: **PNW**

- Manufactured in a range of coloured Neoprene.
- Contrasting colourful bar casing and hook/loop fastening
- 15 cm length
- Sold in single units

 Aqua = AQ	 Black = BK
 Navy = NY	 Pink = PK
 Purple = PU	 Royal Blue = RB
 Red = RD	

Size	Colour	Code (Left)	Code (Right)	Wrist Circumference
XXS	Aqua	PNW/11/AQ	PNW/21/AQ	9 cm
XS	Aqua	PNW/12/AQ	PNW/22/AQ	10.5
S	Aqua	PNW/13/AQ	PNW/23/AQ	11.5 cm
M	Aqua	PNW/15/AQ	PNW/25/AQ	14.5 cm
L	Aqua	PNW/17/AQ	PNW/27/AQ	16 cm
XXS	Black	PNW/11/BK	PNW/21/BK	9 cm
XS	Black	PNW/12/BK	PNW/22/BK	10.5
S	Black	PNW/13/BK	PNW/23/BK	11.5 cm
M	Black	PNW/15/BK	PNW/25/BK	14.5 cm
L	Black	PNW/17/BK	PNW/27/BK	16 cm
XXS	Navy	PNW/11/NY	PNW/21/NY	9 cm
XS	Navy	PNW/12/NY	PNW/22/NY	10.5
S	Navy	PNW/13/NY	PNW/23/NY	11.5 cm
M	Navy	PNW/15/NY	PNW/25/NY	14.5 cm
L	Navy	PNW/17/NY	PNW/27/NY	16 cm
XXS	Pink	PNW/11/PK	PNW/21/PK	9 cm
XS	Pink	PNW/12/PK	PNW/22/PK	10.5
S	Pink	PNW/13/PK	PNW/23/PK	11.5 cm
M	Pink	PNW/15/PK	PNW/25/PK	14.5 cm
L	Pink	PNW/17/PK	PNW/27/PK	16 cm
XXS	Purple	PNW/11/PU	PNW/21/PU	9 cm
XS	Purple	PNW/12/PU	PNW/22/PU	10.5
S	Purple	PNW/13/PU	PNW/23/PU	11.5 cm
M	Purple	PNW/15/PU	PNW/25/PU	14.5 cm
L	Purple	PNW/17/PU	PNW/27/PU	16 cm
XXS	Royal Blue	PNW/11/RB	PNW/21/RB	9 cm
XS	Royal Blue	PNW/12/RB	PNW/22/RB	10.5
S	Royal Blue	PNW/13/RB	PNW/23/RB	11.5 cm
M	Royal Blue	PNW/15/RB	PNW/25/RB	14.5 cm
L	Royal Blue	PNW/17/RB	PNW/27/RB	16 cm
XXS	Red	PNW/11/RD	PNW/21/RD	9 cm
XS	Red	PNW/12/RD	PNW/22/RD	10.5
S	Red	PNW/13/RD	PNW/23/RD	11.5 cm
M	Red	PNW/15/RD	PNW/25/RD	14.5 cm
L	Red	PNW/17/RD	PNW/27/RD	16 cm



Neoprene Wrist/ Thumb Brace

Product Code: **PNWT**

- Manufactured in a range of coloured Neoprene.
- Contrasting colourful bar casing and hook / loop fastening
- 15 cm length
- Sold in single units

Aqua = **AQ** Black = **BK** Navy = **NY** Pink = **PK** Purple = **PU** Royal Blue = **RB** Red = **RD**

Size	Colour	Code (Left)	Code (Right)	Wrist Circumference
1	Aqua	PNWT/13/AQ	PNWT/23/AQ	11.5 cm
2	Aqua	PNWT/15/AQ	PNWT/25/AQ	14.5 cm
3	Aqua	PNWT/17/AQ	PNWT/27/AQ	16 cm
1	Black	PNWT/13/BK	PNWT/23/BK	11.5 cm
2	Black	PNWT/15/BK	PNWT/25/BK	14.5 cm
3	Black	PNWT/17/BK	PNWT/27/BK	16 cm
1	Navy	PNWT/13/NY	PNWT/23/NY	11.5 cm
2	Navy	PNWT/15/NY	PNWT/25/NY	14.5 cm
3	Navy	PNWT/17/NY	PNWT/27/NY	16 cm
1	Pink	PNWT/13/PK	PNWT/23/PK	11.5 cm
2	Pink	PNWT/15/PK	PNWT/25/PK	14.5 cm
3	Pink	PNWT/17/PK	PNWT/27/PK	16 cm
1	Purple	PNWT/13/PU	PNWT/23/PU	11.5 cm
2	Purple	PNWT/15/PU	PNWT/25/PU	14.5 cm
3	Purple	PNWT/17/PU	PNWT/27/PU	16 cm
1	Royal Blue	PNWT/13/RB	PNWT/23/RB	11.5 cm
2	Royal Blue	PNWT/15/RB	PNWT/25/RB	14.5 cm
3	Royal Blue	PNWT/17/RB	PNWT/27/RB	16 cm
1	Red	PNWT/13/RD	PNWT/23/RD	11.5 cm
2	Red	PNWT/15/RD	PNWT/25/RD	14.5 cm
3	Red	PNWT/17/RD	PNWT/27/RD	16 cm



Junior Wrist/Thumb Brace

Product Code: **JWTO**

A wrist brace made of Neoprene which can be cut without fraying to fit around swollen joints, giving compression and excellent conformity.

- Incorporating a malleable metal thumb stay to control an unstable CMC joint; the thumb section can be trimmed to allow IP flexion
- Use for: Juvenile chronic arthritis, soft tissue injury.

Size	Code (Left)	Code (Right)	Wrist Circumference	Palmar Bar Length
S	*/13	*/23	12 cm	14 cm
M	*/15	*/25	14 cm	14 cm
L	*/17	*/27	16 cm	16 cm

* Insert product code



Tracy Brace

Product Code: **TWB**

A wrist brace made of quality woven elastic for comfort and durability, enclosing an aluminium palmar bar which can be easily removed.

- The design is cut below the palmar crease and will not impinge on the MCP joints allowing for a fully functional hand
- Use for: Colle's fracture, wrist fracture and arthritis
- Length of wrist brace is 14-15cm

Size	Code (Left)	Code (Right)	Wrist Circumference
S	*/13	*/23	12 cm
M	*/15	*/25	14 cm
L	*/17	*/27	16 cm

Contains Latex

* Insert product code



Paediatric Ulnar Deviation Support

Product Code: **PUDS**

- Manufactured in colourful Neoprene.
- Contrasting multi-coloured hook/loop closure straps.
- Hard wearing palmar bar casing with adjustable aluminium palmar bar.
- Adjustable ulnar support can be shaped to help control ulnar drift.

Size	Code (Left)	Code (Right)	Wrist Circumference
Baby	*/11	*/21	6.5 cm
XS	*/12	*/22	10 cm
S	*/13	*/23	12 cm
M	*/15	*/25	14 cm
L	*/17	*/27	16 cm

* Insert product code



Paediatric Lancaster Shoulder Brace

Product Code: **LSSB/10** (Baby) , **LSSP/10** (Paediatric)

A universal adjustable sling with shoulder strap and body wrap for shoulder immobilisation.

- Arm envelope can be easily trimmed for custom-fit.
- Fully adjustable shoulder strap with padded sleeve for comfort around neck.
- This sling can be easily adjusted to fit most paediatric patients.

Size	Code	Age
Baby	LSSB/10	0 - 18 Months
Paediatric	LSSP/10	2 - 10 Years



Paediatric Trauma Sling

Product Code: **TAU/P**

The Trauma Sling is a high quality universal sling that has multiple uses within the healthcare field. The Trauma Sling can be used as a traditional sling or for arm elevation following trauma, burns or post-operative management of the shoulder or upper arm.

The design and quality of materials used offer protection, comfort and ease of access to the damaged limb. The Trauma Sling is easy to apply and is supplied vacuum packed to reduce storage volume.



Size
One size

NEW

KIDS



Paediatric Walker

Product Code: **PJWF**

A comfortable paediatric walker boot with adjustable uprights for custom fit on lower limb, providing improved stability and immobilisation.

New lightweight design sole unit to greatly improve shock absorption and gait pattern. Robust foot section for improved wear.

Simple, secure hook / loop fastening. Comfortable detachable liners with replacements available.

Product Features

- Lightweight & durable design
- Improved stability and support
- No-slip, rocker sole



No-slip tread

Size	Code	Child UK Shoe Size
S	PJWF/03	2 - 6
M	PJWF/05	7 - 10
L	PJWF/07	11 - 13



Rocker Sole



Protem Belt

Product Code: **TBP**

This paediatric lumbar support has been used effectively in the management of abnormal tone by providing proximal stability and improving proprioceptive feedback.

Contains Latex

Size	Code	Waist Circumference
XXS	*/01	46 - 54 cm
XS	*/02	56 - 64 cm

* Insert product code



Pavlik Harness

Product Code: **PAV**

To control hip flexion and abduction in the primary treatment of congenital dislocation of the hips. Colourful strap identification with adjustable hook and loop closures, a padded harness with a unique reversible toe cover which can be worn open or closed toe.

Size	Code	Chest Circumference	Age of Baby
XXS	*/01	28 - 34 cm	Prem
XS	*/02	30 - 36 cm	0 - 3 Months
S	*/03	36 - 42 cm	3 - 6 Months
M	*/05	39 - 45 cm	6 - 9 Months
L	*/07	46 - 55 cm	9 - 12 Months
XL	*/08	56 - 59 cm	12 Months +

* Insert product code



CDH Splint

A lightweight, rigid polypropylene frame is lined with a hygienic, non-absorbant, closed cell foam padding.

Use to maintain legs in 90 degrees hip flexion and up to 60 degrees of abduction.

Size	Code	Age of Baby
XS	FCDH/02	0 - 3 Months
S	FCDH/03	3 - 6 Months
M	FCDH/05	6 - 12 Months
L	FCDH/07	12 - 18 Months
XL	FCDH/08	18 + Months



Fun Pattern Arm & Knee Immobiliser

Product Code: **FPI**

- Hard wearing fun print material
- Wrap-around with hook / loop closures
- Three rigid stays for effective immobilisation
- Please quote required length of brace after product code
- After each code please specify choice of fabric from below.
- Patterns may vary, but will be age/gender appropriate.

Size	Code	Length	Cicumference (Top)	Cicumference (Bottom)
1	FPI8/02/*	20 cm	17 - 24 cm	12 - 19 cm
2	FPI8/03/*	20 cm	23 - 30 cm	14 - 21 cm
3	FPI10/02/*	20 cm	17 - 24 cm	12 - 19 cm
4	FPI10/03/*	25 cm	23 - 30 cm	14 - 21 cm
5	FPI12/03/*	30 cm	23 - 30 cm	14 - 21 cm
6	FPI12/05/*	30 cm	30 - 37 cm	17 - 24 cm
7	FPI14/05/*	35 cm	30 - 37 cm	17 - 24 cm
8	FPI16/07/*	40 cm	37 - 44 cm	20 - 27 cm
9	FPI18/07/*	45 cm	37 - 44 cm	20 - 27 cm
10	FPI18/08/*	45 cm	44 - 51 cm	23 - 30 cm
11	FPI20/07/*	50 cm	37 - 44 cm	20 - 27 cm
12	FPI20/08/*	50 cm	44 - 51 cm	23 - 30 cm

* Add one of the codes below for fabric selection:
 (e.g. Size 2, with fabric code 'K' = FPI8/03/K)

Fabrics available:

Code E	Code F	Code G	Code H
Code I	Code J	Code K	Code L
Code C			

NEW



Wee Bow Post-Op

Product Code: **AE028100**

Locked or limited motion control of the elbow during rehabilitation, after operative procedures or injury to the elbow. Malleable arm cuffs allow for adjustability and improved fit.

Common Examples of Use

- Ligament injuries of the elbow (UCL, RCL, Annular Ligament)
- Strains and sprains

Product Features

- Universal sizing
- Telescoping bars adjust from 10.25" – 13.25" (26 – 34 cm)
- Full range-of-motion hinge adjustable from 10° – 110°

Code	Description
AE028100	Left, Universal
AE028200	Right, Universal

NEW



Wee ROM Post-Op Knee

Product Code: **EK090000**

A post-op knee brace with telescoping bars to accommodate varying leg lengths.

Common Examples of Use

- Non-operative or post-op procedures to the lower leg
- ACL, PCL, MCL and LCL injuries
- Meniscus
- Cartilage
- Juvenile Osteochondritis Dissecans
- Patella injuries

Product Features

- Universal sizing
- Telescoping bars adjust from 16" – 22" (41 – 56 cm)
- Full range-of-motion hinge adjustable from 10° – 110°
- Extends to malleolus

Code	Description
EK090000	Universal



Paediatric Hi Sided Cast Sandal

Product Code: **HSCSL**

Designed with hard wearing, non-stretch nylon faced padded upper for long term wear. Loop-lock D-ring fastening for secure and easy fastening. Non-slip rocker sole.



Size	Code	Length	Width
S	*/CS	14 cm	5.5 cm
M	*/CM	17 cm	6.5 cm
L	*/CL	20 cm	7.5 cm

* Insert product code

Paediatric Cast Sandal

Product Code: **CS**

- Strong and durable
- Canvas upper with jolly print
- Grooved non-slip rocker sole unit
- Hook / loop fastening



Size	Code	Length	Width
XXXS	*/98	12.5 cm	9 cm
XXS	*/99	14 cm	9 cm

* Insert product code

Contents

Hand Therapy

Code	Name	Pg.	Code	Name	Pg.
ATW	Abducted Thumb Wrap.....	19	PFLB	Poroflex Lace	9
ASX	Air SX Wrist Brace	30	PFWB	Poroflex Pull.....	9
LPS	Air X Low Profile Spica.....	29	PFTB	Poroflex Pull & Thumb	9
AWB	Air X One Piece Wrist.....	13	PRHS	Pro-rheuma Thumb Spica	21
ATRS	Air X Thumb Restriction Splint.....	12	PRHW	Pro-rheuma Wrist Brace	10
ATS	Air X Thumb Spica.....	13	PRHT	Pro-rheuma Wrist Thumb Brace	11
AFS	Air X Two Piece Wrist.....	13	PDR	Procool D-ring Wrist Splint	25
AWT	Air X Wrist Thumb.....	12	PDWT	Procool Deluxe Wrist Thumb Brace	25
BWTB	Bedford Wrist Thumb Brace.....	18	PTRS	Procool Thumb Restriction Splint.....	24
BWB	Beta Wrist Brace	6	PWS	Procool Wrist Support.....	24
SSX	Black Seal SX Wrist Brace.....	30	PWTRS	Procool Wrist Thumb Restriction.....	25
CAP/5	Cambridge Finger Splint.....	32	PWTS	Procool Wrist Thumb Support.....	30
CSX	Carbon SX Wrist Brace	30	TS	Promedics Thumb Shield.....	20
DFP	Deltaform PLUS Wrist Brace.....	7	31011	Push CMC	31
DFS	Deltaform Wrist Brace	7	RWB	Reversible Wrist.....	27
CAP/3	Double Armchair	32	CAP/4	Simpson Modified Armchair.....	33
DBF	Double Bedford Finger Stall.....	33	CAP/2	Single Armchair	32
ESNB	Economy Standard Neobrace.....	11	SWB	Standard Wrist Brace.....	6
ETS	Elastic Thumb Spica.....	20	TUJ	Thumbs Up.....	15
ETSM	Elastic Thumb Spica - Modified.....	20	CAP/8	Triple Armchair.....	32
WTE	Elastic Wrist Thumb Brace	17	TPNS	Two Piece Neoprene Spica	28
FSX	Freedom SX Wrist Brace.....	8	TRSV	Ventilated Thumb Restriction Splint.....	27
FWB	Freedom Wrist Brace.....	8	VNB	Ventoprene Neobrace.....	16
JRS	Functional Resting Splint	31	VWNB	Ventoprene Wrist Brace	16
IWB	Intrinsic Wrist Brace.....	15	WFS	Whitechapel Finger Splints	33
JRWTB	Jura Black Wrist / Thumb Brace.....	26	WNB	Wrightington Neobrace.....	16
JRWB	Jura Black Wrist Brace	26	WWB	Wrightington Wrist Brace	6
JLFW	Jura Latex Free Wrist Brace.....	8	WTN	Wrist Thumb Brace.....	18
JPC	Jura Palm Protector	15	WTW	Wrist Thumb Wrap	19
JUJDS	Jura Ulnar Deviation Support	22	WW	Wrist Wrap	19
JUWB	Jura Universal Wrist Brace	28	UDSA	Wrist/Ulnar Deviation Support.....	17
JWB	Jura Wrist Brace.....	14	WL	Wristlet.....	17
JWT	Jura Wrist Thumb.....	14			
JWTL	Jura Wrist Thumb (Long).....	14			
JLWW	Juraprene Long Wrist Wrap - Sizes.....	23			
JSWWU	Juraprene Short Wrist Wrap - Universal.....	23			
JTWWS	Juraprene Wrist Thumb Wrap.....	22			
JTWWL	Juraprene Wrist Thumb Wrap - Long.....	23			
BLTS	Lycrafleece Thumb Spica.....	10			
BLWB	Lycrafleece Wrist Brace	10			
MTW	Mini Thumb Wrap	28			
NNWB	Neobrace PLUS Wrist Brace	29			
SNB	Neobrace Wrist Brace	11			
NTS	Neoprene Thumb Spica	22			
WTO	Neoprene Wrist Thumb Brace	18			
NUF	Nuform Wrist Brace.....	7			
P-1008	Oval - 8 Finger Splint	34			
MFS	Plastic Mallet Splint.....	33			
LPF	Pollo Form Spica - Long.....	21			
SPF	Pollo Form Spica - Standard	21			

Cervical & Upper Limb

Code	Name	Pg.
AC	AC Collar	36
00052X	Atlas Universal Shoulder Brace	42
BOS	Boscombe Sling	39
AESB	Bradford Arm Elevator Sling	38
BC	Budget Collar	36
CLS	Clavicle Splint	44
CTS/10	Cubital Tunnel Splint	45
DSI	Deluxe Shoulder Immobiliser	38
ESI	Economy Shoulder Immobiliser	38
EAS	Evans Abduction Sling	43
FC	Halton Collar	36
HEW	Hinged Elbow Wrap	45
JASS	Jura Abduction Sling	43
JASSP/10	Jura Abduction Sling PLUS	41
JUC	Jura Clasp	46
JCS/10	Jura Clavicle Splint	44

JSI	Jura Shoulder Immobiliser	41
JSL/10	Jura Shoulder Lock	41
NSB	Neoprene Shoulder Brace	44
OES	Orthoprene Elbow Wrap	45
TAU/P	Paediatric Trauma Sling	40
PCC	Philly Collar	36
RVP	Replay Tennis Elbow Support	46
SSW	Sling and Swathe	39
004X	SlingShot™ 3 Shoulder Brace	42
SCF	Subluxation Cuff	44
0725x	T Scope Elbow Premier	45
TEP	Tennis Elbow Strap	45
TAU/10	Trauma Sling	40
UCS	Universal Combo Sling	43
LSS	Universal Lancaster Shoulder Sling	37
UPS	Universal Pro Sling	39
UTP/10	Universal Two Piece Collar	36
VPS	Vari Pad Tennis Elbow Support	46

Lower Limb

Code	Name	Pg.
ARK	Air X Long ROM Knee Brace.....	58
AMP	Air X Polycentric Knee Brace	58
ASK	Air X SK Knee Brace	58
AVK	Air X VK Long Knee Brace.....	59
AUT	Autolock II	61
PK428	Axiom-D Elite Ligament Knee Brace.....	49
BN	Buttress Knee Sleeve.....	63
CNI	Canvas Knee Immobiliser	57
ED312000	Dynamic Unloading Osteoarthritis (DUO) Knee Brace	52
EN	Economy Knee Sleeve.....	62
100379-1XX	Freerunner Knee Brace	51
1171X	Freestyle OA Knee Brace	50
1209	Fusion Knee Brace	53
1208	Fusion OA Osteoarthritis Knee Brace	50
GN/08	Genurange Lite.....	55
GN/0510	Genurange Universal	54
GN/04	Genuranger	54
GN/0310	Genurehab	55
HN	Hinged Knee Sleeve.....	63
8910	Hipo.....	64
BTPK	Jura Knee Immobiliser	56
OPK	Jura One Piece Immobiliser	57
022XX	LPR Knee Brace	51
BLPL	Lycrafleece Patella Lux	61
NS	Neosplint.....	57
NSK	Nuffield SK Knee Brace	60
PKS/20	Post-Op Knee Splint.....	56
M.2122X	Promedics Hip Joint.....	64
SK	SK Knee Brace.....	60
SP	Spiral Knee Sleeve.....	62
ST	Stabilised Knee Sleeve	63
SNO	Standard Knee Sleeve - Open Patella.....	62
8910	T Scope Hip Brace.....	64
08814	T Scope™ Premier Post-op Knee Brace	48
TPK	Three Panel Knee Immobiliser	56
AG061012	Thruster RLF.....	53
JBB/10	Universal Hip Abduction Belt with Shoulder Strap.....	64
PHJ/10	Universal Hip Abduction Waist Belt	64
VSKP	Ventilated SK Knee Plus	60
VK	VK Knee Brace.....	61

Ankle & Footwear

Code	Name	Pg.
AGA/10	Air/Gel Ankle.....	75
LEB	Ankle Contracture AFO	72
ANK	Anklet.....	77
L-478	Anti-slip Pad	88

EWB	Barnsley Ankle Brace.....	75
BPO	Basic Post-Op Shoe	83
BFS	Boot Foot Support	71
L-45	Cap.....	86
1144	Dorsal Night Splint.....	70
L-DS	Dress Shoe Comfort Inserts.....	84
L-46	Elastic Sheath	87
EPCB	Epoda Dual Heel Cups	85
EPB	Epoda Dual Heel Pads	85
F8LA	Figure of 8 Ankle Brace.....	76
821	Fixed Ankle.....	74
BL5150	Genesis Full Shell Walker	67
HVN	Hallux Valgus Night Splint	76
HUP/10	Heal-up™ Protection Boot	73
ANTFORTH	Heel Relief Orthosis.....	71
HRS	Heel Relief Shoe.....	80
HSCS	Hi Sided Cast Sandal.....	82
JWW	Jura Achilles Wedges.....	68
JWA	Jura Air Shield Walker	68
JWA-ES	Jura Air Shield Walker Extension Straps.....	68
JWF	Jura Fixed Walker.....	69
HP-OP	Jura High Post-Op Shoe	78
P-OP	Jura Post-Op Shoe.....	78
825	Leaf Spring Orthosis	74
LAB	Lincoln Ankle Brace.....	77
M.A1028	Malleostrom.....	76
MCD	McDavid Ankle.....	77
811	Medi Shoe	82
LEM	Multi Use Boot.....	72
OWH	Off-loader Heel Shoe.....	79
-	Off-loading Insoles	79
CS/08	Open Toe Cast Sandal	83
CS/30	Open Toe Cast Sandal with Buckles	82
CS/48	Open Toe Cast Sandal with Buckles	83
1136	Plantar Fasciitis Wrap	75
PB	Podus Boot	71
L-PQ6	PQ - 6mm Heel Lift.....	85
L-PQ9	PQ - 9mm Heel Lift.....	85
ASW	Promedics Airstep Walker	66
445	Rigid Post Op Shoe.....	84
JWR	ROM Walker	70
L-495	Silicone Gel Discs.....	87
L-4**	Silicone Gel Hallux Valgus Protector.....	88
L-47	Silicone Gel Metatarsal Pad with Elastic Band.....	86
L-41*	Silicone Gel Retractor / Protector	88
L-40*	Silicone Heel Cup	86
L-40*	Silicone Heel Cup with Soft Heel Spot	86
L-SHP	Star Heel Pad.....	85
L-SSP	Star Spur Pad.....	85
442	Supershoe	84
RPO	Superstep Post Op Shoe	81
M.9050	Swedish Ankle Foot Orthosis.....	74
L-43	Toe Retractor	87
L-44*	Tubular Protection A	87
L-44*	Tubular Protection B	87
UFO	Universal Plantar Faciitis Orthosis.....	72
JSHV	Universal Soft Hallux Valgus.....	76

Spinal Support

Code	Name	Pg.
600XX	BOA with Chairback Back Brace	93
CBB	Chairback Brace.....	90
CFB	Crossform Belt.....	92
LB*	Lumbarform Belt.....	94
MP45	MP45 Plus 3-Point Hyperextension Brace	90
NHYCO	Nhyco Cruciform Hyperextension Orthosis	90
ORT-TLS	Ortholux + Extender	91
PF*	Proflex Belt.....	94
PFM	Proform Belt.....	92
TB	Protem Belt.....	91
SLS	Sacro Lumbar Support.....	92
60002	TLSO Back Brace	93

Maternity Support

Code	Name	Pg.
DMB	Deluxe Maternity Belt	98
DUB	Duofit Maternity Belt	98
LPM	Lo-pro Maternity Belt.....	98
MBAS	Maternity Belt Style A	97
MBB	Maternity Belt Style B.....	97
MBS	Maternity Belt Style C.....	97
SPD	Pelvic Instability Belt	96
TROCH	Trochanteric Belt	96
WMB	Wiltshire Maternity Belt.....	96

Fracture Bracing & Splinting

Code	Name	Pg.
AJ/10	Alignment Jig for Genurange Hinges	102
APS	Alphastrap	105
DHC	Derby Heel Cup.....	103
EHC	Extended Humeral Fracture Brace.....	100
FEM	Femoral Sock	103
7130	Free Knee Polycentric Hinges.....	103
GNF	Genurange Fracture Bracing Hinges	102
GON/10	Goniometer	105
HS/10	Humeral Sock	103
HSO/10	Humeral Stabilising Orthosis.....	101
JPR	Juraprene Splinting Material	108
7135	Lerman Polycentric Hinges.....	103
LS	Loop Strap.....	105
	- Neofrakt.....	104
	- Orfit	107
PFP	Proform Plus	106
PFX	Proform XL.....	107
PHJ/10	Promedics Hip Joint.....	101
GNP/10	Replacement Headplates.....	102
GEN/CS	Replacement Hinge Cover Screws	102
GNS/10	Replacement Screws.....	102
SA	Self Adhesive Hook / Loop.....	108
SS	Sling strap.....	100
PSH	Soft Humeral	101
HAL	Standard Hook / Loop.....	108
SHC	Standard Humeral Fracture Brace.....	100
TS	Terry Strap	105
BK/10	Tibial Sock	103
JBB/10	Universal Hip Abduction Waist Belt	101
JBBS/10	Universal Hip Abduction Waist Belt with Shoulder Strap.....	101

Cold Therapy

Code	Name	Pg.
100291	Polar Care Kodiak Cold Therapy.....	110
CTKW	Universal Cold Therapy Wraps	110

Paediatric

Code	Name	Pg.
FCDH	CDH Splint.....	117
FPI	Fun Pattern Arm & Knee Immobiliser.....	118
JWTO	Junior Wrist/Thumb Brace.....	114
PNWT	Neoprene Wrist/Thumb Brace	114
CS	Paediatric Cast Sandal.....	120
PEW	Paediatric Elastic Wrist Brace.....	113
PEWT	Paediatric Elastic Wrist/Thumb Brace	113
HSCSL	Paediatric Hi Sided Cast Sandal.....	120
LSSP	Paediatric Lancaster Shoulder Brace	115
PNW	Paediatric Neoprene Wrist Brace.....	113
PLTS	Paediatric Thumb Spica.....	112
TAU/P	Paediatric Trauma Sling	115
PUDS	Paediatric Ulnar Deviation Support	115
PJWF	Paediatric Walker	116

PLWB	Paediatric Wrist Splint.....	112
PLWTB	Paediatric Wrist/Thumb Splint	112
PAV	Pavlik Harness.....	117
TBP	Protem Belt	117
TWB	Tracy Brace.....	114
AE028100	Wee Bow Post-Op	119
EK090000	Wee ROM Post-Op Knee.....	119

Terms & Conditions

These terms and conditions apply to all orders placed. Please read all information carefully. By placing an order with us you are accepting the Terms and Conditions set out herein. If you have any queries regarding these terms and conditions, please contact us by e mail at orders@promedics.co.uk or telephone on +44 (0)1475 746400 or fax 01475 746410 or in writing to Promedics Orthopaedic, Port Glasgow Industrial Estate, Gareloch Road, Port Glasgow, Scotland, PA14 5XH

"Conditions" means the terms and conditions as set out in the document.

"Product" means any Product supplied by Us in accordance with these Conditions.

"United Kingdom" means England, Wales, Scotland, Northern Ireland and the Channel Islands

"We/Us" means Promedics Limited and "Our" shall be construed accordingly.

"You" means any person ordering Products from Us.

1. Payment

1.1 All payment on approved credit accounts must be made within 30 days of the date of the invoice and are made subject to acceptance by You of Our terms and conditions.

Fulfillment of the order is subject to stock availability and acceptance of Your order and We are under no obligation to fulfill Your order.

1.2 Payments can be made with most major credit and debit cards.

1.3 If You experience any ordering difficulties or have any concerns please phone or e mail us. Please ensure that You leave the correct details.

2. Price

2.1 Prices shown by Us are accurate on the date they were published but are subject to change.

3. Currency

3.1 All prices quoted by Us are in pounds sterling and exclude VAT at the current rate.

4. Delivery

4.1 Standard timescales for delivery of Products to UK addresses by a nominated carrier chosen by Us are as follows:

4.1.1 1 -5 working days from receipt of order and payment of order for manufactured Products.

4.1.2 10 working days for any made to measure or specifically requested Products.

4.1.3 10 working days for any Breg, North Coast and Orfit Industries Products.

4.2 We will not be liable for late delivery of Products and time shall not be of the essence in these Conditions and We reserve the right to make partial deliveries.

4.3 Standard delivery charges (which include packaging and inspection of the Products prior to dispatch) is £5.95 but We reserve the right to increase this delivery charge where such an increase is beyond Our control. Contact Us if You have any special delivery requirements.

4.4 The title of the goods remains with Promedics Orthopaedic Ltd. until paid for in full. We retain the right to recover goods that have not been paid for in full and will charge for the cost of doing so.

5. Cancellation procedure for Consumers

5.1 You have the right to cancel the order except for made to measure Products within 7 working days of receipt of the Products provided that the Products are not damaged, opened or used.

5.2 All cancellations must be made in writing to Us with a completed Promedics Returns Note detailing the reason for return and received at Our address within the 7 days cooling off period.

5.3 If You do cancel within 7 days and in accordance with the procedure stated in 5.2 then You will be entitled to a full refund by Us of all monies debited by Us within 30 days of the date the notice of cancellation was given together with the reasonable costs to You of the return of the Products provided such return is a result of the Products being defective.

5.4 If You cancel Your order but do not comply with clause 5.2 then You will be entitled to a refund of the cost of the Products minus the cost for Us to collect the Products from You.

6. Returns

6.1 Except as provided for in Clause 5 above relating to consumers all cancellations or arrangements for the return of Products are at Your own cost.

6.2 If You receive Products that are defective please contact Us by e mail or telephone within 21 working days of receipt of the Products. In this case You will be entitled to a refund from Us limited to the cost of repair or (where appropriate in our opinion) replacement of the Product.

6.3 If You wish to return Products please e mail us at Our address listed above or the telephone listed at the top of these Conditions stating the reason for return.

We will then provide you with a Returns Note for You to complete and return with the packaging of all returned Products.

6.4 Where Products are returned for Credit We will charge a restocking fee as follows; 21 - 60 days 20%, 60 days - 6 months 40 %, 6 - 12 months 60%, Over 12 months - 100%.

7. Refunds

7.1 If You feel that You have been incorrectly charged for the Products that You have ordered or that Your goods are different to Your original order (except where it has been agreed in advance by You and Us) please contact Us by e mail or telephone.

7.2 If You comply with clause 7.1 You will be entitled to a full refund of incorrect amounts debited to Us or all monies debited by Us should You cancel the order within 30 days receipt by Us of Your complaint.

8. Liability

8.1 In relation to the purchase of Products We accept no liability for indirect or consequential loss or damage or for any loss of data, profit or revenue or business (whether direct or indirect) however caused.

8.2 We accept no liability for the loss however caused whether due to the inaccuracy, error, omission or any other cause or whether on Our part or that of Our servants, agents or any other person.

9. Jurisdiction

The Conditions shall be governed by and construed in accordance with the laws of England and You irrevocably submit to the exclusive jurisdiction of the courts of England.

10. Force Majeure

We shall not be liable to You where performance of any of Our obligations to You is prevented or impeded by any act of God, war and other hostilities, civil commotion, accident, strikes, lock outs, trade disputes, acts or restraints of Government or any other cause that is not reasonably within Our control.

11. Severance

If any part of the Conditions shall be deemed unlawful, void or for any reason unenforceable, then that provision shall be deemed to be severable from the Conditions and shall not effect the validity and enforceability of any of the remaining provisions of the Conditions.

12. Waiver

No waiver by Us of any provision of the Conditions shall be construed as a waiver of any proceeding, succeeding or continuing breach of any provision of the conditions.

13. Entire Agreement

These Conditions govern Our relationship with You. Any changes to these Conditions must be in writing and signed by both parties. You confirm that in agreeing to accept the Conditions You have not relied on any representation save insofar as the same has expressly been made a term of these Conditions and You agree that You shall have no remedy in respect of any representation.

14. Statutory Rights

Where You are a consumer as defined by the Distance Selling Regulations 2000 Your statutory rights are not affected by these Conditions.

Promedics Orthopaedics Ltd.

Port Glasgow Industrial Estate
Gareloch Road, Port Glasgow
Scotland PA14 5XH

Solutions for today and tomorrow

T: 01475 746 400
F: 01475 746 410
E: enquiries@promedics.co.uk

twitter.com/promedics 

facebook.com/promedicsorthopaedic 

youtube.com/promedicsorthopaedic 

www.promedics.co.uk

