

Lacer Ankle

The Lacer Ankle is a low-profile, comfortable brace, offering stability in a discreet, yet supportive design that is suitable for sports. The brace features conformable stays for stability, especially for patients who have previously suffered ankle injury. The laced enclosure design is easy to apply with figure of 8 strapping and an elastic cuff for adjustable compression. Provides ankle support for an array of injuries and conditions, such as;

- + Sprains
- + Medial / lateral instability
- + Subtalar joint instability



The hook/loop figure of 8 closure system allows for easy fastening



Low profile ankle orthosis made from a comfortable, breathable material



Elastic cuff for adjustable compression around the leg



Flexible plastic stays conform to the malleoli for stability

Size	Code	Ankle Circumference	UK Shoe Size
XXS	LCRA/01	22cm - 25cm	3
XS	LCRA/02	25cm - 28cm	4 - 6
S	LCRA/03	28cm - 30cm	6.5 - 8
M	LCRA/05	30cm - 33cm	8.5 - 10
L	LCRA/07	33cm - 35cm	10.5 - 12
XL	LCRA/08	35cm - 38cm	12.5 - 14
XXL	LCRA/09	38cm - 40cm	14+