

Philly Collar

A lightweight, two-piece collar in Plastazote, with rigid occipital and mandibular posts. Strong comfortable support for the head and neck, incorporating a tracheal access hole.

The Philly Collar is designed to provide support and relief for;

- + C-Spine precautions for trauma patients
- + Immobilisation for pre c-spine surgery
- + Other conditions which require gross immobilisation of the mid cervical spine



Size	Code	Neck Circumference	Collar Depth
S/Shallow	PCC/33	25.5cm - 30.5cm	8.5cm
M/Shallow	PCC/53	33cm - 38cm	8.5cm
L/Shallow	PCC/73	41cm - 48cm	8.5cm
S/Mid	PCC/34	25.5cm - 30.5cm	11cm
M/Mid	PCC/54	33cm - 38cm	11cm
L/Mid	PCC/74	41cm - 48cm	11cm
S/Deep	PCC/35	25.5cm - 30.5cm	13.5cm
M/Deep	PCC/55	33cm - 38cm	13.5cm
L/Deep	PCC/75	41cm - 48cm	13.5cm

5.5cm depths available on request.