

Paediatric Fixed Walker

A comfortable, one-piece injection moulded walker boot with adjustable 'paddle' uprights to reduce pressure on lower limb and provide improved stability and immobilisation.

The lightweight design sole unit greatly improves shock absorption and the gait pattern. Robust foot section improves wear. Simple, secure hook / loop fastening. Comfortable detachable liners with replacements available.

- + Foot and ankle conditions that benefit from motion restriction
- + Ankle sprains
- + Trauma and rehabilitation
- + Stable foot fracture



Rocker sole



No-slip tread sole



Size	Code	Children's Shoe Size (UK)
S	PJWF/03	2 - 6
M	PJWF/05	7 - 10
L	PJWF/07	11 - 13